

ToxCup[®] Drug Screen Cup Distributed by CLIAwaived[™] Inc.



ToxCup[®] Drug Screen Cup

The ToxCup[®] Drug Screen Cup is a one-step test for the rapid detection of up to seven drugs of abuse or five drugs and four adulteration tests. Its design integrates sample collection, drug/adulteration testing and shipment for confirmation in one device.

BENEFITS:

- Simultaneously tests for drugs and adulterants
- Rapid one-step procedure
- No sample exposure
- Integrates sample collection, testing and shipment for confirmation in a single container
- Easy to interpret visual results
- Room temperature storage
- Built-in temperature strip
- FDA Cleared and CLIA waived

Ordering Information:	
Product Description	Catalog #
ToxCup [®] Drug Screen Cup – 5 Panel (25 Cups) (THC, COC, OPI, AMP, PCP)	BRN-PT-11
ToxCup [®] Drug Screen Cup – 5 Panel (25 Cups) (THC, COC, OPI, AMP, PCP + CR/NI/PH/OX)	BRN-PT-11A
ToxCup [®] Drug Screen Cup – 5 Panel (25 Cups) (THC, COC, OPI, PCP, MET)	BRN-PT-12
ToxCup [®] Drug Screen Cup – 5 Panel (25 Cups) (THC, COC, OPI, PCP, MET + CR/NI/PH/OX)	BRN-PT-12A
ToxCup [®] Drug Screen Cup – 5 Panel (25 Cups) (THC, COC, OPI, MET, AMP)	BRN-PT-15
ToxCup [®] Drug Screen Cup – 5 Panel (25 Cups) (THC, COC, OPI, MET, AMP + CR/NI/PH/OX)	BRN-PT-15A
ToxCup [®] Drug Screen Cup – 5 Panel (25 Cups) (THC, COC, OPI, MET, BZO + CR/NI/PH/OX)	BRN-PT-23A
ToxCup [®] Drug Screen Cup – 5 Panel (24 Cups) (THC, COC, OPI, MET, OXY)	BRN-PT-24
ToxCup [®] Drug Screen Cup – 6 Panel (25 Cups) (THC, COC, MET, OPI, OXY, MTD)	BRN-HT-14
ToxCup [®] Drug Screen Cup – 7 Panel (25 Cups) (THC, COC, OPI, MET, BZO, BAR, OXY)	BRN-HT-21
Reimbursement Information:	
CPT Code #	National Limit Amount
80101QW (Single Class of Drug Qualitative)	\$19.24

ToxCup® Drug Screen Cup Test Instructions

The ToxCup® Drug Screen Cup is a one-step test for the rapid qualitative detection of drugs of abuse and adulterants.

ToxCup® Procedures -- Refer to package insert for detailed instructions and technical information.

1. Remove the test lid from the sealed pouch.




2. Twist the test lid securely onto the specimen cup after collection.




3. Lean the cup on its side to activate testing.



4. Adulteration Pads: After one minute, compare the colors of the adulterant tests with the corresponding color blocks on the color chart included in the kit. **Do not read after 2 minutes as the color may change.**



5. Drugs of abuse test: Read results once the **control band (C)** appears and the membrane clears (in approximately 5 minutes or less). Results are stable for up to an hour after the control bands form. Interpret according to instructions below.



Interpreting ToxCup® Test Results

Negative Result

For each test, two colored bands should be observed:

- **One in the CONTROL (C) region**
- **One in the specific TEST region**

The color of the test band may be slightly darker or lighter than the control band. Any band that can be seen visually is a **negative** result.




Example Interpretation:

AMP: Negative
COC: Negative
MET: Negative
OPI: Negative
THC: Negative

Presumptive Positive Result

A colored band at the **CONTROL (C) region** should be observed. When there is no colored band at the specific **TEST region**, the test is presumptive positive for that particular drug.



Example Interpretation:

AMP: Negative
COC: **Presumptive Positive**
MET: Negative
OPI: Negative
THC: **Presumptive Positive**

Invalid Result

When no band appears in the **CONTROL (C) region**, the **TEST is invalid** even if there is a band in the test region. Check test procedures, and samples. Repeat the test with a new device.



Example Interpretation:

AMP: **Invalid**
COC: Negative
MET: Negative
OPI: **Invalid**
THC: Negative

Rev. 10.10.09