Eurotrol GlucoTrol-NG

Intended purpose/Intended use

Eurotrol GlucoTrol-NG is an assayed glucose control intended for professional use in the verification of the precision and accuracy of HemoCue® Glucose systems.

IVD Medical Device

Eurotrol GlucoTrol-NG complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark.

Eurotrol GlucoTrol-NG complies with the following US codes of Federal Regulations (CFR): 42 CFR part 72 and 21 CFR parts 606, 640 and 820.

Eurotrol GlucoTrol-NG is for in vitro diagnostic use only.

Summary and Principle

Eurotrol GlucoTrol-NG is a quality control material with a known glucose concentration. Eurotrol GlucoTrol-NG may be used in the periodic verification of the precision and accuracy of the HemoCue Glucose systems. Please refer to local guidelines for recommended frequency of use.

Reagents

Eurotrol GlucoTrol-NG is available at different physiologically relevant levels, each vial containing 1.0 mL of purified bovine plasma with glucose.

Storage and stability

Eurotrol GlucoTrol-NG must be stored in the refrigerator at 2–8 °C (35–46 °F). Stored unopened at this temperature it is guaranteed stable as indicated until the expiry date on the ampoules and outer box. After opening the vial, it is stable for 30 days when properly recapped and stored at 2–8 °C (35–46 °F).

Procedures/instructions for use

1. Allow the vial to stand for 15 minutes at room temperature 15–30°C (59–86°F).
2. Gently mix the vial 8 – 10 times before sampling.
3. Do not fill the cuvette directly from the vial. Dispense a drop of the control material onto a hydrophobic surface, for example a plastic film. Fill the cuvette according to the manufacturer's instructions. Place the cuvette in the analyzer for immediate analysis.
4. Wipe any excess material from the vial and the cap with a clean tissue. Recap the vial tightly.
5. Store in the refrigerator at 2–8 °C (35–46 °F)
6. If the controls do not perform as expected, review the instructions for use to see if the test was performed correctly. Check the expiration date and storage conditions for the control and the cuvettes. Repeat the test. If the controls still doesn't perform as expected, contact technical assistance.

Precautions

1. For in vitro diagnostic use only.
2. Caution: Eurotrol GlucoTrol-NG is an Animal Blood Product. Bovine based materials do not carry biohazards for man, such as Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). This product is free from TSE. However, no test method

can offer complete assurance that products derived from human source material are free from infectious agents.

3. It is good laboratory practice to follow the same precautions as with patient samples when handling or disposing of vials.
4. This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
5. Eurotrol GlucoTrol-NG is not to be used as a calibrator.

Expected results

The mean assay value provided on the enclosed assay sheet is derived from replicate analyses on factory calibrated HemoCue Glucose systems. The values are only valid on the HemoCue Glucose systems as referred to in the assay sheet. Values for "HemoCue whole blood systems" are applicable for HemoCue Glucose systems reporting whole blood values. Values for "HemoCue plasma equivalent systems" are applicable for HemoCue Glucose systems reporting plasma equivalent values. Results obtained should fall within the expected range specified for the control. Not all HemoCue Glucose systems are available world wide.

For further information or technical assistance, please contact:

Eurotrol B.V., The Netherlands	Eurotrol Inc., USA
Phone +31 318 695 777	Toll free phone 1-866-234-5754
Fax +31 318 695 770	Main phone 1-978-598-3779
E-mail: office@eurotrol.com	Fax 1-978-598-3780
www.eurotrol.com	E-mail: officeUSA@eurotrol.com
	www.eurotrol.com

 Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, The Netherlands

Symbols used:

Attention, see instructions for use

Use by

In Vitro Diagnostic Medical Device

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Eurotrol B.V., Pays-Bas
Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, Pays-Bas