Intended purpose/Intended use

Eurotrol AlbuTrol is an assayed albumin control intended for professional use in the verification of the precision and accuracy of the HemoCue® Albumin 201 system.

IVD Medical Device

Eurotrol AlbuTrol complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark.

Eurotrol AlbuTrol complies with the following US codes of Federal Regulations (CFR): 42 CFR part 72 and 21 CFR parts 606, 640 and 820.

Eurotrol AlbuTrol is for in vitro diagnostic use only.

Summary and Principle

Eurotrol AlbuTrol quality control material is available at two levels, each with a known albumin concentration. Eurotrol AlbuTrol may be used in the periodic verification of the precision and accuracy of the HemoCue Albumin 201 system. Please refer to local guidelines for recommended frequency of use.

Reagents

Eurotrol AlbuTrol is provided at two physiologically relevant levels. Each vial contains 1.0 mL of human urine containing albumin.

Storage and stability

Eurotrol AlbuTrol must be stored in the refrigerator at 2–8 °C (35–46 °F). Stored unopened at this temperature it is guaranteed stable until the expiration date as indicated on the vial and the outer box. After opening the vial, it is stable for 30 days stored at 2–8 °C (35–46 °F).

Procedures/instructions for use

1. Allow the vial to stand for 15 minutes at room temperature 15–30°C (59–86°F).
2. Gently mix the vial 8 – 10 times before sampling.
3. Do not fill the cuvette from the vial. Dispense a drop of the control material onto a hydrophobic surface, for example a plastic film. Fill the cuvette according to the manufacturer's instructions. Place the cuvette in the analyzer for immediate analysis.
4. Wipe any excess material from the vial and the cap with a clean tissue. Recap the vial tightly.
5. Store in the refrigerator at 2–8 °C (35–46 °F)
6. If the controls do not perform as expected, review the instructions for use to see if the test was performed correctly. Check the expiration date and storage conditions for the control and the cuvettes. Repeat the test. If the controls still doesn't perform as expected, contact technical assistance.

Precautions

1. For in vitro diagnostic use only.

2. Caution: Eurotrol AlbuTrol is prepared from human urine and contains other human source material. Components of the control, which are derived from human source material, have been tested and found to be non-reactive for HBsAg, anti-HCV and anti-HIV.
3. However, no test method can offer complete assurance that products derived from human source material are free from infectious agents.
4. This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
5. Eurotrol AlbuTrol is not to be used as a calibrator.

Expected results


The mean assay value provided on the outer box is derived from replicate analyses on factory calibrated HemoCue Albumin 201 systems. Results obtained should fall within the expected range specified for the control. The "Alb-201" abbreviations used on the outer box and vial define the results for the HemoCue Albumin 201 system.

For further information or technical assistance, please contact:

Eurotrol B.V.	Eurotrol Inc.
Postbus 722	563 Main Street
6710 BS Ede	Bolton, MA 01740
The Netherlands	USA
Phone +31 318 695 777	Toll free phone 1-866-234-5754
Fax +31 318 695 770	Main phone 1-978-598-3779
E-mail: office@eurotrol.com	Fax 1-978-598-3780
www.eurotrol.com	E-mail: officeUSA@eurotrol.com
	www.eurotrol.com

 Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, The Netherlands

Symbols used:

- | | | | |
|---|-------------------------------------|---|------------------------|
|  | Attention, see instructions for use |  | Batch number |
|  | Use by |  | Temperature limitation |
|  | In Vitro Diagnostic Medical Device |  | Reference number |
|  | Manufacturer |  | CE Mark |

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Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, Pays-Bas