

SureStep[®] Technology

LifeScan makes getting accurate glucose results perfectly easy.

The **excellent performance** of LifeScan's SureStep Technology* is a reflection of our **commitment to accuracy**, the most important feature of any blood glucose monitoring device.

SureStep[®] Test Strips and Meters are painstakingly engineered to be **accurate** and **precise**, even at low glucose levels, and to **avoid interferences** from other substances.

Our innovative test strip holder not only helps to **promote proper strip insertion** and accurate readings, it also **repels blood** to make cleaning easier.

We recognize that great products must be **easy to use**, so we created a **touchable** test strip with an easy-to-see **pink test square** and a **confirmation dot**. And since SureStep[®] Meters perform an internal check each time they are turned on, **no check strip** is required.

* The SureStep[®] Brand is marketed as Novo Assist[™] or GlucoTouch[®] in certain international locations.

SureStep® Technology for the entire continuum of care



SureStep®

COMPLETE BLOOD GLUCOSE MONITORING SYSTEM

- Simple, single-button testing
- Results in as little as 15 seconds
- Confirms enough blood applied
- 150-test memory with date and time
- Compatible with IN TOUCH® Diabetes Management Software



SureStep®Pro

PROFESSIONAL BLOOD GLUCOSE MANAGEMENT SYSTEM

- Capture and control point-of-care testing data with confidence
- Promotes caregiver and patient safety
- Helps make it easier to comply with stringent regulatory requirements
- Integrates clinical, regulatory, and utilization data from bedside to LIS
- DataLink™ enhancements provide options for remote data transfer via modem or network and a LIS interface via script or EDI



SureStep®Flexx

PROFESSIONAL BLOOD GLUCOSE MANAGEMENT SYSTEM

- Small, lightweight ergonomic design
- Enhanced QC and operator compliance
- Comprehensive data management and report capabilities
- Available with all DataLink connectivity options

LifeScan's SureStep Meters provide innovative blood glucose management solutions for the entire continuum of care, from self-monitoring for individuals to integrated professional systems for clinical and institutional settings.

SureStep® Meters offer key advantages


- Innovative SureStep® Test Strip Technology
- Touchable test strip for true off-meter dosing
- Plasma-calibrated glucose results
- Photoreflectance optical system
- Glucose oxidase chemistry for minimal interference

The advanced design of SureStep® Technology provides accurate, dependable glucose results in three simple steps

ONE Sample is applied to the test strip

Absorbent pad


- Touchable
- Special hydrophilic polymer foam promotes quick, convenient, and safe sample application
- Blood is retained and not transferred to other surfaces, thereby reducing the potential for contamination
- Sample flows one way through the porous pad to the reagent membrane, where the reaction occurs




Asymmetric black tip

- Asymmetric shape helps with the proper orientation of test strip into test strip holder
- Conductive black ink provides an indication that the test strip has been inserted and remains inserted by completing an electrical circuit

Test strip at a glance



Blood sample flow



TWO Glucose reacts with reagents in the test strip


Reagent membrane

- Reagent-impregnated microporous membrane provides the reaction site
- Filters out red blood cells while allowing plasma to move through
- Hides red blood cells and simplifies glucose measurement

Enzymatic reaction

- Glucose in the sample is oxidized by glucose oxidase (GO) in the presence of atmospheric oxygen, forming hydrogen peroxide (H₂O₂)
- H₂O₂ reacts with indicator dyes using horseradish peroxidase (HRP), forming a chromophore or light-absorbing dye
- The intensity of the color formed at the end of the reaction is proportional to the glucose present in the sample

Reaction chemistry



Glucose oxidase chemistry

- High degree of specificity for D-glucose
- Very stable
- Efficiently couples to visible dye systems

Patented indicator dye system

- Tailored specifically to absorb a particular wavelength of light where hemoglobin (red blood) does not interfere
- Stable over time
- High level of performance for accuracy and precision is achieved, especially at low glucose concentrations


Confirmation dot

- Turns blue after the blood sample has been applied
- Allows user to see if a sufficient amount of blood has been applied correctly
- Intensity of the blue color can be compared to the optional color chart for quick reference of blood glucose concentration

Endpoint reaction

- Glucose concentration is determined using the amount of color formed at the end of the reaction: the point at which all the glucose has reacted with the enzyme and the color stops changing
- Off-meter dosing is made possible, thereby improving infection control
- User has up to two minutes to insert the test strip once sample has been applied

Test Strip vial color charts



Vial color charts are to be used as an additional check to confirm a meter reading. They are NOT a replacement for a meter test.

THREE Measurement of blood glucose concentration with SureStep® Meters




No check strip is required

- SureStep Meters automatically perform numerous internal checks each time they are turned on

Test strip holder

- Patented design with serpentine path, which in conjunction with the unique asymmetric test strip tip, helps guide proper test strip insertion
- Teflon®-impregnated plastic repels blood



Measurement of glucose level using reflectance photometry

Reflectance photometry quantifies the intensity of the colored product generated by the enzymatic reaction

- Light-emitting diode (LED) emits a specific wavelength of light onto the test strip
- Colored product absorbs the light—the more glucose in a sample (and thus the more colored product on a test strip), the less reflected light
- A detector captures the reflected light, converts it into an electronic signal, and translates it into a corresponding glucose concentration
- System is calibrated to give plasma glucose values
- In general, reflectance photometry methods are less sensitive to the interferences common to electrochemical technologies

Performance data

EXCELLENT ACCURACY

Accuracy is a measure of how close the results are to target reference values.


Testing by lay user¹

Linear regression

$$y = 1.01x - 1.02 \text{ (mg/dL)}$$

$$1.01x - 0.56 \text{ (mmol/L)}$$

$r^2 = 0.969$
 $r = 0.984$
 $n = 335$
 $Sy.x = 13.46 \text{ mg/dL}$
 0.748 mmol/L




Testing by laboratory technician¹

Linear regression

$$y = 1.03x - 4.44 \text{ (mg/dL)}$$

$$1.03x - 0.56 \text{ (mmol/L)}$$

$r^2 = 0.967$
 $r = 0.983$
 $n = 335$
 $Sy.x = 14.20 \text{ mg/dL}$
 0.789 mmol/L



EXCEPTIONAL PRECISION

Precision is a measure of reproducibility.

Within-run precision of SureStepFlexx using whole blood²

Glucose concentration	Mean		SD		CV %
	mg/dL	mmol/L	mg/dL	mmol/L	
Site 1 (n=20)					
Normal	78.4	4.36	1.1	0.061	1.4
High	339.3	18.85	11.2	0.622	3.3
Site 2 (n=20)					
Low	47.8	2.66	1.9	0.106	4.1
High	260.7	14.48	4.9	0.272	1.9

LACK OF INTERFERENCES

Interferences from endogenous substances³

Substances produced or originating from within the body.

NO INTERFERENCE was observed in the following substances at physiological levels up to test concentration:

Substance	Test concentration		Physiological level	
	mg/dL	mmol/L	mg/dL	mmol/L
Bilirubin	20	0.342	1.2	0.021
Cholesterol	500	12.93	300	7.76
Creatinine	30	2.65	1.5	0.133
Triglycerides	3000	33.88	190	2.15
Uric Acid	20	1.19	7	0.42

Interferences from exogenous substances³

Substances originating from outside the body such as therapeutic agents or vitamins.

NO INTERFERENCE was observed in the following substances at therapeutic levels up to test concentration:

Substance	Test concentration		Therapeutic level (or range)	
	mg/dL	mmol/L	mg/dL	mmol/L
Acetaminophen	20	0.342	1–2	0.066–0.132
Ascorbic Acid	3	0.17	0.8–1.2	0.045–0.068
Ephedrine	10	0.61	0.005–0.01	0.0003–0.0006
Ibuprofen	40	1.94	0.5–4.2	0.024–0.204
Methyl Dopa	2.5	0.12	0.1–0.5	0.005–0.024
Salicylate (aspirin)	50	3.62	15–30	1.09–2.17
Tetracycline	4	0.09	0.4	0.009
Tolazamide (Tolamide, Tolinase)	100	3.21	2.0–2.5	0.06–0.08
Tolbutamide (Oramide, Orinase)	100	3.70	5.3–10	0.20–0.37
Sodium Citrate	500	19.4	Blood drawing tube	
EDTA (K+)	400	9.89	Blood drawing tube	

The following substances produced NO INTERFERENCE AT THERAPEUTIC LEVELS up to indicated threshold (i.e., exhibited suprathreshold interference only):

Substance	Threshold Level		Therapeutic level (or range)	
	mg/dL	mmol/L	mg/dL	mmol/L
Dopamine	6	0.392	0.4–1.6	0.026–0.104
Gentisic Acid	10	0.649	3.5–5.0	0.23–0.32
L-Dopa (Levo-Dopa)	20	1.014	0.02–0.3	0.001–0.015

References

1. Stover HM, et al. *Clinical evaluation of the enhanced SureStep System for self-monitoring of blood glucose (SMBG)*. Presented at the AADE 26th Annual Meeting, August 1999. Reprint available as LifeScan AW 056-796.
2. Weber S, Torio M, Correll G, et al. *Multicenter evaluation of a new blood glucose monitoring system for point-of-care testing*. Presented at the 52nd Annual AACC Meeting, July 2000. Reprint available as LifeScan AW 057-263.
3. *Interferences from endogenous and exogenous substances*. LifeScan Technical Bulletin, AW 056-498, 1999.

For more information, check with your sales representative or call the LifeScan Healthcare Professional Line at 1 800 524-7226.

We designed SureStep® Products with you in mind

Each feature of SureStep® Technology is engineered to provide you with benefits in performance, ease of use or infection control.

Performance
Ease of use
Infection control

Technology	Feature	Benefits		
Absorbent pad	• Touchable		✓	
	• Easy-to-see pink sample target		✓	
	• Absorbs, retains, and transfers blood one-way to reagent membrane	✓		✓
Microporous reagent membrane	• Filters out and hides red blood cells to allow for visual confirmation	✓		
	• Glucose oxidase-based chemistry provides excellent accuracy and precision without interferences from other substances <ul style="list-style-type: none"> – very specific for D-glucose – very stable – efficiently couples to visible dye systems 	✓		
	• Patented indicator dye system provides excellent accuracy and precision, especially at low glucose concentrations <ul style="list-style-type: none"> – dye absorbs a specific wavelength of light where hemoglobin does not interfere – dye absorbs light in proportion to amount of color present in the sample – dye is relatively stable over time 	✓		
	• Not affected by PO_2	✓		
	• Off-meter sample application is possible since the glucose concentration is determined by the amount of color formed at the <i>end</i> of the reaction		✓	✓
Confirmation dot	• Turns blue almost immediately after proper application of blood		✓	
	• Provides a visual estimate of blood glucose level		✓	
Asymmetric black tip	• Asymmetric shape helps with the proper orientation of the test strip (pink side up) upon insertion into the test strip holder	✓	✓	
	• Completes an electrical circuit, providing an indication that a test strip has been inserted and remains inserted	✓	✓	
Test strip holder	• Teflon®-impregnated plastic repels blood to make cleaning easier			✓
	• Patented serpentine shape maintains the test strip against the optics, promoting a constant optical reading	✓		
	• If necessary, can be easily removed for cleaning		✓	✓
Calibration	• SureStep® Systems are calibrated to give plasma glucose values, enabling a more direct comparison to laboratory results		✓	
Check strip	• No check strip is required—SureStep® Meters perform an automatic check every time they are turned on	✓	✓	

Distributed by: **CLIAwaived.com™**
 San Diego, CA 92121
 tel 858-481-5031
 toll free 888-882-7739
 www.cliawaived.com