ENGLISH INSTRUCTIONS

A rapid, visual immunoassay for the qualitative detection of drug and drug metabolites in human urine. This product is used to obtain a visual, qualitative result and is intended for professional use. This booklet may be used for any single drug screen assay listed on page 1. Refer to the information specific to the drug screen assay of the device used.

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For In Vitro Diagnostic Use Only

INTENDED USE

The SureStep™ Drug Screen Tests are an *in vitro*, lateral flow, one-step competitive immunoassay for the qualitative detection of drug and drug metabolites in human urine.

Abbreviation		Cut-off
AMP	Amphetamine	1,000 ng/ml
AMP II	Amphetamine II	300 ng/ml
BAR	Barbiturates	300 ng/ml
BZO	Benzodiazepines	300 ng/ml
COC	Benzoylcegonine	300 ng/ml
COC II	Benzoylcegonine II	150 ng/ml
MDMA	Ecstasy	500 ng/ml
M-AMP	Methamphetamine	1,000 ng/ml
MTD	Methadone	300 ng/ml
MOR	Morphine	300 ng/ml
MOR II	Morphine II	2,000 ng/ml
PCP	Phencyclidine	25 ng/ml
PPX	Propoxyphene	300 ng/ml
TCA	Tricyclic Antidepressants	1,000 ng/ml
THC	11-nor- Δ^9 -THC-9-COOH	50 ng/ml

These products are designed to obtain rapid, visual, qualitative results and are intended for professional use. The assay should not be used without proper supervision and is not intended for over-the-counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by Substance Abuse and Mental Health Services Administration (SAMHSA). Clinical consideration and professional judgment should be applied to any drugs of abuse test result particularly when preliminary positive results are indicated.

SUMMARY

Urine based screening tests for drugs of abuse range from simple immunoassays to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method for screening urine for drugs of abuse. The SureStep™ Drug Screen Tests are based on the principle of highly specific immunochemical reactions of antigens and antibodies which are used for the analysis of specific compounds in biological fluids.

Amphetamine: Amphetamine is a sympathomimetic amine whose biological effects include potent central nervous system stimulation, appetite reduction, hyperthermia, insomnia and euphoria. Pharmaceutical names of amphetamine include Adderall® and Dexedrine®. Common street names include "speed" and "uppers". The half-life of the drug is approximately 12 hours. The detection window for amphetamine in urine is 1 to 2 days after use.

Barbiturates: Barbiturates are a class of central nervous system depressants used medically as sedatives and anticonvulsants. Members of the barbiturate class include phenobarbital, pentobarbital, secobarbital, butalbital and amobarbital. Common physical effects of barbiturates include impaired motor coordination, anesthesia, sedation, cardiovascular and respiratory depression. The half lives of barbiturate range from 2 to 40 hours, depending on the duration of use. Short-acting barbiturates, such as secobarbital and butalbital, can be detected in urine up to 4 days after use. Long-acting barbiturates, such as phenobarbital, can be detected in urine up to 3 weeks.

Benzodiazepines: Benzodiazepines are a class of drugs frequently prescribed for the treatment of anxiety, sleep disorders and some seizure conditions. Members of this drug class include diazepam, chlordiazepoxide, alprazolam and clonazepam. Physical and psychological effects of benzodiazepines include lethargy, reduced motor coordination and drowsiness. Physical dependence can occur and withdrawal symptoms may appear when the medication is discontinued. The half-lives of benzodiazepines range from 2 to 40 hours, depending on the duration of use. Benzodiazepines may be detected in the urine up to 10 days after use.

Cocaine: Derived from the leaves of the coca plant, cocaine is a potent central nervous system stimulant and local anesthetic. Physical and psychological effects of cocaine use include increased heart rate, fever, pupil dilation, diaphoresis, euphoria and increased energy. Biologically, cocaine is rapidly metabolized to benzoylecgonine. The half-life of benzoylecgonine (5 to 8 hours) is much longer than that of the parent compound cocaine (0.5 to 1.5 hours). Benzoylecgonine can be detected in urine up to 3 days after cocaine use.

Ecstasy: 3,4-methylenedioxy-n-methamphetamine (MDMA) is a synthetic psychoactive drug chemically related to methamphetamine but with greater hallucinogenic properties. MDMA is taken orally, usually in tablets or capsules of 80 -150 mg. The effects of the drug last 4 to 6 hours following oral administration. Like methamphetamine, adverse effects of MDMA use include elevated blood pressure, hyperthermia, anxiety, paranoia and insomnia. MDMA overdoses can be fatal, often the result of heart failure or extreme heat stroke. Following a typical dose, 65% of MDMA is excreted in the urine as unchanged drug and up to 7% is demethylated and eliminated in the urine as methylenedioxyamphetamine (MDA). MDMA is detectable in the urine for up to 3 days after use.

Methamphetamine: Methamphetamine is a potent sympathomimetic amine with therapeutic applications. High doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, and a sense of increased energy. Other physical responses to methamphetamine use include cardiac dysrhythmias, paranoia and psychotic behaviors virtually indistinguishable from schizophrenia. Methamphetamine has a half-life of about 15 hours. In the urine, both methamphetamine and its metabolite amphetamine can be detected up to 3 days after use.

Methadone: Methadone is a synthetic analgesic drug primarily used in the treatment of narcotic addiction. It is administered orally or intravenously and is metabolized in the liver to EDDP, EMDP and methadol. Physiological effects of methadone use include respiratory depression, sedation and analgesia. Methadone overdoses are characterized by stupor, hypotension, pupillary constriction and clammy skin. Left untreated, methadone overdoses can lead to coma or death. Methadone has a half-life of 15 to 60 hours and is detectable in the urine (as parent methadone, EDDP, EMDP and methadol) for up to 3 days after use.

Morphine: Opiates, including heroin, morphine and codeine, are derived from the opium poppy. Medicinally, opiates are used for pain management and cough reduction. Heroin is an illicit analog of morphine that breaks down in the body to 6-acetylmorphine, morphine and morphine glucuronide. The biological half-lives for opiates range from 3 to 4 hours. In the urine, opiates are detectable up to 3 days after use. Morphine is the metabolite of both heroin and codeine; morphine (or morphine glucuronide) in the urine could indicate morphine, heroin and/or codeine use.

Phencyclidine: Phencyclidine is an arylcyclohexylamine that was originally developed for use as an anesthetic. Although used in veterinary medicine as a tranquilizer, the use of the drug in humans was discontinued due to the negative side effects. PCP is manufactured illegally and sold on the streets under the names "angel dust," "ozone," and "crystal cyclone." Common physical effects of the drug include hallucinations, disorientation, muscular incoordination and delusions. The half-life of PCP is about 12 hours. The detection window in urine is up to 5 days for infrequent use and up to 20 days for chronic use.

Propoxyphene: Propoxyphene or dextropropoxyphene is a narcotic analgesic compound with a structural similarity to methadone. It is prescribed in the United States for the relief of moderate pain. Darvocet™, one of the most common brand names for the drug, contains 50-100 mg of propoxyphene napsylate and 325-650 mg of acetaminophen. Physiological effects of propoxyphene include respiratory depression. Propoxyphene is metabolized in the liver to yield norpropoxyphene. Norpropoxyphene has a longer half-life (30 to 36 hours) than that of propoxyphene (6 to 12 hours). Norpropoxyphene demonstrates substantially less central nervous system depression than propoxyphene but shows a greater local anesthetic effect.

Tricyclic Antidepressants: Tricyclic Antidepressants are used in the treatment of depressive disorders and are often implicated in accidental and intentional overdoses. Members of this drug class include amitriptyline, nortriptyline, imipramine, desipramine and clomipramine. TCA overdoses can result in profound central nervous system depression, cardiotoxicity and anticholinergic effects. Half-lives for these drugs range from 2 to 70 hours. Detection times in urine vary depending on duration of drug use.

THC: Tetrahydrocannabinol (THC) is considered the most psychoactive of the more than 400 chemicals found in the marijuana plant (*cannabis sativa*). Physical and psychological effects of marijuana use include elevated mood, increased appetite, increased cardiac output, apathy and altered perception. The ability of the drug to enhance the appetite has made marijuana attractive for use with AIDS and cancer patients suffering from nausea and vomiting. THC is extensively metabolized to 11-nor- Δ^9 -THC-9-COOH which has a half-life of 24 hours. Urine detection limits vary widely depending on the frequency of drug use. In chronic users, 11-nor- Δ^9 -THC-9-COOH can remain detectable in the urine for up to 28 days.

PRINCIPLE

Each SureStep™ Drug Screen Test is a one-step immunoassay in which a chemically labeled drug (drug-protein conjugate) competes with the drug/metabolites that may be present in urine for limited antibody binding sites. The test device contains a membrane strip that has been pre-coated with drug-protein conjugate on the test line region. The colored antibody-colloidal gold conjugate moves along with urine, chromatographically by capillary action, across the membrane. In the absence of drug/metabolites in the urine, the colored antibody colloidal gold conjugate attaches to the drug-protein conjugate on the test region to form a visible colored line as the antibody complexes with the drug protein conjugate. Therefore, the formation of a visible colored line in the test region occurs when the urine is **negative** for the drug. When the drug/metabolite is present in the urine, the drug/metabolite antigen competes with the drug-protein conjugate on the test region for the limited antibody binding sites. When a sufficient amount of drug/metabolite is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug-protein conjugate on the test region. Therefore, absence of the colored line on the test region indicates a **positive** result.

A control line with a different antigen/antibody reaction is also added to the immunochromatographic membrane strip at the control region to indicate that the test is performed properly. A colored line should always appear in the control region regardless of the presence of the drug or metabolites. Therefore, **negative** urine will produce **two** colored lines, and **positive** urine will produce only **one** line. The presence of this colored line in the control region serves as verification that 1) a sufficient volume of urine has been added, and 2) that proper flow was obtained.

STORAGE AND STABILITY

The test device should be stored refrigerated or at room temperature 2-30°C (36-86°F). Each device should remain in its sealed pouch until use or stated expiration date.

PRECAUTIONS

- FOR IN VITRO DIAGNOSTIC USE.
- For professional use only.
- The pouch containing the device should be sealed. Discard the device if the package is ripped or torn.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Do not reuse.

MATERIALS SUPPLIED

- Twenty-five (25) individually wrapped test devices each containing one disposable pipette.
- One instructional booklet.

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Timer.
- External positive and negative controls.

SPECIMEN COLLECTION AND HANDLING


The SureStep™ Drug Screen Tests are designed for use with urine specimens. Fresh urine does not require any special handling or pretreatment. Urine samples should be collected such that testing may be performed as soon as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8°C for 2 days or frozen at -20°C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing. Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.

Note: Urine specimens, and all materials coming in contact with them, should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire.


TEST PROCEDURE

Review "Specimen Collection" instructions. The test device, patient samples and controls should be brought to room temperature (20-30°C) prior to testing. Do not open pouches until ready to perform the assay.

1. Remove the test device from its protective pouch. Label the device with patient or control identification.
2. Using the pipette with each device, draw the urine sample up to marked line (approximately 0.2 ml) dispense entire contents into the sample well. Avoid adding drops that contain air since air bubbles in the well may cause uneven flow or prevent the flow of the sample onto the test strip.
3. **Results may appear as early as 5 minutes. Verify positive results at a read time of 8 minutes.** Do not interpret result after 8 minutes.

INTERPRETATION OF RESULTS

Negative (-): Two colored lines appear adjacent to each other in the result window. This indicates that no drug above the cut-off level has been detected. **The test line intensity may be any shade of pink and may be weaker or stronger than the control line.**

Positive (+): A colored line appears in the control region (C) and no colored line appears in the drug test region. A sample is positive for the drug if no colored line appears in the drug test region. A positive result indicates that the drug is present in the specimen at or above the stated cut-off level. A positive result should not be considered conclusive and should be confirmed with a more specific alternative method such as GC/MS.

Invalid: A result is invalid if no colored line appears in the control region (C). Under no circumstances should a positive sample be identified until the control line forms in the viewing area. If the control line does not form, the test result is invalid, improper testing procedure, sample tampering, or test deterioration probably occurred. If an invalid test occurs, retest with a new test device. If problem persists, call Applied Biotech, Inc. at (888) 578-7956 for calls within the U.S. or (858) 713-9668 for international calls for technical assistance. Calls are answered by English speaking personnel. Other languages should email marketing@abiinverness.com.

LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical and/or procedural errors as well as other substances or factors not listed may interfere with the test and cause incorrect results. See the specificity section for lists of substances that will either produce positive results, or that do not interfere with the test performance.
- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected and the test should be repeated.

QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. It is recommended to use external quality control materials with each new lot, new shipment and every 30 days. When testing the external positive and negative controls, use the same test procedure as with a urine specimen. Quality control standards are available from commercial sources. Each laboratory should follow their state and local quality control requirements.

PERFORMANCE CHARACTERISTICS

A. Accuracy

For each individual drug screen assay, approximately sixty (60) urine samples containing specific drugs were obtained from clinical laboratories. The concentration of drug present in the samples was determined by either GC/MS, HPLC or commercial kits. Each sample was then tested with a SureStep™ Drug Screen Test. Additionally, at least one hundred (100) urine samples, collected from presumed non-user volunteers, were tested. The results are listed in Table 1 on page 47.

B. Precision

The precision of the SureStep™ Drug Screen Tests were determined by conducting testing at four clinical sites with spiked urine controls containing known amounts of each drug found on the panel assay. In clinical testing, the precision of the SureStep™ Drug Screen Test was determined to be an average of 99% correlation with specimens containing drugs targeted at 50% of the cut-off value for each drug and greater than 99% correlation with specimens containing drugs targeted at 200% of the cut-off value for each drug.

C. Specificity

Specificity for each of the SureStep™ Drug Screen Tests was evaluated by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal urine. The following structurally related compounds produced positive results only with the specified drug screen assay but not with the other assays when tested at levels equal to or greater than the concentrations listed below:

Compound	Concentration (ng/ml)	Compound	Concentration (ng/ml)
Amphetamine			
d-Amphetamine	1,000	Methamphetamine (con't.)	
l-Amphetamine	10,000	Chloroquine	50,000
MDA	5,000	(+/-)-Ephedrine	50,000
Phentermine	1,000	(-)-Methamphetamine	25,000
Tyramine	75,000	Mephentermine	50,000
Amphetamine II			
d-Amphetamine	300	MDMA	2,000
l-Amphetamine	2,500	β-Phenylethylamine	50,000
MDA	500	Ranitidine	50,000
β-Phenylethylamine	100,000	Trimethobenzamide	10,000
d-Methamphetamine	100,000	Methadone	
l-Methamphetamine	100,000	Methadone	300
MDMA	100,000	Doxylamine	50,000
Barbiturates			
Secobarbital	300	EDDP	50,000
Allobarbital	1,000	Methadol	25,000
Alphenal	300	Morphine	
Amobarbital	1,000	Morphine	300
Aprobarbital	300	Codeine	300
Barbital	300	Ethyl Morphine	300
Butobarbital	300	Hydrocodone	5,000
Butalbital	2,000	Hydromorphone	5,000
Butethal	300	Morphine-3-β-d-glucuronide	1,000
Pentobarbital	300	Thebaine	30,000
Phenobarbital	300	Morphine II	
Benzodiazepines			
Oxazepam	300	Morphine	2,000
Alprazolam	600	Codeine	2,000
Bromazepam	100	Ethyl Morphine	2,000
Chlordiazepoxide	300	Hydrocodone	25,000
Clobazam	300	Hydromorphone	25,000
Clonazepam	300	Morphine-3-β-d-glucuronide	5,000
Chlorazepate	200	Phencyclidine	
Delorazepam	3,000	Phencyclidine	25
Diazepam	300	EDDP	50,000
Estazolam	300	Methylphenidate	25,000
Flunitrazepam	300	Tenocyclidine	2,000
Flurazepam	150	Propoxyphene	
Lorazepam	500	Propoxyphene	300
Lormetazepam	500	Norpropoxyphene	1,000
Medazepam	2,000	Tricyclic Antidepressants	
Nitrazepam	250	Amitriptyline	1,000
Nordiazepam	150	Cyclobenzaprine	1,500
Prazepam	1,500	Clomipramine	5,000
Temazepam	150	Desipramine	600
Triazolam	200	Doxepin	3,000
Cocaine			
Benzoyllecgonine	300	Imipramine	600
Cocaine	300	Nortriptyline	1,000
Cocaine II			
Benzoyllecgonine	150	Nordoxepin	1,000
Cocaine	150	Perphenazine	75,000
Ecstasy			
MDMA	500	Promazine	15,000
MDA	2,000	Protriptyline	2,000
MDEA	1,000	Trimipramine	2,000
Methamphetamine			
(+)-Methamphetamine	1,000	THC	
d-Amphetamine	50,000	11-nor-Δ ⁹ -THC-9-COOH	50
		11-nor-Δ ⁸ -THC-9-COOH	50
		11-hydroxy-Δ ⁹ -Tetrahydrocannabinol	2,500
		Δ ⁹ -Tetrahydrocannabinol	7,500
		Δ ⁹ -Tetrahydrocannabinol	10,000
		Cannabinol	10,000
		Cannabidiol	100,000

MDA=(+/-)3,4-methylenedioxyamphetamine
MDMA=(+/-)3,4-methylenedioxy-n-methamphetamine

MDEA=Methylenedioxyethylamphetamine
EDDP=2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine

The following compounds listed were found not to cross react with any of the SureStep™ Drug Tests. Each compound was found not to cross react when tested at concentrations up to 100 µg/ml.

Acetaminophen	Diphenhydramine	β-Naphthalene Acetic Acid
Acetone	5,5-Diphenylhydantoin	(+/-)Naproxen
Albumin	Dopamine	(-)-Nicotine
Amoxapine	Ecgonine	Nicotinic Acid
Ampicillin	Ecgonine Methyl Ester	(+/-)-Norephedrine
l-Ascorbic Acid	(-)-ψ-Ephedrine	Normetanephrine
Aspartame	(+)-ψ-Ephedrine	Noscapine Hydrochloride
Aspirin	(+/-)-Epinephrine	Orphenadrine
Atropine	(+)-Epinephrine	Oxalic Acid
Baclofen	Erythromycin	Penicillin-G
Benzocaine	Ethanol	Phenothiazine
Bilirubin	Fenofibrate	Phenelzine
(+)-Brompheniramine	Fentanyl	l-Phenylephrine
Buprenorphine	Fluoxetine	(+/-)-Phenylpropranolamine
Caffeine	Furosemide	Primidone
Carbamate	Gemfibrozil	Procaine
Carbamazepine	Glucose	Promethazine
Carisoprodol	Guaiacol Glyceryl Ether	2-Propylpentanoic acid
(+)-Chlorpheniramine	Hemoglobin	Pseudoephedrine
(+/-)-Chlorpheniramine	DL-Homatropine	d-Propoxyphene
Chlorpromazine	Hydrochlorothiazide	Quinine
Chlorprothixene	3-Hydroxytyramide	Quinidine
Chlorthalidone	Ibuprofen	Riboflavin
Clofibrate	(+/-)- Isoproterenol	Salicylic acid
(-)-Cotinine	Lidocaine	Sodium Chloride
Creatine	Maprotiline	Sulindac
Creatinine	Meperidine	Theophylline
r-cyclodextrin	Methanol	Thioridazine
Cyproheptadine	Methapyrilene	cis-Thiothixene
Dantrolene	Methaqualone	d(+)-Trehalose
Dexamethasone	(1R,2S)-(-)-n-Methyl-Ephedrine	Trifluoperazine
Dexbrompheniramine	2-(N-morpholine)-Ethane	Vitamin C
Dextromethorphan	Sulphonic Acid	
4-Dimethylaminoantipyrine	Naltrexone	

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1. A. Shulgin. The background and chemistry of MDMA. *Journal of Psychoactive Drugs*, 18(4), 291-304, 1986.
2. Aniline O., Pittes, F. N., Phencyclidine (PCP): A Review and Perspectives. *CRC Critical Reviews in Toxicology*, 1982, 10, 145-177.
3. Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.
4. Ellenhorn, M.J. and Barceloux, D.G. *Medical Toxicology*. Elsevier Science Publishing Company, Inc., New York, 1988.
5. Federal Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69,11970-11979, 1988.
6. Gilman, A. G., and Goodman, L. S. The Pharmacological Basis of Therapeutics, eds., MacMillan Publishing, New York NY, 1980.
7. Gorodetzky, C. W., Detection of Drugs of Abuse in Biological Fluids, Martin WR(ed): Drug Addiction I, New York, Spring - Verlag, 1977.
8. Greenblatt, D.J., Shader, R.I. Benzodiazepines in Clinical Practice. New York: Raven Press, 1974.
9. H-J. Helmin et al. Analysis of 3,4-methylenedioxymethamphetamine (MDMA) and its metabolites in plasma and urine by HPLC-DAD and GC/MS. *Journal of Analytical Toxicology*, 20, 432-438, 1996.
10. H. Maurer et al. Toxicokinetics and Analytical Toxicology of Amphetamine-derived Designer Drugs ("Ecstasy").
11. Harvey, R.A., Champe, P.C. *Lippincotts Illustrated Reviews*. Pharmacology. 91-95, 1992.
12. Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.
13. Horsmans Y, Saliez A, van den Berge V, Desager J.P., Geubel A.P., Pauwels S, and Lambotte L. 14C-Propoxyphene Demethylation in the Rat. An Example of Differences Between Liver and Intestinal Drug-presystolic Metabolism. Drug Metabolism and Disposition: The Biological Fate of Chemicals, 1997, 25(11): 1257-9.
14. Jonasson U, Jonasson B, and Saldeen T. Middle-aged Men - A Risk Category Regarding Fatal Poisoning Due to Dextropropoxyphene and Alcohol in Combination. *Preventative Medicine*. 2000, 31 (2 Pt 1): 103-06.
15. McBay, A.J., *Clinical Chemistry*, 33, 33B-40B, 1987.
16. Miller, R.R. Propoxyphene: A Review. *American Journal of Hospital Pharmacology*, 1977, 34 (4): 413-23.
17. Physicians Desk Reference. 53rd. Medical Economics Company. 1999, pp 1567-8.
18. R. Torre et al. Non-linear Pharmacokinetics of MDMA ("ecstasy") in Humans. *Journal of Clinical Pharmacology*, 49, 104-109, 1999.
19. S-J. Peroutka ed. Ecstasy: The Clinical, Pharmacological and Neurotoxicological Effects of the Drug MDMA. Kluwer Academic Publishers, 1990.
20. Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA), Research Monograph, 73, 1986.
21. Young R. J. Dextropropoxyphene Over Dosage. Pharmacological Considerations and Clinical Management. *Drugs*, 1983, 26 (1):70-79.

REF

5108 KAB AMP
 5608 KAB BAR
 5708 KAB BZO
 5208 KAB COC
 5110 KAB MDMA
 5308 KAB M-AMP
 5908 KAB MTD

REF

5408 KAB MOR
 5428 KAB MORII
 5808 KAB PCP
 5974 KAB PPX
 5008 KAB TCA
 5508 KAB THC

**Export Only/Nur für den Export/Sólo para exportación/Réservé à l'exportation/
 Solo esportazione/Apenas para exportação/Endast export/Kun eksport**

REF

5128 KAB AMPII
 5228 DAB COCII

<p>CE CE Mark CE-Kennzeichen Marca CE Marque CE Marchio CE Marca de Conformidade CE CE-märke CE märke</p>	<p>i Consult Insert Packungsbeilage beachten Consulta el prospecto Consulter le texte inséré Consultare il foglietto illustrativo Consultar o Folheto Informativo Se Produktbladet Se indlægsseddel</p>	<p>ⓧ Do not Reuse Nicht wiederverwenden No reutilizar Ne pas réutiliser Non riutilizzare Não Reutilizar Återanvänd ej Má ikke genanvendes</p>
<p>EC REP Authorized Representative Autorisierter Vertreter Representante autorizado Représentant agréé Rivenditore autorizzato Representante Autorizado Autoriseerad Representant Autoriseert representant</p>	<p>IVD In vitro diagnostic device In-vitro Diagnostikum Dispositivo para diagnóstico Dispositif de diagnostic in vitro Dispositivo ad uso diagnostico in vitro Dispositivo para diagnóstico in vitro Enhet för in vitro-diagnostik In vitro diagnostisk udstyr</p>	<p>LOT Lot number/Batch code Chargenbez./Loscode Numero de lote/Código de remesa Número de lote/Code du lot Numero lotto/codice batch Número de Lote/Código do Lote Satsnummer/Batchkod Lot-nummer/Batch-kode</p>
<p>Manufactured by Hersteller Fabricado por Fabriqué par Produttore Fabricado Por Tillverkad av Fremstillet af</p>	<p>REF Catalog Number Best.-Nr. N. de referencia Référence Catalogue N. di catalogo Número de Catálogo Katalognummer Katalognummer</p>	<p>20°C Temperature Limitation Temperaturbeschränkung Limite de temperatura Limites de température Limiti di temperatura Limitação de temperatura Temperaturbegränsning Temperaturbegrensning</p>
<p>Use By/Expiration Date Verwendbar bis/Verfalldatum Utilizar antes de/Fecha de caducidad A utiliser jusque/Date d'expiration Utilizzare entro/Data di scadenza Utilizar Até/Prazo de Validade Använd före/Utgångsdatum Anvendes inden/Udløbsdato</p>	<p>Σ Quantity Quantität Cantidad Quantité Quantità Quantidade Antal</p>	<p>! Attention Achtung Atención Attention Attenção OBS! Bemærk</p>

Applied Biotech, Inc.
 10237 Flanders Ct.
 San Diego, CA 92121 USA
 (858) 587-6771 • Fax (858) 587-6757
 email: marketing@abinverness.com

Dr. Roland Seidel GmbH
 Paul-Ehrlich-Str. 32-H8
 D63322 Roedermark/Frankfurt
 Germany
 Tel: +49 (0) 6074/97725
 Fax: +49 (0) 6074/881446