## (KDTC) Training Quiz

* Name:	*Organization:	•	•	<u> </u>
*Address:				
*City:	*State:	*Zip:		
*E-mail:				
*Phone:  *Denotes required fields. Your test will not be evaluated until these fields.  Instructions: This is a multiple-choice test. Read the choosing the best answer. Click on the button that co You must complete the test with a score of 80 percent	questions completely before rresponds to the best answer.			
<ol> <li>The Know Drug Test Cup (KDTC)</li> <li>( ) a. Preliminary drug test results</li> <li>( ) b. quantitative drug test results</li> <li>( ) c. qualitative drug test results</li> <li>( ) d. confirmed test results</li> <li>( ) e. both a and c are correct</li> </ol>	provides			
2. The presence of a colored line in (T1 and T2) indicates that ( ) a. the test is invalid ( ) b. the test is negative ( ) c. the test is presumptive positive	the control region a	nd a colored line	in both test region	1S
3. The absence of a colored line in to ( ) a. the test is invalid ( ) b. the test result is negative ( ) c. the test result is presumptive po ( ) d. the test should be repeated usin ( ) e. both a and d are correct	ositive	;) indicate		
4. The Know Drug Test Cup (KDTC)	should be stored			
( ) a. in the freezer ( ) b. at 2-30 C (36-86 F) (either roon ( ) c. only in the refrigerator	n temperature or the	refrigerator)		
5. Test results can be read up to: ( ) a. 2 days ( ) b. 1 hour ( ) c. 15 minutes ( ) d. Between 5-60 minutes				
6. Provided the control line appears ( ) a. Yes ( ) b. No	s, it is ok to use test	cups past their st	ated expiration da	ite

## (KDTC) Training Quiz Continued

7. Test results can be read up to: The color intensity of a line at the "T" region must be as dark as the color intensity of the line at the "C" region.	
( ) a. True ( ) b. False	
<ul><li>8. You should wait 5 minutes before reading the test results</li><li>( ) Yes</li><li>( ) No</li></ul>	
<ul><li>9. A "ghost" line with valid control lines should be interpreted as</li><li>( ) a. Presumptive positive</li><li>( ) b. Negative</li></ul>	
<ul><li>10. External quality controls are recommended to ensure proper kit performance</li><li>( ) a. True</li><li>( ) b. False</li></ul>	

PRINT AND FAX THIS TEST TO 801.720.7568