

# HIV 1/2 STAT-PAK<sup>®</sup> Assay

A rapid, qualitative assay for the detection of antibodies to HIV-1 and HIV-2 in fingerstick whole blood, venous whole blood, serum, or plasma.

## Convenient, Reliable and Cost Effective

- Ideally suited for point-of-care testing
- Rapid reporting: Results in 15 minutes
- 99.7% sensitivity (HIV-1), 99.9% specificity<sup>1</sup>
- Built-in IgG procedural control
- Easy to perform
- No special preparation of reagents required
- Room temperature storage 8° to 30°C (46° to 86°F)
- Long shelf life
- No special equipment required
- Each control pack contains sufficient volume to run 50 tests (available separately)

**Minimal sample size: 5 µL**



Ideal for: Community Based Organizations, AIDS Service Organizations, Planned Parenthood Affiliates, Family Planning Councils, Physician Offices and Outreach Centers


3661 Horseblock Road, Medford, NY 11763 USA  
844.CHEMBIO (844.243.6246)  
CustomerService@chembio.com

 **CHEMBIO**  
DIAGNOSTIC SYSTEMS, INC.  
*Rapid Tests for Earlier Treatment™*


<sup>1</sup> See product insert for complete performance details

## Test Procedure\*


**1** TRANSFER loop sample to Sample Well (S)



**2** ADD three drops of Running Buffer



**3** READ results in 15 minutes



\* Before testing, please refer to the complete product insert provided in the test kit

## Ordering Information

Reimbursement Code: 86703-QW, G0433-QW

| Description   | Contents   | Catalog Number |
|---|--|----------------|
| HIV 1/2 STAT-PAK® Assay<br>FDA Approved, CLIA Waived/Non-Waived | 20 Individually Pouched Test Devices<br>20 Disposable 5µL Sample Loops<br>1 HIV Running Buffer (3.5 mL)<br>20 Copies of Subject Information Notice<br>1 Quick Reference Instructions (for CLIA Waived)<br>1 Product Insert | 60-9505-1      |
| HIV Rapid Test Control Pack<br>FDA Approved                     | One Pack Contains:<br>- 1 HIV-1 Reactive Control 0.25 mL<br>- 1 HIV-2 Reactive Control 0.25 mL<br>- 1 HIV Nonreactive Control 0.25 mL<br>Each kit control vial contains sufficient volume to run 50 tests.                 | 60-9549-0      |