

# CLIAwaived™ Inc. Rapid Drug Test Cup

## Step-by Step Instructions

The CLIAwaived Inc. Rapid Drug Test Cup is an *in vitro* screening test for the rapid detection of multiple drugs in human urine at or above the following cutoff concentration:

THC	11-nor- $\Delta^9$ -Tetrahydrocannabinol-9-carboxylic acid	50 ng/ml
COC	Benzoyllecgonine	150 ng/ml
OPI	Morphine	300 ng/ml
MET	d-Methamphetamine	500 ng/ml
AMP	d-Amphetamine	500 ng/ml
BZO	Oxazepam	300 ng/ml
BAR	Secobarbital	300 ng/ml
MTD	Methadone	300 ng/ml
BUPG	Buprenorphine Glucuronide	10 ng/ml
TCA	Nortriptyline	1000 ng/ml
MDMA	3,4-Methylenedioxymethamphetamine	500 ng/ml
OXY	Oxycodone	100 ng/ml
PCP	Phencyclidine	25 ng/ml
PPX	Propoxyphene	300 ng/ml

These tests provide visual qualitative results and are intended for *in vitro* diagnostic use only. The CLIAwaived Inc. Rapid Drug Test Cup is available in double drug analyte cassette dip format. It is intended for prescription point-of-care use and over-the-counter consumer use.

These tests provide only a preliminary test result and are the first step in a two-step process for detecting drugs of abuse in urine. The second step is confirming the results in a certified laboratory. For a quantitative result or to confirm preliminary positive results obtained by the CLIAwaived Inc. Rapid Drug Test Cup, a more specific alternative method such as Gas Chromatography/Mass Spectrometry (GC/MS) must be used. Clinical consideration and professional judgment must be applied to any drug of abuse test result, particularly when a preliminary positive result is indicated.

**This is a preliminary screening test that detects drug-of-abuse in urine at specified detection levels. To confirm preliminary positive results, a more specific method such as Gas Chromatography/Mass Spectrometry (GC/MS) must be used.**


### CONTENTS OF KIT

<p><b><u>For Testing:</u></b></p> <ul style="list-style-type: none"> <li>✓ 1 Step-by-Step Test Instructions</li> <li>✓ 25 Individually Wrapped Test Lids</li> <li>✓ 25 Specimen Cups</li> </ul>	 <p>Specimen Collection Cup</p>	 <p>Individually Wrapped Test Lid</p>
---	--	--

# CLIAwaived™ Inc. Rapid Drug Test Cup

## Step-by Step Instructions

### Items for Confirmation Testing (Optional items for professional point-of-care use):

<ul style="list-style-type: none"> <li>✓ 1 Specimen Bag</li> <li>✓ 2 Identification Labels</li> <li>✓ 1 Mailing Label</li> </ul>	 <p style="text-align: center;">Specimen Bag</p>	<table border="1" style="width: 100%; text-align: center;"> <tr> <td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td> </tr> <tr> <td>THC</td><td>OPI</td><td>AMP</td><td>BAR</td><td>BUPG</td><td>MDMA</td><td>PCP</td> </tr> <tr> <td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td> </tr> <tr> <td>COC</td><td>MET</td><td>BZO</td><td>MTD</td><td>TCA</td><td>OXY</td><td>PPX</td> </tr> </table> <p style="text-align: center;">ID#: DT14-XXXXX-XXXX</p> <p style="text-align: center;">Identification Label</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	THC	OPI	AMP	BAR	BUPG	MDMA	PCP	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	COC	MET	BZO	MTD	TCA	OXY	PPX	<div style="border: 1px solid black; padding: 5px; text-align: center;"> <p>CLIAwaived Inc. 11578 Sorrento Valley Rd. San Diego, CA 92121 USA</p> </div> <p style="text-align: center;">Mailing Label</p>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																									
THC	OPI	AMP	BAR	BUPG	MDMA	PCP																									
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																									
COC	MET	BZO	MTD	TCA	OXY	PPX																									


### **STORAGE**

Store the CLIAwaived Inc. Rapid Drug Test Cup at room temperature 59°F to 86°F (15°C to 30°C).


### **INSTRUCTION**

#### Step 1.

Collect fresh urine in the specimen cup. Make sure the urine is above **the minimum line.**



Open foil pouch. Remove test lid from pouch. Discard desiccant(s).



#### Step 2.

Twist the lid onto the cup. The cup lid must be closed tightly.

#### Step 3.

Tilt the cup on its side to activate the test. **Read test results at 5 minutes. Do not read after 8 minutes.**



5 minutes

# CLIAwaived™ Inc. Rapid Drug Test Cup

## Step-by Step Instructions

### INTERPRETATION OF RESULTS

Each strip contains two drug tests. C region shows validity of a test result. T1 region shows result for Test 1. T2 region shows result for Test 2.

#### For C region:

The appearance of a line indicates a valid result.

No line means an **Invalid** result. If a test strip does not have a line in the C region, test results are **Invalid** for both T1 and T2 on that strip.



#### For T1 and T2 regions:

The appearance of a line indicates a **Negative** result.


Note: **Any test line, even a very faint test line, is considered a negative result.**

No line indicates a **Preliminary Positive** result.


Note: **Any urine with preliminary positive results should be sent to a laboratory for confirmation.**



**Example #1:** There is a line appearing in both T1 and T2 regions on all test strips. Therefore, it is **Negative** for all tests.



**Example #2:** There is no line appearing in the T2 region on the third test strip. Therefore, it is **Preliminary Positive** for BZO test. All other tests are **Negative**.



**Example #3:** There is no line appearing in the C region on the seventh test strip. Therefore, it is **Invalid** for both PCP and PPX tests. All other tests are **Negative**.

# CLIAwaived™ Inc. Rapid Drug Test Cup

## Step-by Step Instructions

### **APPROXIMATE DRUG DETECTION TIME TABLE**

<u>Illicit Drug</u>	<u>Identifier</u>	<u>Cut-off Level<sup>1</sup></u>	<u>Minimum Detection Time<sup>2</sup></u>	<u>Maximum Detection Time<sup>2</sup></u>
Marijuana	THC	50 ng/ml	1-3 hours	1-7 days
Cocaine	COC	150 ng/ml	2-6 hours	2-3 days
Opiates	OPI	300 ng/ml	2-6 hours	1-3 days
Methamphetamine	MET	500 ng/ml	4-6 hours	2-3 days
Amphetamine	AMP	500 ng/ml	4-6 hours	2-3 days
Ecstasy	MDMA	500 ng/ml	2-7 hours	2-4 days
Phencyclidine	PCP	25 ng/ml	4-6 hours	7-14 days
Propoxyphene	PPX	300 ng/ml	2-8 hours	2-7 days

<u>Prescription Drug</u>	<u>Identifier</u>	<u>Cut-off Level<sup>1</sup></u>	<u>Minimum Detection Time<sup>2</sup></u>	<u>Maximum Detection Time<sup>2</sup></u>
Benzodiazepines	BZO	300 ng/ml	2-7 hours	1-4 days
Barbiturates	BAR	300 ng/ml	2-4 hours	1-3 weeks
Methadone	MTD	300 ng/ml	3-8 hours	1-3 days
Buprenorphine	BUPG	10 ng/ml	2-7 hours	1-6 days
Tricyclic Antidepressants	TCA	1000 ng/ml	8-12 hours	2-7 days
Oxycodone	OXY	100 ng/ml	1-3 hours	1-2 days

<sup>1</sup> Cut-off level is the lowest drug concentration in the urine that can be detected by the CLIAwaived Inc. Rapid Drug Test Cup.

<sup>2</sup> Drug clearance rates are dependent on many factors such as frequency of drug use, the amount of drug taken, metabolism rates, and even body fat content.

### **RECOMMENDED PROCEDURE FOR CONFIRMATION**

The following is a recommended procedure for confirming preliminary positive results:

#### **Step 1:**

Make sure the lid is twisted tightly. Place a check mark in the box on both Identification Labels for the drug(s) with preliminary positive result. Place one Identification Label onto the CLIAwaived Inc. Rapid Drug Test Cup. Place the other Identification Label below for your record:

*Place Identification Label Here*

#### **Step 2:**

Place the labeled CLIAwaived Inc. Rapid Drug Test Cup into the Specimen Bag and seal the bag.

#### **Step 3:**

Place specimen into a shipping box and affix the Mailing Label. Drop the shipping box with POSTAGE into any mailbox. Specimen should be mailed within 24 hours of collection.

**NOTE: Please call 1-888-882-7739 M-F 8am-4pm Pacific Time to request additional Identification Labels and Mailing Labels if needed.**

# CLIAwaived™ Inc. Rapid Drug Test Cup

## Step-by Step Instructions

### **OBTAIN CONFIRMATION RESULTS**

- Results will be ready 5 to 7 days after the sample is received in our laboratory.
- Dial 1-888-882-7739 M-F 8am-4pm Pacific Time to obtain confirmation result.
- You will need the identification number to access the confirmation result.

### **WARNINGS AND PRECAUTIONS**

- ❖ For *in vitro* diagnostic use only (not for internal use).
- ❖ The test is for one time use only. It is not reusable.
- ❖ Do not use the CLIAwaived Inc. Rapid Drug Test Cup after the expiration date printed on the pouch.
- ❖ Keep the CLIAwaived Inc. Rapid Drug Test Cup in its original sealed pouch until ready for use. Do not use the test if the pouch is ripped or torn.
- ❖ Certain foods or medications may cause the test to give false results.
- ❖ Contaminated or tainted urine sample may give false results.
- ❖ Send specimen with preliminary positive or uncertain results to a laboratory for confirmation.
- ❖ Urine may contain infectious diseases. Always wear gloves and wash hands with soap after handling.
- ❖ Do not use this test if you are color-blind.

### **UNDERSTAND THE TEST RESULTS & FOLLOW UP**

#### **What is a Preliminary Positive Result?**

The CLIAwaived Inc. Rapid Drug Test Cup is a screening test to detect the presence of drugs in human urine. This means that if a drug is present, you will usually get a preliminary positive test result. It is important to send out the specimen to confirm the preliminary positive result. This is because certain foods, supplements, beverages, or medicines can affect the results of the CLIAwaived Inc. Rapid Drug Test Cup.

#### **If the test results are preliminary positive, does it mean that you found drugs of abuse?**

No. Take no serious actions until you get the laboratory's result. Many factors may cause a false positive result in the home test. A positive test for a prescription drug does not mean that a person is abusing the drug. This is because the test does not indicate acceptable levels compared to abusive levels of prescribed drugs.

This test cannot be used for legal purposes. If you get a preliminary positive test result, you should send the urine to a laboratory to confirm the test result. The confirmation test is called gas chromatography/mass spectrometry or GC/MS.

#### **Many things can affect the accuracy of the tests, including (but not limited to):**

- The way you did the test.
- The way you stored the test or urine.
- What the person ate or drank before taking the test.
- Any other prescription or over-the-counter drugs the person may have taken before the test.

*Note: Some over-the-counter medications will produce the same test results as illegally-abused amphetamines. Please call 1-888-882-7739 for a complete list of substances or over-the-counter medications which may cross-react with this test.*

#### **Where can I get some help?**

Consult with a counselor, doctor, or a qualified professional to help you address drug abuse problems. The following organizations provide helpful resources on drug abuse prevention and recovery programs. These resources are for information purposes only.

Alcohol and Drug Information <http://store.samhsa.gov/home>

Community Anti-Drug Coalitions of America <http://www.cadca.org>

National Institute on Drug Abuse [www.drugabuse.gov](http://www.drugabuse.gov)

#### **QUESTIONS?**

Call 1-888-882-7739 M-F 8am-4pm Pacific Time

# For Professional Point-of-Care Use

## QUALITY CONTROL

**Internal control:** The *CLIAwaived Inc.* Rapid Drug Test Cup test device has built-in internal procedural controls. The appearance of the control band © is considered an internal procedural control. This band should always appear if adequate sample volume is used and the testing procedure is followed. Additionally, the background color should become clear and provide distinct test result. If the control band © does not appear then the test is invalid. The test should be repeated using a new device.

**External control:** It is recommended that negative and positive urine controls be used to initially test each new lot of product to ensure proper kit performance. The same assay procedure should be followed with external control materials as with a urine specimen. If external controls do not produce the expected results, do not run test specimens. Follow the proper federal, state and local guidelines when running external controls.

**Quality control testing at regular intervals is a good laboratory practice and may be required by federal, state or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure that the quality program employed meets the established standards.**

## PERFORMANCE CHARACTERISTICS

### PRECISION

A study was conducted at two laboratory and one physician offices in an effort to determine the precision of *CLIAwaived Inc.* Rapid Drug Test Cup over 12 or more consecutive days. Testing was conducted on the Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Marijuana, Methamphetamine, Methylendioxyamphetamine, Methadone, Opiates, Oxycodone, Phencyclidine, Propoxyphene, and Tricyclic Antidepressants assays by operators using three different lots of product to demonstrate the within-run, between-run and between-operator precision. An identical panel of coded samples, containing drugs at the concentration of ± 50% cut-off level was labeled as a blind and tested at each site. The correlation with expected results was >99% across all lots and sites (with a 95% confidence interval).

### ACCURACY

The accuracy of the *CLIAwaived Inc.* Rapid Drug Test Cup was evaluated in comparison to the results from GC/MS or LC/MS analysis. Thirty-six (36) negative drug-free urine samples were collected from volunteer donors and tested with both the *CLIAwaived Inc.* Rapid Drug Test Cup and the GC/MS or LC/MS method. Of the 36 negative urine samples tested, all were found negative by both methods. Additionally, for each drug test, a minimum of 40 clinical urine samples previously analyzed by GC/MS or LC/MS method with known concentration(s) of drug(s) values were blind labeled and evaluated. The results are summarized below:

Drug Test	GC/MS Neg.	GC/MS < -50%	GC/MS -50% to Cutoff	GC/MS Cutoff to +50%	GC/MS > +50%	% Agreement w/ GC/MS		
						Neg (-)	Pos (+)	
THC 50	Pos. (+)	0	0	1	6	35	97.7%	100%
	Neg. (-)	36	2	4	0	0		
COC 150	Pos. (+)	0	0	3	3	37	92.7%	97.6%
	Neg. (-)	36	0	2	1	0		
OPI 300	Pos. (+)	0	0	3	7	34	92.5%	100%
	Neg. (-)	36	0	1	0	0		
MET 500	Pos. (+)	0	0	0	5	67	100%	96.0%
	Neg. (-)	36	2	4	3	0		
AMP 500	Pos. (+)	0	0	2	5	36	95.1%	100%
	Neg. (-)	36	1	2	0	0		
BZO 300	Pos. (+)	0	0	3	4	39	92.5%	100%
	Neg. (-)	36	0	1	0	0		

Drug Test	GC/MS Neg.	GC/MS < -50%	GC/MS -50% to Cutoff	GC/MS Cutoff to +50%	GC/MS > +50%	% Agreement w/ GC/MS		
						Neg (-)	Pos (+)	
BAR 300	Pos. (+)	0	0	1	6	33	97.5%	95.1%
	Neg. (-)	36	0	3	2	0		
MTD 300	Pos. (+)	0	0	0	3	36	100%	97.5%
	Neg. (-)	36	0	4	1	0		
BUPG 10	Pos. (+)	0	0	1	4	38	97.5%	97.7%
	Neg. (-)	36	0	3	1	0		
TCA 1000	Pos. (+)	0	0	0	27	11	100%	92.7%
	Neg. (-)	36	0	4	3	0		
MDMA 500	Pos. (+)	0	0	1	3	40	97.5%	97.7%
	Neg. (-)	36	0	3	1	0		
OXY 100	Pos. (+)	0	0	2	6	38	95.2%	100%
	Neg. (-)	36	0	4	0	0		
PCP 25	Pos. (+)	0	0	0	3	36	100%	95.1%
	Neg. (-)	36	0	4	2	0		
PPX 300	Pos. (+)	0	0	2	4	36	95.0%	100%
	Neg. (-)	36	0	2	0	0		

## SPECIFICITY

The specificity for the *CLIAwaived Inc.* Rapid Drug Test Cup was determined by testing various drugs, drug metabolites, structurally related compounds, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine. The effect of specimens with various pH (4.5–9) and specific gravity (1.005–1.030) ranges was also evaluated and found not to interfere with *CLIAwaived Inc.* Rapid Drug Test Cup.

The following compounds produced positive results when tested at or above the concentrations listed below.

### AMP 500 ng/ml

Compound	ng/ml	Compound	ng/ml
d-Amphetamine	500	Phentermine	1,000
l-Amphetamine	20,000	β-Phenylethylamine	80,000
d,l-3,4-MDA	1,500		

### BAR 300 ng/ml

Compound	ng/ml	Compound	ng/ml
Allobarbitol	1,500	Butalbital	300
Alphenal	400	Butethal	400
Amobarbital	1,500	Pentobarbital	400
Aprobarbital	400	Phenobarbital	400
Barbital	400	Secobarbital	300
Butobarbital	400		

### BZO 300 ng/ml

Compound	ng/ml	Compound	ng/ml
A-Hydroxy Alprazolam	50	Lorazepam	1,500
Alprazolam	150	Lormetazepam	1,000
Bromazepam	800	Medazepam	2,000
Chlordiazepoxide	2,000	Nitrazepam	1,000
Clobazam	200	Nordiazepam	100
Clonazepam	4,000	Oxazepam	300
Delorazepam	6,000	Phenazepam	1,000
Diazepam	150	Prazepam	1,000
Estazolam	300	Temazepam	150
Flunitrazepam	1,000	Triazolam	1,500
Flurazepam	300		

### BUPG 10ng/ml

Compound	ng/ml	Compound	ng/ml
Buprenorphine	100	Norbuprenorphine	100
Buprenorphine Glucuronide	10	Norbuprenorphine Glucuronide	100

### COC 150 ng/ml

Compound	ng/ml	Compound	ng/ml
Benzoylcegonine	150	Ecgonine	65,000

### MDMA 500 ng/ml

Compound	ng/ml	Compound	ng/ml
d,l-3,4-MDA	2,000	d,l-3,4-MDMA	500
d,l-3,4-MDEA	250	d-Methamphetamine	50,000

**MET 500 ng/ml**

Compound	ng/ml
Ephedrine	10,000
p-Hydroxymethamphetamine	1,750
d,l-3,4-MDMA	1,000
d,l-3,4-MDEA	20,000

**MTD 300 ng/ml**

Compound	ng/ml
Doxylamine	50,000
2-Ethylidene-1,5-Dimethyl-1-3,3-Diphenylpyrrolidine	50,000

**OPI 300 ng/ml**

Compound	ng/ml
6-Acetylmorphine	500
6-Acetylcodeine	600
Codeine	300
Dihydrocodeine	500
Ethyl morphine	300
Heroin	100

**OXY 100 ng/ml**

Compound	ng/ml
6-Acetylcodeine	15,000
Codeine	5,000
Dihydrocodeine	2,000
Hydrocodone	300

**PCP 25 ng/ml**

Compound	ng/ml
4-Hydroxy Phencyclidine	500
Metaphit	500

**PPX 300ng/ml**

Compound	ng/ml
Propoxyphene	300

**TCA 1000 ng/ml**

Compound	ng/ml
Amitriptyline	1,000
Clomipramine	7,500
Cyclobenzaprine	1,500
Desipramine	750
Doxepin	1,000
Imipramine	750

**THC 50 ng/ml**

Compound	ng/ml
Cannabidiol	100,000
Cannabinol	50,000
11-nor- $\Delta$ 8-THC-9-COOH	50
11-nor- $\Delta$ 9-THC-9-COOH	50

**CONSUMER STUDY**

A consumer study was conducted to determine the performance of the device when used by untrained, laypersons following only the instructions in the product labeling. A total of 153 participants read a total of 5460 assays during the study and 5228 of those 5460 assays (95.8%) was interpreted correctly. Each assay was tested by these participants using spiked solutions targeted to 0%, 25%, 50%, 75%, 125%, 150%, and 175% of the assay cutoff level.

**INTERFERENCE**

The following compounds were found not to cross-react when tested at concentrations up to 100  $\mu$ g/ml (100,000 ng/ml).

Acetaminophen	Amitriptyline (except TCA assay)
Acetone	Amobarbital (except BAR assay)
Acetylsalicylic acid (Aspirin)	Amoxapine
6-Acetylcodeine (except OPI & OXY assay)	Amoxicillin
6-Acetylmorphine (except OPI assay)	Aprobarbital (except BAR assay)
Albumin	d-Amphetamine (except AMP assay)
Allobarbitol (except BAR assay)	l-Amphetamine (except AMP assay)
Alphenal (except BAR assay)	Ampicillin
Alprazolam (except BZO assay)	Apomorphine
	l-Ascorbic Acid (Vitamin C)

Aspartame	$\alpha$ -Hydroxy Alprazolam (except BZO assay)
Atropine	4-Hydroxy Phencyclidine (except PCP assay)
Barbital (except BAR assay)	p-Hydroxymethamphetamine (except MET assay)
Benzilic acid	11-Hydroxy- $\Delta$ 9-THC (except THC assay)
Benzocaine (Ethyl p-Aminobenzoate)	Ibuprofen
Benzoic acid	Imipramine (except TCA assay)
Benzoyllecgonine (except COC assay)	d,l-Isoproterenol
Benzphetamine	Ketamine
Bilirubin	Lidocaine
Bromazepam (except BZO assay)	Lorazepam (except BZO assay)
d-Brompheniramine	Lormetazepam (except BZO assay)
Buprenorphine (except BUPG assay)	Medazepam (except BZO assay)
Butabarbital (except BAR assay)	Meperidine
Butalbital (except BAR assay)	Metaphit (except PCP assay)
Butethal (except BAR assay)	Methadone (except MTD assay)
Caffeine	d-Methamphetamine (except MET & MDMA assay)
Cannabidiol (except THC assay)	l-Methamphetamine (except MET assay)
Cannabinol (except THC assay)	Methaqualone
Chlordiazepoxide (except BZO assay)	Methoxyphenamine (1R,2S) N-Methyl-Ephedrine
Chloroquine	2-Methylamine-Propiophenone
d,l-Chlorpheniramine	d,l-3,4-Methylenedioxyamphetamine (except AMP & MDMA assays)
Chlorpromazine	d,l-3,4-methylenedioxyethylamphet (except MET & MDMA assays)
Cholesterol	d,l-3,4-Methylenedioxyamphetamine (except MET & MDMA assays)
Clobazam (except BZO assay)	Methylphenidate
Clomipramine (except TCA assay)	Morphine (except OPI assay)
Clonazepam (except BZO assay)	Morphine-3- $\beta$ -D-Glucuronide (except OPI assay)
Cocaine	Nalidixic acid
Codeine (except OPI & OXY assays)	Nalorphine (except for OPI assay)
Cortisone	Naloxone
l-Cotinine	d-Naproxen
Creatine	Niacinamide
Creatinine	Nitrazepam (except BZO assay)
Cyclobenzaprine (except TCA assay)	Nordiazepam (except BZO assay)
Delorazepam (except BZO assay)	Nordoxepin (except TCA assay)
Deoxycorticosterone	Nicotine, (S)-
Desipramine (except TCA assay)	Norepinephrine
Dextromethorphan	Norethindrone
Diazepam (except BZO assay)	Norpropoxyphene (except PPX assay)
Dihydrocodeine (except OPI & OXY assay)	Nortriptyline (except TCA assay)
4-Dimethylaminoantipyrine	Oxalic Acid
Diphenhydramine	Oxazepam (except BZO assay)
Dopamine (3-Hydroxytyramine)	Oxolinic acid
Doxepin (except TCA assay)	Oxycodone (except OXY assay)
Doxylamine (except MTD assay)	Oxymorphone (except OXY assay)
Ecgonine (except COC assay)	Papaverine
Ecgonine Methyl Ester	Penicillin-G (Benzylpenicillin)
l-Epinephrine	Pentazocine
d,l-Ephedrine (except MET assay)	Pentobarbital (except BAR assay)
Erythromycin	Perphenazine (except TCA assay)
Estazolam (except BZO assay)	Phenazepam (except BZO assay)
$\beta$ -Estradiol	Phencyclidine (except PCP assay)
Estrone-3-Sulfate	Phencyclidine Morpholine (except PCP assay)
Ethanol	Pheniramine (except MTD assay)
Ethyl Morphine (except OPI & OXY assay)	Phenobarbital (except BAR assay)
Ethyl-p-aminobenzoate	Phenothiazine (Thiodiphenylamine)
2-Ethylidene-1,5-Dimethyl-1-3,3-Diphenylpyrrolidone (except MTD assay)	Phentermine (except AMP assay)
Flunitrazepam (except BZO assay)	Phenylephrine
Flurazepam (except BZO assay)	$\beta$ -Phenylethylamine (except AMP assay)
Furosemide	Prednisolone
Glucose	Prazepam (except BZO assay)
Gentisic acid	Procaine (except MET assay)
Glutethimide	Promazine (except TCA assay)
Guaiaicol Glyceryl Ether	Promethazine
Hemoglobin	Propoxyphene (except PPX assay)
Heroin (except OPI assay)	Protriptyline (except TCA assay)
Hippuric acid	11-nor- $\Delta$ 9-THC-9-Carboxylic Acid
Hydrochlorothiazide	
Hydrocodone (except OPI & OXY assays)	
Hydrocortisone	
Hydromorphone (except OPI & OXY assays)	

d-Pseudoephedrine	( <i>except THC assay</i> )
Pyrrrolidine	Thiamine
Quinidine	Thioridazine
Quinine	Triazolam ( <i>except BZO assay</i> )
Ranitidine	Trifluoperazine
Riboflavin	Trimethobenzamide ( <i>except MET assay</i> )
Salicylic acid	Trimipramine ( <i>except TCA assay</i> )
Secobarbital ( <i>except BAR assay</i> )	Tryptamine
Serotonin	d,l-Tryptophan
Sertraline	Tyramine
Sodium Chloride	d,l-Tyrosine
Sulfamethazine	Uric Acid
Sulindac	Verapamil
Temazepam ( <i>except BZO assay</i> )	Zomepirac
Tetracycline	
$\Delta$ 8-THC ( <i>except THC assay</i> )	
$\Delta$ 9-THC ( <i>except THC assay</i> )	
11-nor- $\Delta$ 8-THC-9-Carboxylic Acid ( <i>except THC assay</i> )	

#### BIBLIOGRAPHY OF SUGGESTED READING

1. Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man, Biomedical Publications, Davis, CA, 1982.
2. Urine testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.
3. Wong, R., The Current Status of Drug Testing in the US Workforce, Am. Clin. Lab., 2002; 21(1): 21-23
4. Wong, R., The Effect of Adulterants on Urine Screen for Drugs of Abuse: Detection by an On-site Dipstick Device, Am. Clin. Lab., 2002; 21(3); 14-18
5. Young, D.S. et. al., Clinical Chemistry, 21 (9), 1975.
6. U.S. Dept. of Transportation, Procedures for Transportation Workplace Drug and Alcohol Testing Programs. Federal Register, 1999 Dec.; 64(236); 69076
7. U.S. Dept. of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs. Federal Register, 2001 Aug.; 66(162): 43876
8. Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970–11979, 1988.
9. Liu, Ray H. and Goldberger, Bruce A., Handbook of Workplace Drug Testing, AACC Press (1995).
10. Gilman, A. G. and Goodman, L. S., The Pharmacological Basis of Therapeutics, eds. MacMillan Publishing, New York, NY, 1980.
11. McBay, A.J. Clin. Chem. 33, 33B-40B, 1987.
12. Ringsrud, K.M and Linne, J.J., Urinalysis and Body Fluids, A color Text and Atlas, Mosby-Year Book, Inc., 1995.
13. Baselt RC. Disposition of toxic Drugs and chemicals in Man. 6<sup>th</sup> Ed. Biomedical Publ., Davis, CA. 2002; 129

**CLIAwaived™ Inc.**  
Rapid Drug Test Cup



Emergo Europe  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands

Tel: 31-0-70-345-8570  
Fax: 31-0-70-346-7299



Manufactured for CLIAwaived, Inc.  
11578 Sorrento Valley Road  
San Diego, CA 92121 USA

Tel: 1-888-882-7739 (toll free)

Tel: 1-858-481-5031

Fax: 1-801-720-7568

Email: [info@cliawaived.com](mailto:info@cliawaived.com)

Web: [www.cliawaived.com](http://www.cliawaived.com)

P/N: PI-X07-CLIA-14-RDTC-A Rev. A, 09/2014

Ref: PI-DTA Rev. A

© 2014 CLIAwaived, Inc. Made in the U.S.A.