

CLIAwaived™ Inc. Rapid Test Cup Drug Training

CLIA-14-RDTC

**Product
Performance / Interpretation**



The *CLIAwaived Inc.* Rapid Drug Test Cup (RDTC)

- ***CLIAwaived Inc.* RDTC Cup**
 - Rapid on-site drug screening test device for the detection of drugs in urine samples.
 - Collect anywhere, anytime.
 - Results within minutes.

- **The *CLIAwaived Inc.* RDTC Cup tests for the following drugs:**

• AMP: Amphetamines	• mAMP: Methamphetamine
• COC: Cocaine	• OPI: Opiates
• PCP: Phencyclidine	• THC: Marijuana
• MTD: Methadone	• BZO: Benzodiazepines
• BAR: Barbiturates	• TCA: Tricyclic Antidepressants
• OXY: Oxycodone	• MDMA: Methylenedioxymethamphetamine (Ecstasy)
• PPX: Propoxyphene	• BUP G: Buprenorphine

- **Gather all necessary testing supplies**
 - *CLIAwaived Inc.* RDTC Test Cup
 - *CLIAwaived Inc.* RDTC specimen lid

The Process: Performing the Drug Screen


- **Have Donor Select Test Kit**

- Allow donor to select sealed *CLIAwaived Inc.* RDTK from test kit box.
- Donor should hand sealed test kit to collector to record lot number and test expiration date.
- Ensure “Exp.” (Expiration date) is within range.

- **Return Collection Device to Donor**

- Instruct donor on proper specimen collection.
- Have donor open foil pouch and remove test lid.
- Have donor provide urine sample in specimen cup.

CUP COMPONENTS



The Process: Performing the Drug Screen (cont' d)

Temperature Verification

- A temperature strip is present on the back of the collection cup serves as an initial specimen validity check.
- The temperature should be checked within **4 minutes** of the donor providing the specimen.
- A freshly voided specimen will be in the range of **90° -100° F.**




Read Green Color

The Process: Performing the Drug Screen (cont' d)

- Have collector screw the *CLIAwaived Inc.* RDTC lid onto the collection cup by twisting a quarter turn **after** lid is snug.
- **Tilt** cup on it's side to **activate** test.
- Read results of drugs of abuse tests in 5 minutes.
- **DO NOT** interpret the test results after 8 minutes for CLIA-14-RDTC and 10 minutes for CLIA-12-RDTC.

Secure test device cap to the specimen cup. NOTE: Screw the cap tightly to ensure adequate seal.




IMPORTANT:
Cup lid must be secured tightly by twisting lid a quarter turn **AFTER** lid is snug

CLIAwaived Inc. RDTC Test Results: Interpretation

- Interpret the drug test results at five **(5) minutes**.
- Drug test results are stable for up to **8 minutes** for CLIA-14-RDTC and **10 minutes** for CLIA-12-RDTC.
- Each drug test strip in the device includes an internal procedural control (C) ensuring proper device function.
- **Control lines** should form next to the “C” or control area on both strips.
- Drug test interpretation:
 - The **presence of a line** (or any indication of a line) indicates a normal or **Negative** result.
 - The **absence of a line** indicates an abnormal or **Presumptive Positive** result.

Read results at 5 minutes.
Do not interpret result after
10 minutes.



CLIAwaived Inc. RDTTC Test Results: Invalid

- Each test strip within the device includes an internal procedural control (C) that ensures proper device function.
- Control lines should form on all drug test strips, indicating the proper functioning of the test device.
- The **absence** of a control line in one or multiple test strips indicates that the test is **Invalid** (even if one or more drug lines are absent).
- **DO NOT** record test results. This test must be administered again, with a new test cup.

Note: It is recommended that drug screen results not be given to the donor



No line next to C =
Invalid Result

CLIAwaived Inc. RDTC Test Results: Negative

- The appearance of a line next to each and every **T1** and **T2 (Test Lines)** corresponding to a specific drug and **ALL** control areas.
- Record **Negative** test results and discard test device.
- **Any indication of a colored line** regardless of color intensity would be correctly interpreted as a **Negative** test result.



All Lines Present =
Negative Test Result

Note: It is recommended that drug screen results not be given to the donor

CLIAwaived Inc. RDTC Test Results: Negative

- The appearance of a line next to each and every **T1 and T2 (Test Lines)** corresponding to a specific drug and **ALL** control areas.
- This shows an example of a **light line** for MET. This should still be interpreted as a **Negative**.
- Record **Negative** test results and discard test device.
- **Any indication of a colored line** regardless of color intensity would be correctly interpreted as a **Negative** test result.

Note: It is recommended that drug screen results not be given to the donor



All Lines Present =
Negative Test Result

CLIAwaived Inc. RDTC Test Results: Negative

- The appearance of a line next to each and every **T1 and T2 (Test Lines)** corresponding to a specific drug and **ALL** control areas.
- This shows an example of a **pencil thin line**. This should still be interpreted as a **Negative**.
- Record **Negative** test results and discard test device.
- **Any indication of a colored line** regardless of color intensity would be correctly interpreted as a **Negative** test result.



All Lines Present =
Negative Test Result

Note: It is recommended that drug screen results not be given to the donor

CLIAwaived Inc. RDTC Test Results: Negative

- The appearance of a line next to each and every **T1 and T2 (Test Lines)** corresponding to a specific drug and **ALL** control areas.
- This shows an example of a **broken line**. This should still be interpreted as a **Negative**.
- Record **Negative** test results and discard test device.
- **Any indication of a colored line** regardless of color intensity would be correctly interpreted as a **Negative** test result.

Note: It is recommended that drug screen results not be given to the donor



All Lines Present =
Negative Test Result

CLIAwaived Inc. RDTC Test Results: Presumptive Positive

- The **absence** of a line next to one or more **T1 and T2 (Test Lines)** corresponding to a specific drug *and* the presence of a line in **ALL** control areas.
- Record **Presumptive Positive** test results.
Note: We suggest all Presumptive Positive test results be confirmed by an alternative method.
- Example shows **Presumptive Positive** test result.

Note: It is recommended that drug screen results not be given to the donor



Presumptive Positive

No line next to drug name =
Presumptive Positive test result

CLIAwaived Inc. RDTC Test Results: Presumptive Positive

- This shows an example of a **“ghost” line** for COC with valid control (C) lines. This should still be interpreted as a **Presumptive Positive**.
- Although there appears to be a change in the test area, the **absence of color** makes it different from a Negative.
- Record **Presumptive Positive** test results.
Note: We suggest all Presumptive Positive test results be confirmed by an alternative method.

Note: It is recommended that drug screen results not be given to the donor



Presumptive Positive

No line next to drug name =
Presumptive Positive test result

CLIAwaived Inc. RDTC Test Results: Presumptive Positive

- This shows an example of **two Presumptive Positive** drugs on one specimen.
- We suggest that **both results** need to be recorded and confirmed by an alternative method.
- Record **Presumptive Positive** test results.
Note: We suggest all Presumptive Positive test results be confirmed by an alternative method.

Note: It is recommended that drug screen results not be given to the donor



Presumptive Positive

No line next to drug name =
Presumptive Positive test result

Thank You

For more information, please visit us on the web:

www.cliawaived.com

Or call (858) 481-5031 *(phone)*

(888) 882-7739 *(toll-free)*

info@cliawaived.com *(e-mail)*