

Surgilance™ Safety Lancets

Distributed by

CLIAwaived™ Inc.



Surgilance™ Safety Lancets

The Surgilance® Safety Lancet is used for convenient micro blood sampling, offering simplicity and single-use safety for healthcare professionals as well as patients. The lancet is one of the easiest to use and safest lancets available in the market for micro blood sampling. Its innovative design protects users from needlestick injury and the unique one-step activation mechanism makes it extremely simple and reliable to use. The device is available with either a needle or blade, each with different depths of penetration, to provide optimal blood flow for different applications and skin types.

BENEFITS:

- Does not require arming
- Color-coded for easy recognition
- Needle is safely concealed before and after use
- Patented design makes blood sampling easy and safe
- Once the lancet is used, it is rendered inoperative, providing added safety for patient and clinician
- High-speed delivery and penetration method minimizes patient pain and clinician error
- Latex-free and hypoallergenic

Ordering Information:

Product Description	Catalog #
Surgilance™ Yellow 1.0mm Needle Lancet, Low Blood Flow 5-10 µl (100/box)	SLN-100
Surgilance™ Gray 1.8mm Needle Lancet, Low Blood Flow 10-20 µl (100/box)	SLN-200
Surgilance™ Orange 2.2mm Needle Lancet, Medium Blood Flow 20-40 µl (100/box)	SLN-240
Surgilance™ Pink 2.8mm Needle Lancet, High Blood Flow 40-60 µl (100/box)	SLN-300
Surgilance™ Green 1.8mm Blade, High Blood Flow 75-100 µl (100/box)	SLB-200
Surgilance™ Blue 2.3mm Blade, High Blood Flow 150-200 µl (100/box)	SLB-250
Reimbursement Information:	
CPT Code #	National Limit Amount
N/A	N/A

Rev. 02.03.10

Surgilance™ Safety Lancets Test Instructions



HOW TO USE THE SURGILANCE SAFETY LANCET:

SurgiLance makes blood sampling easier for the clinician and more comfortable for the patient.

1. Simply clean test site with isopropyl alcohol and wipe dry.
2. Remove protective cap from the safety lancet.
3. Place red raised platform end of safety lancet on the test site.
4. Gently push the safety lancet down against test site to activate the lancet mechanism.

IT'S THAT EASY! AND SAFE!

WHY THE SURGILANCE SAFETY LANCET IS SUPERIOR:

Simple. No arming is required. Lancets are color-coded for easy recognition.

Safe. The needle or blade lancet is safely concealed before and after use.

Innovative. SurgiLance safety lancet makes blood sampling the easiest and safest method. Once the lancet is used, it is rendered inoperable providing added safety for patient and clinician.

Efficient. The high-speed delivery and penetration method minimizes patient pain and clinician error.

Latex Free. SurgiLance Safety lancet is hypoallergenic.

Rev. 02.03.10