DRAFT – PENDING BAEBIES APPROVAL

FINDER® Multifunctional Near-Patient Diagnostic Platform

Multifunctional Testing – One Platform, Multiple Types of Tests, Delivering Rapid Results for Near-Patient Use.

Baebies FINDER is a platform for near-patient testing in multiple settings. Powered by digital microfluidics, FINDER addresses the critical need for tests requiring low sample volumes, delivering rapid results (positive results in as fast as ~15 minutes) after sample introduction. With its small footprint, cost-effective design, and intuitive workflow, FINDER delivers efficient testing without compromising accuracy.

Rapid testing • Accurate results • Portable platform • Robust pipeline

FINDER Features

- < 1-hour for lab installation
- Toaster-sized instrument
- Fully automated sample preparation
- Preloaded reagents on cartridge
- Gold standard RT-PCR
- Intuitive 3-step workflow
- Positive results in as fast as ~15 minutes
- Built-in Internal Quality Control (IQC)
- Low maintenance – no pumps, tubing, or valves



Built for Your Practice

Designed for urgent care centers, primary care and pediatric clinics, outpatient and occupational health clinics, and independent physician groups – lab-quality molecular answers during a single patient visit. FINDER is a multifunctional platform, and additional assays are in development to expand your point-of-care menu without adding instruments.

DRAFT – PENDING BAEBIES APPROVAL

FINDER® Flu A&B/SARS-CoV-2 Test

FDA 510(k) Cleared and CLIA-Waived

The FINDER® Flu A&B/SARS-CoV-2 Test delivers gold-standard RT-PCR testing for Influenza A, Influenza B, and SARS-CoV-2 right at the point of care. With a compact, automated system powered by digital microfluidics, FINDER® is an ideal solution for urgent care centers, clinics, and other CLIA-waived settings.

Results: positives in as fast as ~15 minutes; negatives in ~20 minutes.



Multiplex Testing
Detects and differentiates Flu A, Flu B, and SARS-CoV-2 from a single nasopharyngeal swab

Gold-Standard RT-PCR
Multiple gene targets ensure detection and robust variant coverage

Ultra-rapid RT-PCR
Positive results in as fast as 15 minutes within a fully enclosed cartridge

Specimen Adequacy Controls
Includes testing for endogenous RPP30 to check specimen quality

Intuitive 3-Step Workflow
No specialized training required and minimal hands-on time

Built-In Quality Controls
Over 100 integrity checkpoints on-cartridge to ensure reliable results

Clinical Performance

Target	Positive Percent Agreement (95% CI)	Negative Percent Agreement (95% CI)
SARS-CoV-2	97.6% (40/41) 87.4–99.6%	99.8% (484/485) 98.8–100.0%
Influenza A	95.4% (62/65) 87.3–98.4%	97.2% (449/462) 95.2–98.3%
Influenza B	96.0% (24/25) 80.5–99.3%	99.4% (499/502) 98.3–99.8%

Prospective clinical study, nasopharyngeal swab specimens, compared with an FDA-cleared RT-PCR method. Source: FDA 510(k) K252269 decision summary.

Coding & Reimbursement	
CPT code	87636-QW
Method	RT-PCR, multiplex amplified probe
CLIA waiver	CW250011 (510(k) K252269)
Medicare national limit	$142.63

Distributed by **CLIAwaived, Inc.** · 2721 Loker Ave W, Carlsbad, CA 92010 888-882-7739 · sales@cliawaived.com · www.cliawaived.com

DRAFT – PENDING BAEBIES APPROVAL

Ordering Information

Everything you need to run the FINDER® Flu A&B/SARS-CoV-2 Test. Contact CLIAwaived for pricing and instrument placement options.

Order #	Description	Pack Size
BAE-100-000082	**FINDER® Instrument** Initial instrument – placement program available, ask us for details	1 instrument
BAE-190-000018	**FINDER® Flu A&B/SARS-CoV-2 Test** (cartridges)	Box of 60
BAE-100-000076	FINDER® RVP Buffer	Box of 60
BAE-100-000077	FINDER® Flu A&B/SARS-CoV-2 Positive Control Swabs	Bag of 2
BAE-100-000078	FINDER® Flu A&B/SARS-CoV-2 Negative Control Swabs	Bag of 2
BAE-100-000061	FINDER® Cleaning Accessory	Box of 10
BAE-999-000032	COPAN FLOQSwabs® (nasopharyngeal)	Box of 100
BAE-100-000014	FINDER® Printer Accessory	1 each
BAE-999-000013	FINDER® Thermal Printer Replacement Paper	Qty of 10
BAE-100-000082-IR	FINDER® Instrument Replacement	1 instrument
BAE-100-000082-EW	FINDER® Instrument Extended Warranty	1 each

Part numbers provided by Baebies, Inc. (manufacturer of FINDER®), October 2026.

Getting Started Is Simple

1
Contact us
Call or email CLIAwaived to discuss your testing volume and placement options.

2
Install & train
Installation in under an hour, with onboarding and training by Baebies.

3
Test & reorder
Reorder cartridges, buffer, and controls through CLIAwaived.

DISTRIBUTED BY

CLIAwaived, Inc.
2721 Loker Ave W
Carlsbad, CA 92010

Toll-free: 888-882-7739
Email: sales@cliawaived.com
Web: www.cliawaived.com

The FINDER Flu A&B/SARS-CoV-2 Test is for in vitro diagnostic use, Rx only, with nasopharyngeal swab specimens from individuals with signs and symptoms of respiratory tract infection. Refer to the package insert for complete intended use, limitations, warnings, and precautions. Manufactured by Baebies, Inc., Durham, NC. Baebies and FINDER are registered trademarks of Baebies, Inc. Assays in development are not available for sale. baebies.com/tests/finder-flu-a-b-rsv-sars-cov-2-test/