



piccolo^{xpress} chemistry analyzer

FAST, EASY, ACCURATE ON-SITE CHEMISTRY TESTING

REAL-TIME COMPREHENSIVE CHEMISTRY RESULTS FOR INTERNAL MEDICINE

MAXIMIZE YOUR PRACTICE THROUGHPUT

- Complete CLIA waived chemistry results on-site in minutes
- Monitor, diagnose and treat in the span of normal patient visit
- See more patients per day at a higher value to the practice
- Lab-accurate liver, kidney, metabolic and lipid function testing
- Increase reimbursement and evaluation revenues


MINIMIZE INEFFICIENCIES AND WASTE

- Reduce time spent calling patients to relay lab results
- Diminish non-billable hours
- Boost staff productivity
- Reduce unnecessary office visits
- Make room for new patients and accompanying revenue

Information Provided By: **CLIAwaived.com**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com

FAST AND DECISIVE

From draw to printed results in 3 easy steps and approximately 12 minutes



1 / Add sample



2 / Insert disc



3 / Read results

SIMPLE AND SMART

Fully automated system requires no special operating skills

- Self contained reagent discs deliver complete results in real time
- Very small sample – only 100 microliters of whole blood, serum or plasma
- As easy to use as a CD player – load sample, insert disc and read results
- Barcoded discs self calibrate with every run and minimize data input errors
- CLIA waived – any of your clinicians can use it
- Comprehensive chemistry results in a single use disc
- Small footprint occupies very little counter space
- 5 bi-directional USB ports for easy connectivity to EMR/LIS or peripherals

ACCURATE AND DEPENDABLE

Amazing technology ensures accurate results every time

- Precision and specificity comparable to expensive high-throughput analyzers
- Long 12-month shelf life ensures minimal waste
- Use discs directly from the refrigerator – no warming required
- iQC™ system provides continuous quality control checks on every run to eliminate the potential for erroneous results
- Panels are fully reimbursable through medicare and private payors



3240 Whipple Road, Union City, CA 94587 | 1-800-822-2947 | www.abaxis.com

©2010 Abaxis, Inc., Piccolo Xpress™ is a trademark of Abaxis, Inc. 888-3222 Rev.A