

NEW

# Control™

Blood Glucose Monitoring System

Biosensor Technology for fast, accurate, and reliable results.

- **Fast Results in 5 seconds**  
*Saves time and effort*
- **Tiny Sample Size (1uL)**  
*Means less pain*
- **Capillary Action**  
*Quickly draws blood into the strip*
- **Immediate on/off with strip insertion**  
*Automatic start*
- **250 - Test Memory**  
*Helps keep track of the results*
- **7-90 - Day Averaging**  
*Allows flexible monitoring*
- **5,000 - Test Battery Life**  
*Makes Control™ extremely reliable*
- **Comprehensive Lifetime Warranty**  
*Assures easy meter replacement*



Product Name	Product Number
Control™ Meter	G3-103MK
Control™ Starter Kit	G3-003SK
Control™ Test Strips - 50/box	G3-203S10
Control™ Control Solution (Level Low)	G3-CSL
Control™ Control Solution (Level Normal)	G3-CSN
Control™ Control Solution (Level High)	G3-CSH

Distributed by: **CLIAwaived.com™**  
San Diego, CA 92121  
tel 858-481-5031  
toll free 888-882-7739  
www.cliawaived.com