

# BACtrack® Select S75 Pro Breathalyzer Test

Distributed by

CLIAwaived™ Inc.



The NEW 2010 BACtrack S75 Pro Breathalyzer Test is the most affordable BACtrack® fuel-cell model ever available! The BACtrack® S75 Pro Breathalyzer Test is a stunning combination of accuracy, affordability, and style. The S75 has been designed to reduce unnecessary size and weight, while maintaining all the accuracy benefits of a full-sized fuel-cell-based breathalyzer. And because it is priced affordably and has obtained FDA 510(k) clearance, this professional technology is now directly available for both consumer and professional use.

#### BENEFITS:

- Simple to use breath intake area
- Versatile design allows for quick use with or without mouthpiece
- 3 digit multifunction display
- Easy to read test results and alerts
- One button operation
- Pocket size compact design



Ordering Information:	
Product Description	Catalog #
BACtrack® Select S75 Pro Breathalyzer (Breathalyzer, 6-mouthpieces, case, 2-“AA” batteries)	BAC-S75
BACtrack® Select Series Mouthpieces (10 pack)	BAC-MPS-10
BACtrack® Select Series Mouthpieces (50 pack)	BAC-MPS-50
Reimbursement Information:	
CPT Code #	National Limit Amount
N/A	N/A

All product information and prices are subject to change without notice at any time. Please call our office before ordering to verify current prices.

# BACtrack® Select S75 Professional Breathalyzer Test Specifications

## SPECIFICATIONS

Dimensions	2.3 x 4.8 x 1 inches (5.8 x 12.2 x 2.5 cm)
Weight	4.4 oz (125 g) with mouthpiece and batteries
Sensor Technology	Fuel cell sensor technology
Detection Range	0.000 – 0.400 %BAC
Test Count	Displays total number of tests performed
Power Supply	Two AA alkaline batteries, included
Battery Life	Approximately 1500 tests
Warm Up Time	10-20 seconds
Response Time	10 seconds
Sensor Accuracy	+/- 0.005 %BAC @ 0.1 %BAC
Operating Temperature	32-104 °F (0-40 °C)

**Warranty:** Limited 1 year warranty on labor and parts

**Return policy:** Due to hygiene purposes, this item is non-returnable.

## PRECAUTIONS

1. Wait at least 20 minutes after drinking, eating or smoking before testing. Not observing this waiting period can cause inaccurate readings and damage the sensor.
2. Do not blow smoke, food, or liquids into the BACTRACK as this will damage the sensor.
3. Do not test in areas with strong winds, smoke, or in an area where a large amount of alcohol is being consumed.
4. Replace the two AA batteries when the battery indicator icon reaches one bar.
5. Send your tester in for periodic calibration service as required. (See calibration).
6. The BACTRACK Breathalyzer is designed to be used in a temperature range of 32-104 °F (0-40 °C).
7. Avoid testing in the presence of any substances that contain methyl alcohol, isopropyl alcohol or acetone. These substances may interfere with the results of the test.
8. You cannot use the results of this product in court.
9. Do not use the BACTRACK as a tool to determine whether you should operate a motor vehicle or equipment, or perform any other dangerous act.
10. Do not drink and drive. Always have a designated driver when alcohol is being consumed.