

The New Next Generation™ Reagents



Unique formulas that get the best performance from a cell analyzer

Boule Medical provides hematology reagents for Medonic CA series and M-series instruments respectively. Due to unique instrument and reagents design, only two kinds of reagents (diluent and lyse) are needed using on Medonic hematology analyzer, which achieve good cost per test scale.

For Human Application

Available for **Medonic CA** analyzers

1504082	Medonic CA Diluent	20 litre
1503971	Medonic CA Diluent	10 litre
1503972	Medonic CA Lyse	10 litre
1504083	Medonic CA Lyse	5 litre

Available for **Medonic M-series**

1504128	M-series Dual pack (Diluent+lyse)
1504122	Medonic M-series Diluent
1504123	Medonic M-series Lyse

For Veterinary Application

1503974	Mediton II vet	10 liter
1503975	Medilyse II vet	5 liter

Distributed by: **CLIAwaived.com™**
 San Diego, CA 92121
 tel 858-481-5031
 toll free 888-882-7739
 www.cliawaived.com