

# CDS 3PD Hematology Control

KIT LOT:

0812-502

EXP:

27-Apr-09

Open Vial Stability:

14 days

LOW CONTROL		Abbott CELL-DYN Series				± Range
		1400	1600	1700	1800	
Lot 0812-011						
Parameters		Mean	Mean	Mean	Mean	± Range
WBC	x10 <sup>9</sup> /L	3.7	3.7	3.7	3.6	0.3
Lymph #	x10 <sup>9</sup> /L	1.1	1.1	1.1	1.1	0.5
Mid #	x10 <sup>9</sup> /L	na	0.2	0.2	0.3	0.3
Gran #	x10 <sup>9</sup> /L	2.6	2.5	2.3	2.3	0.8
Lymph	%	29.1	29.1	31.0	29.3	5.0
Mid	%	na	4.6	6.2	8.2	4.0
Gran	%	70.9	66.3	62.8	62.5	5.0
RBC	x10 <sup>12</sup> /L	2.58	2.58	2.58	2.58	0.12
Hgb	g/dL	6.4	6.4	6.4	6.5	0.3
Hct	%	18.6	18.6	18.6	18.6	2.0
MCV	fL	72	72	72	72	4.0
MCH	pg	24.8	24.8	24.8	25.2	1.5
MCHC	g/dL	34.5	34.5	34.5	35.0	2.2
RDW	%	na	20.2	20.2	19.0	3.0
PLT	x10 <sup>9</sup> /L	80	80	80	79	15
MPV	fL	na	na	na	na	1.5
NORMAL CONTROL		1400	1600	1700	1800	± Range
Lot 0812-012						
Parameters		Mean	Mean	Mean	Mean	± Range
WBC	x10 <sup>9</sup> /L	9.0	9.0	9.0	8.8	0.6
Lymph #	x10 <sup>9</sup> /L	3.6	3.6	3.9	3.6	0.6
Mid #	x10 <sup>9</sup> /L	na	0.5	0.5	0.7	0.5
Gran #	x10 <sup>9</sup> /L	5.4	4.9	4.6	4.5	1.0
Lymph	%	40.2	40.2	43.6	41.0	5.0
Mid	%	na	5.2	5.2	7.9	4.0
Gran	%	59.8	54.6	51.2	51.1	5.0
RBC	x10 <sup>12</sup> /L	4.23	4.23	4.23	4.17	0.18
Hgb	g/dL	12.5	12.5	12.5	12.7	0.4
Hct	%	36.0	36.0	36.0	35.4	2.5
MCV	fL	85	85	85	85	4.0
MCH	pg	29.6	29.6	29.6	30.5	1.5
MCHC	g/dL	34.8	34.8	34.8	35.8	2.2
RDW	%	na	17.1	17.1	15.7	3.0
PLT	x10 <sup>9</sup> /L	245	245	245	232	30
MPV	fL	na	na	na	na	1.5
HIGH CONTROL		1400	1600	1700	1800	± Range
Lot 0812-013						
Parameters		Mean	Mean	Mean	Mean	± Range
WBC	x10 <sup>9</sup> /L	22.1	22.1	22.1	21.5	1.8
Lymph #	x10 <sup>9</sup> /L	11.3	11.3	12.0	11.1	1.5
Mid #	x10 <sup>9</sup> /L	na	1.6	1.5	2.1	1.0
Gran #	x10 <sup>9</sup> /L	10.8	9.2	8.6	8.3	2.5
Lymph	%	51.2	51.1	54.3	51.6	5.0
Mid	%	na	7.4	6.7	9.6	4.0
Gran	%	48.8	41.5	39.0	38.8	5.0
RBC	x10 <sup>12</sup> /L	5.29	5.29	5.29	5.17	0.22
Hgb	g/dL	17.8	17.8	17.8	18.0	0.5
Hct	%	51.3	51.3	51.3	49.6	3.0
MCV	fL	97	97	97	96	4.0
MCH	pg	33.6	33.6	33.6	34.8	1.5
MCHC	g/dL	34.7	34.7	34.7	36.3	2.2
RDW	%	na	16.3	16.3	15.1	3.0
PLT	x10 <sup>9</sup> /L	599	599	599	568	60
MPV	fL	na	na	na	na	1.5

Note: Platelet/RBC - causing alarms, and WBC differential "R" flags should be disregarded for Abbott Series.



Distributed by: **CLIWAIVED.com™**  
 San Diego, CA 92121  
 tel 858-481-5031  
 toll free 888-882-7739  
 www.cliwaived.com

PN 201052B R04.24.07

