

CDS Hematology Calibrator

KIT LOT: 1008-731

EXP: 19-Oct-10

Open Vial Stability: 5 days

Analyzer	WBC x10 ⁹ /L	RBC x10 ¹² /L	Hgb g/dL	Hct %	MCV fL	PLT x10 ⁹ /L	MPV fL
Beckman Coulter Large Systems S Plus Series, JT Series, STKR, STKS™, MAXM™	7.9	4.25	12.6	35.5	83.6	188	11.2
Beckman Coulter Small Systems MD, MD II Ac-T™, Ac-T diff™ Series Onyx and T Series	7.9 7.9	4.28 4.28	12.7 12.7	36.0 36.0	84.0 84.0	188 188	11.4 11.4
Abbott Cell-Dyn® 1400, 1500, 1600, 1700 Cell-Dyn® 1800	7.9 7.4	4.42 4.41	12.6 12.7	36.7 36.6	83 83	239 222	--- ---
ABX Micros 45/60	7.9	4.34	12.5	35.7	82.3	237	11.4
Range Limits:	± 0.3	± 0.10	± 0.2	± 1.0	± 2.0	± 15	± 1.5

Range Limits are provided for evaluating analyzer data, and for determination of analyzer reproducibility and operator skill.

All trademarked analyzer names belong to the respective trademark holders.



Distributed by: **CLIAwaived.com™**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com

201044F R12.12.08

