

CLIAwaived™ Inc. Thyroid Test Rapid TSH Cassette

NEW!



CLIAwaived Inc. Thyroid Test For Elevated TSH

The CLIAwaived Inc. Thyroid Test for Elevated TSH was developed by doctors, to test for elevated thyroid stimulating hormone (TSH) in whole blood. Elevated TSH is the first clinical indicator of hypothyroidism. The CLIAwaived Inc. Thyroid Test is a rapid immunoassay that is easy to use, cost effective and accurate. It is ideal for physician office labs, pharmacy programs and health promotion. Knowing who's not hypothyroid lets you make timely decisions in your office, saving you and your patients time and money. The test requires a drop of blood from a finger stick or venous draw, with results in minutes.

BENEFITS:

- Test patient on first visit
- Staff members perform test
- Results in 10 minutes
- High precision, sensitivity, specificity at 5 mIU/L
- Store at room temperature
- FDA Cleared and CLIA waived

Ordering Information:

Product Description	Catalog #
CLIAwaived™ Inc. Thyroid Test - Rapid TSH Cassette	CLIA-TSH-I20-A
CLIAwaived™ Inc. Thyroid Test Controls	CLIA-TSH-I20-C

Reimbursement Information:

CPT Code #	National Limit Amount
84443QW Hyperthyroidism (TSH)	\$ 24.53

CLIAwaived™ Inc. Thyroid Test - Rapid TSH Cassette Instructions



Two pink lines at C and T indicate a positive result. TSH greater than 5 mIU/L, the first clinical indicator of hypothyroidism, the most prevalent disease of the thyroid.

US Patent 5821073

Procedure:



Place one drop of blood from a finger-stick or venous draw in S1



Wait 90 seconds



Place 4 drops of buffer in S2



Read results within 10 minutes and before 15 minutes

Time saving:

Qualify your patient's thyroid health at the first visit.

Cost effective:

Priced to make it worthwhile to test for TSH first.

Easy to Use:

One drop of blood from a finger-stick or venous draw.

Rapid Results:

Know in 10 minutes or less if TSH is elevated.

Precision:

100% agreement with lab, 98.5% sensitivity, 96.9 specificity at 5mIU/L +/- 10%.

CPT Code:

84443QW

Reimbursement:

\$ 23.47 in most States.

Storage:

Room temperature. Refrigeration not required.

Shelf life:

18 months on all components.