

Now test for all the 'usual suspects' and alcohol

with

DRUGCHECK[®]
NxStep OnSite

Shown with adulterant and alcohol test strips



NEW!
1st Urine Alcohol Test Cup!

Available with test strips for alcohol and adulterants, or alcohol only



Alcohol only

NEW!

1st Urine Alcohol Test Cup!

DRUGCHECK[®]

NxStep OnSite

15 tests in one



- ✓ 15-panel test
- ✓ Available with alcohol and/or adulterant tests
- ✓ 8 oz. cup
- ✓ 9% larger opening
- ✓ Unbreakable cup surface

- ✓ Leak-proof transport technology
- ✓ 15 individual chambers eliminate drug interpretation confusion
- ✓ Tamper-proof strip insert
- ✓ Made in U.S.A.
- ✓ Results in 5 minutes




More great innovations from the leader

The manufacturer of the DrugCheck[®] Drug Test Cup and Dip Drug Test proudly introduces the first 15-panel drug test cup on the market. With the new NxStep OnSite Drug Test Cup, you can now test for up to 15 drug types in a single cup with **results in 5 minutes**.

Designed for the person giving the test

With its superior design, the DrugCheck[®] NxStep OnSite Drug Test Cup requires no handling or manipulation to activate testing. The NxStep OnSite is **FDA cleared** and has the same accuracy when compared to certified lab testing.

- ✓ NO ASSEMBLY
- ✓ NO TILTING TO ACTIVATE
- ✓ NO SPLIT SAMPLES OR CHAMBER TO RELY ON
- ✓ NO CARD DIPPING
- ✓ NO CROSS CONTAMINATION
- ✓ NO MOVING PARTS



For use with alcohol test strips



For use with adulterant test strips

AMP AMP AMP	Amphetamine	1000 ng/ml
BAR BAR BAR	Barbiturate	300 ng/ml
BZD BZD BZD	Benzodiazepine	300 ng/ml
COC COC COC	Cocaine	300 ng/ml
MET MET MET	Methamphetamine	1000 ng/ml
MTD MTD MTD	Methadone	300 ng/ml
OPI OPI OPI	Opiates	300/2000 ng/ml
OXY OXY OXY	Oxycodone	100 ng/ml
PCP PCP PCP	Phencyclidine (PCP)	25 ng/ml
THC THC THC	Marijuana	50 ng/ml
MDMA MDMA MDMA	Ecstasy	500 ng/ml
PPX PPX PPX	Propoxyphene	300 ng/ml
TCA TCA TCA	Tricyclic Antidepressants	1000 ng/ml



Distributed by: **CLIAwaived.com[™]**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com