

## Hyperbaric QC

For Validation of elevated  $pO_2$  Blood Gas Measurements



### Ideal for Analytical Measurement Range (AMR) Validation:

- Access accuracy and precision of elevated  $pO_2$  measurements
- Extend your lab's reportable range
- Allows reporting of actual results--simplify recordkeeping
- Satisfy federal and accrediting agency AMR requirements
- Ideal for cardiac catheterization labs & OR suites
- Three year shelf-life at room temperature (15 - 30°C)

<u>Level</u>	<u>Target</u>	<u>Packaging</u>
High	460mmHg	4 amps/box
X-High	710mmHg	4 amps/box