

## Glucotrol-NG

### Superior QC for HemoCue B-Glucose & Glucose 201 Systems



No other QC offers so many benefits:

- **15-month refrigerated shelf-life** – ensures long term QC continuity and virtually eliminates product waste due to outdating
- **Better consistency & performance** – plasma-like matrix outperforms ‘pickled’ or stabilized RBC formulations
- **Hypo- and Hyperglycemic target values** – select from 5 levels to validate your lab’s entire reportable range
- **Optimized fill volume** – 20 tests per vial reduces risk of using expired product
- **Recommended and distributed by HemoCue**

Best of all, Glucotrol-NG is competitively priced and available from stock for immediate delivery!

<u>Level</u>	<u>Target</u>	<u>Packaging</u>
Level 1	~38 mg/dL	2 vials/box
Level 2	~105 mg/dL	2 vials/box
Level 3	~168 mg/dL	2 vials/box
Level 4	~289 mg/dL	2 vials/box
Level 5	~421 mg/dL	2 vials/box