

Is she Pregnant?



- Accurate
- Quick Results
- Easy to Use



CLIAwaivedTM.com

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121
Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739

E-mail: info@cliawaived.com

Visit us on the web at: www.cliawaived.com

Distributed by: CLIAwaived, Inc.



Global Diagnostics for Local Markets

ACON One-Step Pregnancy tests provide accurate results in minutes

Accurate

- >99% compared to other lab tests

Quick Results

- Some positive results can be seen as soon as 1 minute



Easy to Use

- Add sample, wait 3-5 minutes and read results



Strip Test



- Immerse the strip into urine or serum
- Read results at 3 minutes for urine or 5 minutes for serum



Positive Negative Invalid

Device Test


- Add 3 drops of urine or serum
- Read results at 3 minutes for urine or 5 minutes for serum



Positive Negative Invalid

Midstream Test

- Urinate on the Absorbent Tip or dip it into urine
- Read results at 3 minutes



Positive Negative Invalid

Product Name	Specimen	Format	Catalog No.	Sensitivity	Standard Size	
					Kit Size	Dimension
hCG Pregnancy Test	Urine	Strip	FHC-101	25 mIU/mL	50 Tests	4.0 mm*
			FHC-U101	10 mIU/mL		
			FHC-111	25 mIU/mL	25 Tests/Canister**	
			FHC-U111	10 mIU/mL		
			FHC-101H	25 mIU/mL	1 Test	
			FHC-U101H	10 mIU/mL		
		Device	FHC-102	25 mIU/mL	40 Tests	4.0 mm
			FHC-U102	10 mIU/mL	1 Test	
			FHC-102H	25 mIU/mL		
	Midstream	FHC-U102H	10 mIU/mL	1 Test	7.2 mm	
		FHC-103	25 mIU/mL			
	Urine/Serum	Strip	FHC-U103	10 mIU/mL	50 Tests	4.0 mm*
			FHC-201	25 mIU/mL		
			FHC-U201	10 mIU/mL	25 Tests/Canister**	
FHC-211			25 mIU/mL			
Device		FHC-U211	10 mIU/mL	40 Tests	4.0 mm	
		FHC-202	25 mIU/mL			
		FHC-U202	10 mIU/mL			

* Also available in dimensions of 2.5, 3.0, 3.5 or 5.0 mm

** Also available in canisters of 20 or 30

All 25 mIU/mL hCG tests have been cleared by the US FDA

All 25 mIU/mL and 10 mIU/mL hCG tests are CE Marked for sale in the European Community



We also offer one-step professional and/or home strip, device, midstream and/or panel tests for: Ovulation (LH), Menopause (FSH), HIV 1/2, HBV, HBsAg, HBsAb, HBeAb, HBeAg, HBcAb, HCV, *H. pylori*, Strep A, Malaria, Chlamydia, Syphilis, PSA, FOB, AFP, Drugs of Abuse and others.
Blood Glucose Monitoring System Now Available !