

GAS-ISE-HCT QC

The *ONE* Quality Control for i-STAT[®] Analyzers



One QC now provides a full critical care analyte menu:

- pH, gases, electrolytes, metabolites, **and** hematocrit in one ampule
- Use *one* i-STAT cartridge for each control check, not two
- Cuts QC time in half, significantly lowering costs
- Optimized for i-STAT with physiological target values
- Linearity panel for Analytical Measuring Range (AMR) validation

Available in three levels and as a linearity panel:

Level	Packaging
1	12 amps/box
2	12 amps/box
3	12 amps/box
1-5 Linearity	5 amps/box