

BIONINE



BIONINE

Distributed by: **CLIAwaived.com™**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com

Rightest™ GM300

Blood Glucose Monitoring System
Noble Metal Electrode Strip



FDA



swissdesign 




- Comfortable vision - Large display
- Simple installation of Smart Code Key without code number input
- Big capacity - 300 test memory
- Plug & Play - Auto start & off
- Less pain - Only 1.4µL blood
- High precision & Accuracy - Noble metal electrode technology
- PC Link Adapter

Blood Glucose Monitoring System
Rightest™ GM300



Your Rightest System Kit includes:

- Rightest Meter
- 10 pcs of Rightest Test Strips
- Rightest Smart Code Key
- Rightest Control Solution
- Rightest Check Key
- User's Manual
- Getting Started Guide
- Test Strip Package Insert
- Lancing Device
- 10 pcs of Lancets
- Two AAA Batteries
- Carrying Case
- Logbook
- Emergency Information Card



Friendly- large strip hand bar



Less pain-1.4µL blood sample only, automatically drawn in



Rapidly- results less than 8 seconds



Safety- Keep from touching blood sample

Distributed by: **CLIAwaived.com™**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com

Noble Metal Electrode Strip
High Precision and Accuracy