

# OSOM<sup>®</sup>

## Strep A Test

FOR THE DETECTION OF *STREPTOCOCCUS A.* ANTIGEN




FOR USE WITH THROAT SWABS

FEATURES:

- ✓ **VERY ACCURATE;** 96% sensitivity; 98% specificity versus culture
- ✓ **FAST;** Results in 5 minutes
- ✓ **QC INSIDE;** Includes 2 additional test sticks for external QC testing
- ✓ **EASY-TO-READ;** Two-color results
- ✓ **WORLD-CLASS TECHNICAL SUPPORT;** Staffed with certified medical technologists and skilled professionals



TEST PROCEDURE / RESULTS



ORDERING INFORMATION

- OSOM<sup>®</sup> STREP A TEST
- OSOM<sup>®</sup> STREP A TEST

CONFIGURATION

- 50 Tests
- 20 Tests

CATALOG NUMBER

- 141 OR 141E (CE marked)
- 141E-20 (CE marked)



Distributed by: **CLIAwaived.com<sup>™</sup>**  
 San Diego, CA 92121  
 tel 858-481-5031  
 toll free 888-882-7739  
 www.cliawaived.com

