



PRODUCT CODE
HANDYNEB-3145H

HCPCS Code
E0570 or 571

This information is included for your convenience and is not intended, nor should it be considered billing advice.

an ultra compact nebulizer with carrying case, car charger and rechargeable battery pack

- Product Overview** CLIawaived.com offers an ultra compact piston compressor nebulizer that gives home care dealers the ability to provide a light, powerful option for mobile patients with oxygen needs. The unit comes with a handy tote for easy traveling and other convenient options. Includes battery pack, car adapter, adult and pediatric masks.
- Features**
- Ultra compact design goes everywhere
 - Dry piston compressor
 - Lightweight—less than 1 pound
 - Continuous or intermittent controls
- Warranty** 3-year against manufacturer defects

METHOD OF NEBULIZATION	dry piston compressor
OPERATION VOLTAGE	12 volts
CLASS	Class II
POWER CONSUPTION	1.100A
SOUND NOISE LEVEL (at distance of 1m)	53 dB
JET NEBULIZER MAXIMUM CAPACITY	7 ml
PRODUCT WEIGHT	.77 lbs. (350 g)
UNIT DIMENSIONS	4.8" x 3.4" x 1.9"
	123 x 86 x 48 mm
MEAN MASS AERODYNAMIC DIAMETER (MICRONS)	less than 5 microns
MEDICATION OUTLET	0.25 ml/min. or greater*
OPERATION HUMIDITY AND TEMPERATURE	41° F – 104° F, 30% to 95%
STORAGE HUMIDITY AND TEMPERATURE	33.8° F – 104° F, 10% to 95%
ATMOSPHERIC PRESSURE RANGE	50 – 106 kPa, (500 –1,060 mbar)
FOR CARRIAGE AND STORAGE	

*Value measured with 0.9% NaCl solution. This specification may vary with the kind of the medication being used and it's viscosity. Electric Power Safety Standards EN 13544, EN 60601.1 / Electromagnetic Compatibility as per EN 60601.2

Distributed by: **CLIawaived.com™**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com