



Clearly different.
Clearly better.

Women's Health | Infectious Disease | Cardiovascular | Drugs of Abuse | GI Disorders

Clearview® Chlamydia

Clearly different.

Clearview Chlamydia uses a patented, innovative technology to provide highly sensitive and specific results.

- **Fast:** test provides positive results in just 15 minutes
- **Easy:** assay can be easily performed and interpreted by all staff members
- **Reliable:** highly sensitive assay with built-in control ensures accuracy
- **Flexible:** test provides results while the patient waits



Clearly better

Diagnosis with Clearview Chlamydia enables immediate and necessary treatment.

- Treatment may eliminate the development of serious, potentially life-threatening health problems
- Detection and proper treatment may prevent transmission to an unborn child as well as to sexual partners



Clearview
Chlamydia
 20 Test Kit
 Order No. 135301

CLIA Non-waived

Chlamydia Facts: ¹

- Approximately 75% of women have no symptoms
- Nearly one half of women with untreated infection will develop pelvic inflammatory disease (PID)
- 20% of women with PID will become infertile
- 9% of women with untreated infection will have a life threatening tubal pregnancy
- Infants born to infected mothers may develop conjunctivitis and/or pneumonia

Clearview Chlamydia is a rapid immunoassay for the detection of the Chlamydia trachomatis antigen. It can be performed alongside the patient, with results available during the same patient visit. **Such a simple diagnostic tool is now an essential tool in the effective management of Chlamydia.**

Rapid, on-the-spot detection of genital chlamydial infection in women

Clearview Chlamydia Procedure

