

**BinaxNOW<sup>®</sup>**

***Accurate, Rapid  
Results Mean  
Focused Therapy...  
NOW***



Respiratory | Women's Health & STD | Enterics & GI | General Health | Instrumentation

## **BinaxNOW<sup>®</sup> Legionella**

The BinaxNOW Legionella Urinary Antigen Test provides accurate, rapid identification of Legionella pneumophila serogroup 1 with an easy-to-use technology, aiding physicians in providing rapid, focused therapy.

- **Highly accurate** – 95% sensitivity and specificity; treat with confidence
- **Easy-to-use** – Easy to manage and train; minimizes operator error
- **Results in 15 minutes or less** – Faster throughput in lab or Emergency Room; aiding in rapid, focused therapy
- **Room temperature storage** – Doesn't take up valuable refrigerator space

## **Cost-effective Identification of Legionella pneumophila for Pneumonia and Legionnaires Disease**

- **Decrease mortality rate<sup>2</sup>**
- **Decrease number of broad-spectrum antibiotics prescribed**
- **Shorten patient wait time/hospital stay**
- **Decrease need for unnecessary ancillary testing**
- **Improve patient outcomes through rapid, focused therapy**

**BinaxNOW Legionella**  
Rapid identification of  
*Legionella pneumophila*  
serogroup 1 in urine



**Urinary Antigen Tests  
Recommended  
as Part of the  
IDSA/ATS CAP  
Guidelines<sup>1</sup>**


**Point of Care Diagnostics**



**BinaxNOW Legionella**  
Order No. 852-000  
22 Test Kit  
CPT Code: 87449  
CLIA: Moderate complexity  
Accessories:  
852-010 Control Swab Pack

## BinaxNOW® Legionella Performance

### Package Insert Data\* - Urine




\*Compared with culture.


See product insert for complete performance characteristics.

### Performing the test is easy.

1. Place sample swab into device.



2. Add reagent to swab, then close device.



3. Read at 15 minutes.



### Reading the test couldn't be easier.

One line = **negative**



Two lines = **positive**



visit

[www.invernessmedicalpd.com](http://www.invernessmedicalpd.com)  
for more information

<sup>1</sup>Mandell, L. et al. *Infectious Diseases Society of America/American Thoracic Society Consensus Guidelines on the Management of Community-Acquired Pneumonia in Adults*. Clinical Infectious Diseases 2007; 44:S27-72.  
<sup>2</sup>Benin, A. et al. *Trends in Legionnaires Disease, 1980-1998: Declining Mortality and New Patterns of Diagnosis*. Clinical Infectious Diseases 2002; 35:1039-46.

BinaxNOW® is a registered trademark of Inverness Medical Professional Diagnostics; Inverness Medical-Binax Inc.

10 Southgate Road 1.800.323.3199  
Scarborough, ME 04074 USA P: 207.730.5700  
[www.invernessmedicalpd.com](http://www.invernessmedicalpd.com) F: 207.730.5710

© 2007 Inverness Medical-Binax Inc.  
MKT000003 Rev. 1 2/21/2007