



Clearly different.
Clearly better.



Women's Health | Infectious Disease | Cardiovascular | Drugs of Abuse | GI Disorders

Clearview® Strep A Exact II

Clearview Strep A Exact II

Rapid technology provides direct detection of Group A Streptococcal Antigen.



Clearly different

Clearview's patented innovative technology provides highly sensitive and specific tests for Group A Strep.

- **Rapid Response:** Positive result in as soon as 1 minute.
- **Simplicity:** CLIA waived, easy-to-use dipstick test.
- **Reliability:** Built-in control ensures accuracy.
- **Flexibility:** Results while patient waits.

Clearly better

Clearview Strep A Exact II provides physicians with fast, accurate results for an early diagnosis than can result in:

- Shorter course of illness.
- Reduced severity of symptoms.
- Reduced spread of bacteria to others.
- Reduced likelihood of complications.




Clearview Strep A Exact II
30 Dipsticks
Product No. 92002
CLIA Waived

Clearview Strep A Exact II
30 Cassettes
Product No. 92003
CLIA Non-waived


Clearview Strep A Comparison Grid

Product Name	Clearview Strep A Exact II Dipstick	Quidel QuickVue® Dipstick Strep A	Beckman Coulter ICON® Strep A	Genzyme OSOM® Ultra Strep A
Format	30 Dipsticks	50 Dipsticks	25 Dipsticks	25 Dipsticks
CLIA Status	Waived	Waived	Waived	Waived
Sample	Swab	Swab	Swab	Swab
Sensitivity	97%	92%	97%	92.6%
Specificity	95%	98%	95%	92.8%
Procedural Steps	5	5	5	10
Time to Result	5 minutes	5 minutes	5 minutes	7 minutes
Shelf Life	24 months	18 months	12 months	12 months
Storage	Room Temperature	Room Temperature	Room Temperature	Room Temperature

Clearview Strep A Exact II Dipstick Procedure



Clearview Strep A Exact II Cassette Procedure



visit

www.invernessmedicalpd.com
for more information

All trademarks referenced are trademarks of their respective owners. All rights reserved.

331 S. 104th Street
Louisville, CO 80027 USA
www.invernessmedicalpd.com

Tel: 800.637.3717
303.530.3888
Fax: 800.637.3827

© 2006 Inverness Medical–BioStar Inc.
12/06 02901.01