

# BinaxNOW® Influenza A/B Test

Distributed by

CLIAwaived™ Inc.

## BinaxNOW® Influenza A/B Test

The BinaxNOW Influenza A & B test identifies A or B type in positive samples in 15 minutes. It facilitates better patient management decisions and enables patients to be diagnosed and treated in one visit. The BinaxNOW Influenza A & B test determines quickly if the influenza is the cause of the outbreak. It allows for rapid initiation of antiviral treatment and limits the spread of nosocomial infection. *(Not approved for home use!)*




### BENEFITS:

- Differentiated results for type A and type B
- Simple 3-step procedure
- Easy to read test results
- Results in 15 minutes
- Room temperature storage
- 15-21 month shelf life
- Built-in procedural control features
- FDA Cleared and CLIA waived



<b>Ordering Information:</b>	
Product Description	Catalog #
BinaxNOW® Influenza A/B Tests (22 Tests)	INV-416-022
BinaxNOW® Influenza A/B Tests (10 Tests)	INV-416-110
<b>Reimbursement Information:</b>	
CPT Code #	National Limit Amount
87804QW*	\$16.76

\*It is always the provider's responsibility to establish appropriate coding & charges. Reimbursement shown indicates National Average. It is recommended provider bill 87804QW for Influenza A and 87804QW for Influenza B. Total combined reimbursement \$33.52.

# BinaxNOW® Influenza A/B Test Instructions



Product Name	BinaxNOW® Influenza A & B
Format	Patented card format
Sample Types	NP swab Nasal wash / aspirate Nasal swab <sup>2</sup>
Differentiated Results	YES
Reagents	0
Timed Steps	1
Tests per Kit	22
Kit Storage	Room temperature
Shelf Life	21 months <sup>3</sup>
Controls	2 controls 1 A/B positive, 1 viral negative
QC Requirements	Once per shipment <sup>6</sup>



Rev. 10.10.09