



Clearly different.
Clearly better.

Women's Health | Infectious Disease | Cardiovascular | Drugs of Abuse | GI Disorders

Clearview® Strep A Exact II

Clearview Strep A Exact II

Rapid technology provides direct detection of Group A Streptococcal Antigen.



Clearly different

Clearview's patented innovative technology provides highly sensitive and specific tests for Group A Strep.

- **Rapid Response:** Positive result in as soon as 1 minute.
- **Simplicity:** CLIA waived, easy-to-use dipstick test.
- **Reliability:** Built-in control ensures accuracy.
- **Flexibility:** Results while patient waits.

Clearly better

Clearview Strep A Exact II provides physicians with fast, accurate results for an early diagnosis than can result in:

- Shorter course of illness.
- Reduced severity of symptoms.
- Reduced spread of bacteria to others.
- Reduced likelihood of complications.

Point of Care Diagnostics

CLIAwaived.com

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121

Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739

E-mail: info@cliawaived.com

Visit us on the web at: www.cliawaived.com

Distributed by: CLIAwaived, Inc.




Clearview Strep A Exact II
30 Dipsticks
Product No. 92002
CLIA Waived

Clearview Strep A Exact II
30 Cassettes
Product No. 92003
CLIA Non-waived


Clearview Strep A Comparison Grid

Product Name	Clearview Strep A Exact II Dipstick	Quidel QuickVue® Dipstick Strep A	Beckman Coulter ICON® Strep A	Genzyme OSOM® Ultra Strep A
Format	30 Dipsticks	50 Dipsticks	25 Dipsticks	25 Dipsticks
CLIA Status	Waived	Waived	Waived	Waived
Sample	Swab	Swab	Swab	Swab
Sensitivity	97%	92%	97%	92.6%
Specificity	95%	98%	95%	92.8%
Procedural Steps	5	5	5	10
Time to Result	5 minutes	5 minutes	5 minutes	7 minutes
Shelf Life	24 months	18 months	12 months	12 months
Storage	Room Temperature	Room Temperature	Room Temperature	Room Temperature

Clearview Strep A Exact II Dipstick Procedure



Clearview Strep A Exact II Cassette Procedure



visit

www.invernessmedicalpd.com
for more information