



Clearly different.
Clearly better.



Women's Health | Infectious Disease | Cardiovascular | Drugs of Abuse | GI Disorders

Clearview® Mono

Clearview Mono

Two-step technology provides a quick, easy solution for diagnosing Infectious Mononucleosis.



Clearly different

Clearview's patented innovative technology provides highly sensitive and specific tests for Infectious Mononucleosis (IM).

- **Rapid Response:** Two-step test, with results in 5 minutes.
- **Simplicity:** CLIA waived*, easy-to-use test requires no training.
- **Reliability:** Built-in control ensures accuracy.
- **Flexibility:** Use with multiple sample types — whole blood, serum or plasma.

Clearly better

Clearview Mono provides a clear solution to the common challenge of diagnosing Infectious Mononucleosis (IM).

- Immediate therapy reduces time of discomfort for patient.
- Minimizes the risk of more serious illness.
- Reduces the spread of virus to others.

CLIAwaived.com

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121

Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739

E-mail: info@cliawaived.com

Visit us on the web at: www.cliawaived.com

Distributed by: CLIAwaived, Inc.

* for whole blood specimens

Point of Care Diagnostics

 **inverness medical**
professional diagnostics



Clearview Mono
30 Test Kit
Order No. 92401

Whole Blood - CLIA Waived
Serum, Plasma -
Non-waived

Clearview Mono Comparison Grid

Product Name	Clearview Mono	Quidel QuickVue® + Infectious Mono	Genzyme OSOM® Mono
Format	30 Cassettes	20 Cassettes	25 Dipsticks
CLIA Status	Waived: Whole Blood Non-waived: Whole Blood, Serum or Plasma	Waived: Whole Blood Non-waived: Whole Blood, Serum or Plasma	Waived: Whole Blood Non-waived: Whole Blood, Serum or Plasma
Sample	Whole Blood (Fingerstick or Venipuncture), Serum or Plasma	Whole Blood (Fingerstick or Venipuncture), Serum or Plasma	Whole Blood (Fingerstick or Venipuncture), Serum or Plasma
Sensitivity	98% *	99.9%	100%
Specificity	>99% **	99.8%	90.3%
Procedural Steps	2	2	3
Time to Result	5 minutes	5 minutes	5 minutes
Shelf Life	18 months	18 months	12 months
Storage	Room Temperature	Room Temperature	Room Temperature

Clearview Mono Procedure

2 drops Venipuncture whole blood + 1 drop Buffer

OR

50µL Fingerstick whole blood + 1 drop Buffer

OR

1 drop Serum or Plasma + 1 drop Buffer

Read in 5 minutes

Positive **Negative** **Invalid**

Interpret results

CLIAwaived.com

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121
Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739
E-mail: info@cliawaived.com
Visit us on the web at: www.cliawaived.com
Distributed by: CLIAwaived, Inc.

* positive agreement with comparative test
** negative agreement with comparative test

All trademarks referenced are trademarks of their respective owners. All rights reserved.

331 S. 104th Street 1.800.637.3717
Louisville, CO 80027 USA P: 303.530.3888
www.invernessmedicalpd.com F: 800.637.3827

© 2006 Inverness Medical–BioStar Inc.
02902.00

visit
www.invernessmedicalpd.com
for more information

