

# Para 4<sup>®</sup>

## A four parameter assayed hematology control



### Assays

Para 4 is a four parameter hematology control designed for small laboratories or Point of Care facilities performing a limited number of hematology tests, utilizing manual or semi-automated methods. Assays include:

- ITC Hgb Pro Professional Hemoglobin Testing System™
- HemoCue<sup>®</sup> B-Hemoglobin Photometer, Hb 201+ Analyzer
- HemataSTAT II<sup>®</sup>

Para 4 is assayed for RBC, WBC, Phase Platelet (Low, Normal) and Hct manual methods. The technologist can use one product, saving time and money. Para 4 is available in plastic dropper vials for easy sample allocation with hemoglobin and hematocrit instruments. Para 4 is also available in glass vials for automated methods.

### STATS<sup>®</sup>

Streck offers a free quality control program to all qualifying customers. STATS provides convenient data entry, quick turnaround time and relevant reports. Data can be submitted by mail, fax or e-mail.

### Typical assay values for Para 4

RBC X 10 <sup>12</sup> /L	2.00 - 5.50
WBC 10 <sup>9</sup> /L	2.0 - 22.0
HGB g/dL	6.1 - 16.5
HCT %	17 - 50
PLT x 10 <sup>9</sup> /L	60 - 255

### Stability for Para 4

Closed-vial stability	110 days
Open-vial stability	14 days

Laboratories using Para 4 to control spun hematocrits, or values on the HemoCue, can store vials at 18°-30°C for the 14-day open-vial stability.

Level	Description	Catalog #
<b>Glass vials</b>		
6 Low	6 x 1.5 ml	215426
6 Normal	6 x 1.5 ml	215416
6 High	6 x 1.5 ml	215436
4 Low, 4 Normal, 4 High	12 x 1.5 ml	215412
6 Low, 6 Normal, 6 High	18 x 1.5 ml	215418
<b>Plastic squeeze dropper vials</b>		
6 Low	6 x 1.5 ml	215106
6 Normal	6 x 1.5 ml	215107
6 High	6 x 1.5 ml	215108
2 Low, 2 Normal, 2 High	6 x 1.5 ml	215110
4 Low, 4 Normal, 4 High	12 x 1.5 ml	215112

**STATS-Link<sup>®</sup>** Streck has added STATS-Link to the current STATS<sup>®</sup> Interlaboratory Quality Control Program. Using a personal security code to access their reports, STATS users now have the ability to view and download their STATS reports online by visiting [www.streck.com](http://www.streck.com). STATS-Link participants will receive e-mail notification that their report is available for access online, eliminating the wait associated with mailing reports. In addition to current reports, archived reports are available.

**CLIAwaived.com**<sup>™</sup>

For more information please contact CLIAwaived<sup>™</sup> Inc. at:

TEL: (858) 481-5031 TOLL FREE: (888) 882-7739

E-MAIL: [info@cliawaived.com](mailto:info@cliawaived.com) WEB: [www.cliawaived.com](http://www.cliawaived.com)

Distributed by: CLIAwaived<sup>™</sup> Inc., San Diego, CA

**Streck**

Clinical Laboratory Products You Can Count On.