

MASIMO RAD-8™ HOMECARE SENSOR & CABLE CATALOG 2007



RAD-8 Homecare Package: 9139

Includes oximeter, LNCS patient cable, power cord, quick start guide, LNCS adhesive sensor sample 4 pack and choice of LNCS DCI Adult reusable sensor or LNCS Y1 Adult/Pediatric reusable Y sensor.



For Ordering Information Call 1(888) 882-7739 or Phone: (858) 481-5031 or visit us on the Web at www.cliawaived.com
Distributed By: www.cliawaived.com

CLIAwaived.com™

CLIAwaived.com - 4332 Corte de la Fonda - San Diego, CA 92130




Toll Free : (888) 882-7739 - Phone (858) 481-5031

FAX: (801) 720-7568 - www.cliawaived.com




This catalog contains Masimo proprietary material and is provided for information purposes only. No part of this catalog may be reproduced or retransmitted without Masimo's prior written consent. Product availability is subject to change without notice.

Masimo, CleanShield, DCI, LNCS, LNOP, Rad-5, Rad-8, Radical, SET, Signal Extraction Technology and Signal I.Q., are registered trademarks of Masimo Corporation. Rad-5v and Rad-57 are trademarks of Masimo Corporation. Rainbow and SpCO are registered trademarks of Masimo Laboratories Inc. SpMet is a trademark of Masimo Laboratories Inc.








Masimo SET® Pulse Oximetry Devices

Masimo Part Number	Description	UM	Product Image
9019	Rad-8®, Horizontal Comes with either a LNCS Patient Cable, a choice in the language of Operating Manual and a choice in power cord.	EA	
9025 (5) 9026 (5v)	Rad-5®/Rad-5v™ Pulse Oximeter Comes with an LNCS Patient Cable, a choice among available colors of Handheld Boots and a choice in the language of Operating Manual.	EA	
9063	Rad-57™ Pulse CO-Oximeter Includes SpCO, SpMet, SpO ₂ , Pulse Rate and Perfusion Index parameters. Also includes either a 3 foot Rainbow Adult reusable sensor or 3 foot Rainbow Pediatric reusable sensor, a choice between a protective carrying case or handheld boot (both available in multiple colors) and a choice in the Operator's Manual language.	EA	






LNCS® SENSORS

Masimo Part Number	Description	UM	Product Image	Substitute For
1859	LNCS® Adtx Adult Adhesive Sensors, 18 in. Single Patient Use 20/box, Latex free, Sterilizable Weight > 30kg Blue Packaging	Box/20		Nellcor D25 (Blue) Adult adhesive sensor
1860	LNCS® Pdtx Pediatric Adhesive Sensors, 18 in. Single Patient Use 20/box, Latex free, Sterilizable Weight 10kg – 50kg Pink Packaging	Box/20		Nellcor D20 (Orange) Pediatric adhesive sensor
1861	LNCS® Inf-L Infant Adhesive Sensors, 3 ft. Single Patient Use 20/box, Latex free, Sterilizable Weight 3 – 20 kg Orange Packaging	Box/20		Nellcor I20 (Yellow) Infant adhesive sensor


LNCS® SENSORS

Masimo Part Number	Description	UM	Product Image	Substitute For
2328	<p>LNCS® Inf Infant Adhesive Sensors, 18 in. Single Patient Use 20/box, Latex free, Sterilizable Weight 3kg – 20kg Use replacement tape 2307 Orange Packaging</p>	Box/20		<p>Nellcor I20 (Yellow)</p> <p>Infant adhesive sensor *LNCS 2328 has longer tape than I20</p>
1862	<p>LNCS® Neo-L Neonatal Adhesive Sensors, 3 ft. Single Patient Use 20/box, Latex free, Sterilizable Weight < 3 or > 40 kg Purple Packaging</p>	Box/20		<p>Nellcor N25 & Nellcor N25L (Purple)</p> <p>Neonatal/adult adhesive sensor</p>
2329	<p>LNCS® Neo Neonatal Adhesive Sensors, 18 in. Single Patient Use 20/box, Latex free, Sterilizable Weight < 3kg or > 40kg Use replacement tape 2308 Purple Packaging</p>	Box/20		<p>Nellcor N25 (Purple)</p> <p>Neonatal/adult adhesive sensor *LNCS 2329 has longer tape than N25</p>
1901	<p>LNCS® NeoPt-L Neonatal Adhesive Sensors, 3 ft. Single Patient Use 20/box, Latex free, Sterilizable Weight < 1 kg Use replacement tape 2309 Purple Packaging</p>	Box/20		
1798	<p>LNOP® Neo-L LNOP Sensor for sensitive skin for the Respironics Smart Monitor. Single Patient Use, 20/box, Latex free, non-sterile Weight < 3kg or > 40kg Use 1799 replacement tape Purple Packaging</p>	Box/20		
1396	<p>LNOP® Direct Connect Sensor Cable DCSC Adult Reusable Sensor Cable, 3 ft. 1/box, Latex free, Non-sterile Weight > 30kg</p>	EA		<p>Nellcor DS100A Reusable direct connect clip type sensor</p>
1969	<p>LNOP® DC-12 Adult Reusable Sensor Cable, 12 ft. 1/bo, Latex free, Non-sterile Weight > 30kg</p>	EA		



LNCS® SENSORS

Masimo Part Number	Description	UM	Product Image	Substitute For
1863	LNCS® DCI Adult Reusable Sensor, 3 ft. 1/box, Latex free, Non-sterile Weight > 30kg	EA		Nellcor DS100A-1 Adult reusable clip type sensor
1864	LNCS® DCIP Pediatric Reusable Sensor, 3 ft. 1/box, Latex free, Non-sterile Weight 10kg – 50kg	EA		Nellcor Oxi P/ I Nellcor DYSPD Pediatric reusable finger probe
1895	LNCS® TC-I Reusable Tip-Clip Ear Sensor, 3 ft. 1/box, Latex free, Non-sterile Weight > 30kg	EA		
2258	LNCS® YI Multisite Reusable Sensor, 3 ft. 1/box, Latex free, Non-sterile Multiple Foam and Adhesive Wraps Weight > 1kg	EA		Nellcor D-YS Nellcor Oxi-PI Nellcor Oxi-Ni Adult, pediatric, neonatal reusable sensor
1896	LNCS® TF-I Adult Transreflectance Forehead Sensor, 3 ft. 1/box, Latex free, Non-sterile Weight > 30kg	EA		



LNCS® Patient Cables

Masimo Part Number	Description	UM	Product Image
2017	LNC-4 LNCS Patient Cable, 4 ft.	EA	
1814	LNC-10 LNCS Patient Cable, 10 ft.		
2013	LNC-14 LNCS Patient Cable, 14 ft.		

LNCS® Patient Cables



Masimo Part Number	Description	UM	Product Image
2021	LNC-4-Ext LNCS Extension Cable, 4 ft	EA	
1005	PC08 8 foot LNOP patient cable for the Respiration Smart Monitor	EA	

LNCS® Adapter Cables








Masimo Part Number	Description	UM	Product Image	Substitute For
2268	LNC MAC-180 Connects Masimo LNCS Sensors to Nellcor N-180 and N-200 Pulse Oximeters Adapter Cable, 10ft.	EA		Required to use Masimo LNCS sensors with Nellcor pulse oximeters
2269	LNC MAC-395 Connects Masimo LNCS Sensors to Nellcor N-395, N3000, N290, N295, pulse oximeters. Adapter Cable, 10ft.	EA		Required to use Masimo LNCS sensors with Nellcor pulse oximeters

Accessories




Sensor Accessories

Masimo Part Number	Description	UM	Product Image	Substitute For
2307	Replacement Tapes for LNCS Inf Series Sensors For use with LNCS Inf 100/box Latex free, Non-sterile Orange Packaging	Box/100		Nellcor ADH-P/I Respiration 9541
2308	Replacement Tapes For LNCS Neo Series Sensors For use with LNCS Neo / Neo-3 100/box Latex free, Non-sterile Purple Packaging	Box/100		




Sensor Accessories

Masimo Part Number	Description	UM	Product Image	Substitute For
2309	Replacement Wraps for LNCS NeoPt Series Sensors For use with LNCS NeoPt 10/pack Latex free, Non-sterile	Pack/10		Nellcor ADH-P/I Respironics 9541
1596	CleanShield® Multisite Wrap For use with LNCS YI Adult/Pediatric/Neonatal 100/box Latex free, Non-sterile	Box/100		Nellcor ADH/PI Datex Ohmeda OXY/RTB OXY/RT
1597	Standard Multisite Wrap For use with /LNCS YI Adult/Pediatric/Neonatal 100/box Latex free	Box/100		Nellcor ADH/PI Datex Ohmeda OXY/RTB OXY/RT
1598	Standard Petite Wrap For use with LNCS® YI Adult slender digit/Ped/Neonatal adhesive 100/box Latex free	Box/100		Nellcor ADH/PI Datex Ohmeda OXY/RTB OXY/RT
1602	Foam Wraps for LNOP®/LNCS® YI Sensors For use with LNOP/LNCS YI 12/pack Washable, Latex free	Pack/12		Nellcor JT Posey 6554
1601	Accessory Clothing Clips For securing reusable sensor cables to patient clothing 5/pack	Pack/15		
1608	Masimo Headband For use with LNCS® TF-I 5/box	Box/5		

Sensor Packs

Masimo Part Number	Description	UM	Product Image
1897	LNCS® Adtx/Pdtx Sample Pack Comes with 3 LNCS Adtx Adult and 1 LNCS Pdtx Pediatric Sensor Blue Packaging	Pack	
1898	LNCS® Neo-L/Pdtx/Inf-L Sample Pack Comes with 2 LNCS Neo-L Neonatal/Adult, 1 LNCS Pdtx Pediatric, and 1 LNCS Inf-L Infant sensor Purple/Peach/Yellow Packaging	Pack	
1899	LNCS® Sensor Sample Pack Comes with 1 LNCS Adtx Adult, 1 LNCS Pdtx Pediatric, 1 LNCS Inf-L Infant and 1 LNCS Neo-L Neonatal/Adult sensor Blue/Peach/Yellow/Purple Packaging	Pack	

Handheld Device Accessories

Masimo Part Number	Description	UM	Product Image
2207	Water Resistant Protective Carrying Case for Rad-57™, Black Shown with Quick Reference Card, Rainbow® Sensor and Rad-57. Not included.	EA	
2208	Water Resistant Protective Carrying Case for Rad-57™, Red Shown with Quick Reference Card, Rainbow® Sensor and Rad-57. Not included.	EA	
13158	Handheld Carrying Case (Shown with Rad-5. Rad-5 not included)	EA	

Handheld Device Accessories

Masimo Part Number	Description	UM	Product Image
1842	Handheld Boot, Gray For use with Rad-5, Rad-5v and all Rad-57 handhelds	EA	
1980	Handheld Boot, Yellow For use with Rad-5, Rad-5v and all Rad-57 handhelds	EA	
2138	Rad-5/5v Handheld Boot, Yellow	EA	
1981	Handheld Boot, Red For use with Rad-5, Rad-5v and all Rad-57 handhelds	EA	
2139	Rad-5/5v Handheld Boot, Red	EA	
1982	Handheld Boot, Orange For use with Rad-5, Rad-5v and all Rad-57 handhelds	EA	
2097	Handheld Boot, Royal Blue For use with Rad-5, Rad-5v and all Rad-57 handhelds	EA	
2098	Handheld Boot, Light Blue For use with Rad-5, Rad-5v and all Rad-57 handhelds	EA	
2099	Handheld Boot, Pink For use with Rad-5, Rad-5v and all Rad-57 handhelds	EA	