


## Specifications for ThyroChek®



### 120-A

Box of 20 tests, pipettes, buffer, inserts



### 120-B

Self-sealing bag of 20 tests, pipettes, buffer, inserts



### 120-C

Positive and negative controls

- **Shipper (Case):** 20 Bags or 12 boxes 17.7 by 43.2 by 28 cm or 12 by 18 by 20"
- **Storage:** Room temperature 15-30°C/60-86°F / Do not freeze.
- **Expiry dating:** minimum 18 months from date shipped.

**Shipping:** FOB San Diego

**Controls:** Available as Product Number 120-C.

**Proficiency Testing:** registered with major services.

**Specificity:** 96.95%

**Sensitivity:** 98.5%

**Precision:** 97.7%

**CLIA waived:** approved October 4, 2006

**Moderate Complexity:** 510 K 990658

**UPC code:** application made

**CE mark:** application made

**US Patent:** 5821073

**US manufacture:** ISO 13485:1996

**CPT Code:** 84443QW

**Reimbursement:** \$27.47 in most States

## TO PLACE ORDER

- **FAX:** 801-720-7568
- **EMAIL:** [info@cliawaived.com](mailto:info@cliawaived.com)
- **TEL:** 858-481-5031

Payments should be made out and sent to our mailing address:

CLIAwaived, Inc.  
4332 Corte de la Fonda  
San Diego, CA 92130

**CLIAwaived.com**<sup>TM</sup>  
CLIAwaived.com - 4332 Corte de la Fonda - San Diego, CA 92130  
Toll Free : (888) 882-7739 - Phone (858) 481-5031  
FAX: (801) 720-7568 - [www.cliawaived.com](http://www.cliawaived.com)