


1 Insert MEMo Chip.
Turn Analyzer On.

Inserte MEMo Chip.
Encienda el instru-
mento.


Introduza o MEMo
Chip. Ligue o
Analizador.



2 Insert the Test
Strip.

Inserte tira.


Introduza a Tira
de Teste.



3 Stick finger.

Hurga su dedo.


Espete o dedo.



4 Collect sample.

Toma muestra.


Recolha a amostra.



5 Apply blood to
strip.

Aplique muestra.

Aplique o sangue
na tira.



6 Read results.

Lea los
resultados.

Leia os
resultados.

PS-000525 Rev.0 02/2008