

# QUICKVUE<sup>®</sup>

One-Step hCG  
Urine and Combo TEST

## Two Simple Tests You Can Depend On For The Early Detection of Pregnancy



25 Test Urine Kit



50 Test Combo Kit

### Rapid Diagnosis

- Results in 5 minutes (serum), 3 minutes (urine).
- Strong positive may become visible sooner.

### Excellent Performance

- >99% sensitivity; specificity; accuracy.
- 25 mIU/mL sensitivity in serum or urine.

### Convenient

- Room temperature storage.
- 18-month shelf life.

### Easy to Use

- Simple, one-step procedure uses only 3 drops of sample.
- Two-color endpoint for easy-to-read and interpret results.

### Reliable

- Highly sensitive assay with built-in controls.

### Lot-to-Lot Consistency


- Highly controlled and monitored manufacturing environment.

**QUIDEL**<sup>®</sup>  
CORPORATION


For more information on QuickVue One-Step hCG Urine and Combo test kits, visit our website at [www.quidel.com](http://www.quidel.com), or call your Quidel Account Manager at 1.800.874.1517 (U.S. only) or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [www.quidel.com](http://www.quidel.com)

## Combo Procedure:



## Urine Procedure:



# CLIAwaived.com<sup>TM</sup>

CLIAwaived.com - 4332 Corte de la Fonda - San Diego, CA 92130

Toll Free : (888) 882-7739 - Phone (858) 481-5031

FAX: (801) 720-7568 - [www.cliawaived.com](http://www.cliawaived.com)

**QUIDEL**<sup>®</sup>  
CORPORATION

For more information on QuickVue One-Step hCG Urine and Combo test kits, visit our website at [www.quidel.com](http://www.quidel.com), or call your Quidel Account Manager at 1.800.874.1517 (U.S. only) or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [www.quidel.com](http://www.quidel.com)

Distributor Information:

# CLIAwaived.com<sup>TM</sup>

CLIAwaived.com - 4332 Corte de la Fonda - San Diego, CA 92130

Toll Free : (888) 882-7739 - Phone (858) 481-5031

FAX: (801) 720-7568 - [www.cliawaived.com](http://www.cliawaived.com)