

## Aerosol Therapy

# AlphaNeb Bear Nebulizer System

The child-friendly Roscoe Bear Nebulizer provides smoother, quieter treatment with a state-of-the-art piston compressor. Weighing only 3.4 pounds, the Bear Nebulizer delivers 43 psi at maximum pressure and 16 psi operating pressure. The nebulizer system comes complete with compressor, nebulizer, child mask, tubing, mouthpiece and five extra filters and a five year warranty.



### Product Features

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a built-in carrying handle and a convenient nebulizer parking post
- Comes complete with compressor, nebulizer, angled mouthpiece, tubing, five extra filters and instruction manual
- Five year warranty
- HCPCS Code: E0570\*
- HCPCS Code: A7003\* for nebulizer kits
- HCPCS Code: A7013\* for foam filters

### Technical Specifications

- Power: 120 VAC/60Hz
- Power Consumption: 1.5 amps
- Sound Level: less than 60dBA at 1 meter
- Maximum Compressor Pressure: 29 to 43 psi
- Typical Operating Pressure: 10 to 16 psi
- Typical Operating Flow Rate: 5 ~ 8 LPM
- Medication Capacity: 5 ml
- Particle Size: Less than 5 microns
- Operating Temperature: 50° - 104°F
- Dimensions: 7 ¾" (L) x 5 ¾" (W) x 5" (H)
- Weight: 3.4 lb.

Item #	Description	Qty
NEB-BEAR	AlphaNeb Bear Nebulizer Compressor System	each

Distributed by: **CLIAwaived.com™**  
 San Diego, CA 92121  
 tel 858-481-5031  
 toll free 888-882-7739  
 www.cliawaived.com



\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.