

# Aerosol Therapy

## Mini Nebulizer Compressor

Weighing only 2.55 pounds, the compact, lightweight Roscoe Medical Mini Compressor Nebulizer is easy to handle and transport, providing reliable treatment anywhere. The powerful, efficient piston compressor delivers 29 psi at maximum pressure and 10 psi operating pressure. The Roscoe Nebulizer includes compressor, nebulizer, angled mouth piece, tubing, five extra filters and instruction manual.



### Product Features

- Powerful piston pump is efficient, dependable and quiet
- Small lightweight portable design for easy transport
- Produces a particle size of less than 5 microns
- Comes complete with compressor, nebulizer, angled mouthpiece, tubing, five extra filters and instruction manual
- Five year warranty
- HCPCS Code: E0570\*
- HCPCS Code: A7003\* for nebulizer kits
- HCPCS Code: A7013\* for foam filters

### Technical Specifications

- Power: 120 VAC/60Hz
- Power Consumption: 1.5 amps
- Sound Level: less than 60dBA at 1 meter
- Maximum Compressor Pressure: 17.4 to 29 psi
- Typical Operating Pressure: 7.25 to 10 psi
- Typical Operating Flow Rate: 4.5 ~ 7.0 LPM
- Medication Capacity: 5 ml
- Particle Size: Less than 5 microns
- Operating Temperature: 50° - 104°F
- Dimensions: 8 1/4" (L) x 5 2/3" (W) x 3 3/5" (H)
- Weight: 2.55 lb.

Item #	Description	Qty
NEB-ROS	Roscoe Mini Nebulizer Compressor System	each

Distributed by: **CLIAwaived.com™**  
 San Diego, CA 92121  
 tel 858-481-5031  
 toll free 888-882-7739  
 www.cliawaived.com

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

