

Distributed by

**CLIAwaived™ Inc.**

## Product Sheet

### SURVILANCE SAFETY LANCETS

SAFETY LANCET    COLOR    DEPTH    GAUGE    BLOOD FLOW



SLN100

Yellow    1.0mm    G26

Low Flow  
5-10 ul



SLN200

Grey    1.8mm    G23

Low Flow  
10-20 ul



SLN240

Orange    2.2mm    G22

Med Flow  
20-40 ul



SLN300

Pink    2.8mm    G21

Med-High Flow  
40-60 ul



SLB200

Green    1.8mm    G19

High Flow  
75-100 ul



SLB250

Blue    2.3mm    G18

High Flow  
150-200 ul