



CLIA Waived.com

# STAT-Site<sup>®</sup> M<sup>Hgb</sup>

Power in the Palm  
of Your Hand—  
Truly Portable  
Hemoglobin Testing

TM

CLIA'88  
WAIVED

# STAT-Site® M<sup>Hgb</sup> Products

Catalog No. 900900

## STAT-Site® M<sup>Hgb</sup> Meter

Includes STAT-Site® M<sup>Hgb</sup> meter, operating manual, quick reference guide and battery (extra battery also included)



Catalog No. 901025

## STAT-Site® M<sup>Hgb</sup> Test Cards

Includes 100 test cards (4 packages of 25 each) and a lot-specific CODE key



Catalog No. 900905

## STAT-Site® M<sup>Hgb</sup> Battery Pack

Includes 5 batteries (approximately 5,000 tests)

Catalog No. 503000

## STAT-Site® Hemoglobin Controls

6 x 1.5 mL (3 low and 3 high level)

# STAT-Site® M<sup>Hgb</sup> Specs

|                               |                              |
|-------------------------------|------------------------------|
| Methodology _____             | Reflectance photometry       |
| Principle _____               | Azidemethemoglobin           |
| Calibration _____             | Automatic (CODE Key)         |
| Sample Type _____             | Whole blood                  |
| Sample Size _____             | 1 drop (approximately 12 µL) |
| Linearity _____               | 6 to 21 g / dL               |
| Within Run Precision CV _____ | 3.5% to 4.9%                 |
| Total Precision CV _____      | 2.9% to 4.2%                 |

# Easy as one, two, three.

**Step 1 INSERT THE TEST CARD**

**Step 2 APPLY THE SAMPLE**

**Step 3 READ THE RESULT**



CPT CODE: 85018QW