

QuikScreen® Multiple Drug Cup Test

Distributed by

CLIAwaived™ Inc.



QuikScreen is the easy-to-use urinalysis drug test that detects drugs and drug metabolites at SAMHSA cutoff levels. Unlike other tests, QuikScreen does not require any handling or manipulation to activate the test. Just dispense urine into the cup and read results. The cup does all the work in only 5 minutes!

BENEFITS:

- Simple Procedure
- Easy-to-interpret color bands
- Built-in test control
- Room temperature storage
- 12-18 months shelf life
- No instruments needed
- FDA-Cleared for professional use only.
Not CLIA Waived.

Ordering Information:	
Product Description	Catalog #
QuikScreen® Multiple Drug Cup Test - 5 Panel (25 tests) (THC, COC, OPI, AMP, mAMP) FDA-Cleared / CLIA Waived	SYN60502
QuikScreen® Multiple Drug Cup Test - 5 Panel (25 tests) (THC, COC, OPI, AMP, PCP) FDA-Cleared / CLIA Waived	SYN60505
QuikScreen® Multiple Drug Cup Test - 12 Panel (25 tests) (THC, COC150, OPI300, AMP500, mAMP500, PCP, BAR, BZD, MTD, OXY, PPX, .BUP) FDA-Cleared for Professional Use Only. NOT CLIA Waived.	SYN60C12
Reimbursement Information:	
CPT Code #	National Limit Amount
Effective 1/1/2011 the Healthcare Common Procedure Coding System (HCPCS) included the following new code: G0434QW – Drug screen, other than chromatographic; any number of drug classes, by CLIA waived test or moderate complexity test, per patient encounter.	\$20.83

DISCLAIMER – The CPT information contained in this document is provided as representative examples of reimbursement in this category. It is intended to assist providers in accurately obtaining reimbursement for health care services. It is not intended to increase or maximize reimbursement by any payor. Providers should consult their payor organizations with regard to local reimbursement policies. The information provided in this document is for information purposes only and represents no statement, promise or guarantee by CLIAwaived™ Inc. All CPT codes are supplied for information purposes only and represent no statement, promise or guarantee by CLIAwaived™ Inc. that these codes will be appropriate or that reimbursement will be made. CPT codes and descriptions are copyrights of the American Medical Associations. CPT does not include fee schedules, relative values or related ratings. The source for this information is the Center for Medicare and Medicaid Services. The content provided by the Center for Medicare and Medicaid Services is updated frequently. It is the responsibility of the service provider to confirm the appropriate coding required by their local Medicare carriers, fiscal intermediaries and commercial payors. All product information and prices are subject to change without notice at any time. Please call our office before ordering to verify current prices.


QuikScreen® Multiple Drug Cup Test Instructions

- ★ **Accuracy:** FDA 510(K) registered detection strips have the same high degree of accuracy when compared to SAMHSA certified laboratory screening methods.
- ★ **Self-Contained:** Each cup includes a temperature strip and identification label. Each cup is sealed in its own bag and has an 18-month shelf life.
- ★ **No Leakage:** The cup does not leak so a positive specimen can be shipped in the same container to a laboratory for GC/MS confirmation.
- ★ **Cost:** In addition to providing immediate results, QuikScreen® device costs less than half the price of laboratory testing.



Test procedure


1. Urinate into the cup, filling the cup with urine above the minimum level line.
2. Wait 5 minutes.
3. Read the results.
4. Refer to product insert for detailed instructions.



Confirm: One line indicates a positive result.

Negative: Two lines indicate a negative result.

Easy to read results



Compare adulteration results to color chart for devices with built-in adulteration check.

(Above chart is for illustration purpose only and should not be used to compare results.)