

Pediatric (Train) Nebulizer

RESPIRATORY
PEDIATRIC NEBULIZER



PRODUCT CODE

THN-PN3008

HCPCS Code

E0570

This information is included for your convenience and is not intended, nor should it be considered billing advice.

the train nebulizer —adding a touch of child's play to pediatric nebulization

Product Overview

CLIAwaived.com distributes a pediatric piston compressor nebulizer with the look of a child's train. The fun design houses practical space to stow the power cord.

Features

- Whimsical pediatric appearance
- Dry piston compressor
- Comes with fun stickers and colorful lights

Warranty

3-year against manufacturer defects

| | |
|---|--|
| METHOD OF NEBULIZATION | dry piston compressor |
| OPERATION VOLTAGE | 110 volts ~ 60Hz |
| CLASS | Class II |
| POWER CONSUPTION | 1.2A |
| SOUND NOISE LEVEL (at distance of 1m) | 52 dB |
| JET NEBULIZER MAXIMUM CAPACITY | 10 ml |
| PRODUCT WEIGHT | 3.8 lbs. (1.73 kg) |
| UNIT DIMENSIONS | 9.4" x 5.5" x 4.9" 240 x 140 x 124 mm |
| MEAN MASS AERODYNAMIC DIAMETER | less than 5 microns |
| MEDICATION OUTLET | 0.33 ml/min. or greater* |
| OPERATION HUMIDITY AND TEMPERATURE | 41° F – 104° F, 30% to 95% |
| STORAGE HUMIDITY AND TEMPERATURE | 33.8° F – 104° F, 10% to 95% |
| ATMOSPHERIC PRESSURE RANGE FOR CARRIAGE AND STORAGE | 50 – 106 kPa, (500 –1,060 mbar) |

*Value measured with 0.9% NaCl solution. This specification may vary with the kind of the medication being used and it's viscosity.

The pediatric nebulizer includes:

- 1 jet nebulizer
- 1 PVC tubing
- 1 bag with 3 air filters
- 1 mouthpiece
- 1 pediatric mask
- 1 adult mask
- 1 instruction manual

Distributed by: **CLIAwaived.com™**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com