

CLIAwaived,™ Inc.

AUTOMATED URINALYSIS TEST SYSTEM

A semi-automated analyzer for urine chemistry analysis, this device reads urine test strips, stores results and prints reports without the need for special training. Using the latest in scientific technology, it can test up to 120 strips per hour.

- Accurate, reliable and convenient
- Automatic self-calibration for easy operation
- Optional barcode reader to record patient ID



PROMO OFFER ONLY \$499*

1 Urine Strip analyzer + 200 free tests
***Price offer expires 3/31/2015**

AUTOMATED URINALYSIS TEST SYSTEM

Catalog #

Analyzer	CLIA-U111-101
Analyzer with Barcode Reader	CLIA-U111-111
<i>Barcode Reader</i>	CLIA-U221-111
<i>Barcode Reader with Serial Splitter Cable</i>	CLIA-U221-113
<i>Thermal Printer Paper Rolls</i>	CLIA-U121-101
<i>Data Transfer Kit</i>	CLIA-U121-131
Auto Urinalysis Reagent Test Strip (100 strips)	CLIA-U111-103-111

The 10SG strip is an easy-to-read, convenient dip-and-read test. It is compatible with the CLIAwaived,™ Inc. Automated Urinalysis Test System.

Tests for glucose, bilirubin, ketone, specific gravity, blood, pH, protein, urobilinogen, leukocytes and nitrite.