


The On-The-Spot Leader in Saliva Alcohol Testing



Q.E.D.[®]
SALIVA ALCOHOL TEST

- Simple
- Accurate
- Reliable
- FDA Cleared
- DOT-Approved
- CLIA-Waived

The Q.E.D. Saliva Alcohol Test is the “on-the-spot leader” in the rapid, accurate and quantitative determination of alcohol in saliva. With a high correlation to blood-alcohol levels, Q.E.D. has a testing range of 0-145 mg/dL and provides quantitative results in minutes.

The Value of a Reliable and Trusted Brand




OraSure Technologies

Q.E.D.[®] Rapid Saliva Alcohol Testing


3 Easy Steps to Quantitative Results




1. Swab mouth to collect saliva.



2. Insert collector into test.



3. Read color bar.



On-Site Simplicity with Clinical Lab Accuracy . . .

The Q.E.D. test can be used as an accurate screen in place of blood tests. In clinical trials, saliva alcohol levels measured with the Q.E.D. tests demonstrated high correlation to blood analyzed by gas chromatography ($r=0.98$)¹ and produced no false positive or false negative results².

Data on file. ¹Study conducted by Thomas Jefferson University Hospital, Department of Emergency Medicine.

²Study conducted by Rutgers University, Center of Alcohol Studies.

It's Reliable, Fast, Simple . . .

- Non-invasive
- Quantitative results in minutes
- Reads like a thermometer
- No instrumentation
- Long shelf life
- Built-in QA Spot™ assures quality control for every test

The Q.E.D. test is used for . . .

DOT and Non-DOT Worksite Testing

- Random testing
- Post-accident
- Reasonable cause
- Follow-up testing

Substance Abuse and Criminal Justice Testing

- Testing upon admission
- Interventions
- Counseling sessions
- Random testing

Emergency/Trauma Departments

- Pre-discharge
- Admission screening
- Routine trauma treatment protocol
- Seizure or unconscious victims

Clinical Trials

- Compliance testing for alcohol use

**For More Product Information, Technical Support or
DOT Screening Test Technician Certification Information, Please Call**

1-800-672-7873



OraSure Technologies

220 East First Street • Bethlehem, PA 18015 • www.orasure.com • 1-800-ORASURE

©2011 OraSure Technologies, Inc., Q.E.D.[®] is a registered trademark of OraSure Technologies, Inc.

QED0001 (Rev 02/11)