

1. Identification

Product identifier	OSOM® hCG Serum Control - Negative
Other means of identification	
Kit number	138
Synonyms	hCG Serum Negative Control
Recommended use	Component of OSOM® hCG Serum Control Set (Catalog # 138). For use in quality control testing with the OSOM® hCG Combo Test. For In Vitro Diagnostic use only.
Recommended restrictions	Use in accordance with supplier's recommendations.
Manufacturer/Importer/Supplier/Distributor information	
Corporate Headquarters	Sekisui Diagnostics, LLC 4 Hartwell Place, Lexington, MA 02421, USA www.sekisuidiagnostics.com Phone: 800-332-1042 Americas 1-760-476-3962
Emergency Telephone Numbers	Europe, Middle East & Africa +1-760-476-3961 Asia Pacific +1-760-476-3960 Access code 333512

2. Hazard(s) identification

Physical hazards	Not classified.
Health Hazards	Not classified.
OSHA defined hazards	Not classified.
Label elements	
Hazard symbol	None.
Signal word	None.
Hazard statement	The mixture does not meet the criteria for classification.
Precautionary statement	
Prevention	None.
Response	None.
Storage	None.
Disposal	None.
Hazard(s) not otherwise classified (HNOC)	None known.

3. Composition/information on ingredients

Mixtures

Chemical name	CAS number	%
Sodium azide	26628-22-8	< 0.2
Horse Serum	N/A	~99.8

Composition comments	All concentrations are in percent by weight unless ingredient is a gas. Gas concentrations are in percent by volume.
-----------------------------	--

4. First-aid measures

Inhalation	Move to fresh air. Call a physician if symptoms develop or persist.
Skin contact	Wash skin thoroughly with soap and water. Get medical attention if irritation develops and persists.
Eye contact	In case of contact, immediately flush eyes with fresh water for at least 15 minutes while holding the eyelids open. Remove contact lenses if worn. Get medical attention if irritation persists.

Ingestion	Rinse mouth thoroughly with water. Never give anything by mouth to an unconscious person. Get medical attention if irritation develops and persists.
Most important symptoms/effects, acute and delayed	Ingestion may cause irritation and malaise.
Indication of immediate medical attention and special treatment needed	Provide general supportive measures and treat symptomatically.
General information	Ensure that medical personnel are aware of the material(s) involved, and take precautions to protect themselves.

5. Fire-fighting measures

Suitable extinguishing media	Extinguish with water spray, carbon dioxide, dry chemical or material appropriate for the surrounding fire.
Unsuitable extinguishing media	None known.
Specific hazards arising from the chemical	When heated to decomposition, may produce hydrazoic acid fumes.
Special protective equipment and precautions for firefighters	Self-contained breathing apparatus and full protective clothing must be worn in case of fire.
Fire fighting equipment/instructions	Use standard firefighting procedures and consider the hazards of other involved materials.
General fire hazards	The product is not flammable.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures	Keep unnecessary personnel away. Wear appropriate protective equipment and clothing during clean-up. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing.
Methods and materials for containment and cleaning up	Absorb spill with vermiculite or other inert material. Dispose of waste in accordance with all applicable federal, state, local and provincial environmental regulations, per Section 13.
Environmental precautions	Do not allow to enter drains, sewers or watercourses. This mixture contains a small amount of sodium azide which can react with copper, lead, brass or solder in plumbing systems and form potentially explosive metal azides. Follow proper disposal procedures.

7. Handling and storage

Precautions for safe handling	Avoid contact with skin and eyes. Wash thoroughly after handling. In case of insufficient ventilation, wear suitable respiratory equipment. Handle and open container with care. Observe good industrial hygiene practices.
Conditions for safe storage, including any incompatibilities	Store at 2-8°C (35-46°F). Store in a closed container away from incompatible materials.

8. Exposure controls/personal protection

Occupational exposure limits

US. ACGIH Threshold Limit Values

Components	Type	Value
Sodium azide (CAS 26628-22-8)	Ceiling	0.29 mg/m3
		0.11 ppm

US. NIOSH: Pocket Guide to Chemical Hazards

Components	Type	Value
Sodium azide (CAS 26628-22-8)	Ceiling	0.3 mg/m3
		0.1 ppm

Biological limit values No biological exposure limits noted for the ingredient(s).

Exposure guidelines Follow standard monitoring procedures.

US - California OELs: Skin designation

Sodium azide (CAS 26628-22-8) Can be absorbed through the skin.

US - Tennessee OELs: Skin designation

Sodium azide (CAS 26628-22-8)

Can be absorbed through the skin.

US. NIOSH: Pocket Guide to Chemical Hazards

Sodium azide (CAS 26628-22-8)

Can be absorbed through the skin.

Appropriate engineering controls

Provide adequate ventilation. Observe Occupational Exposure Limits and minimize the risk of inhalation of vapors.

Individual protection measures, such as personal protective equipment**Eye/face protection**

Wear approved safety glasses or goggles.

Skin protection**Hand protection**

Wear appropriate chemical resistant gloves.

Other

Remove contaminated clothing promptly.

Respiratory protection

In case of inadequate ventilation or risk of inhalation of vapors, use suitable respiratory equipment.

Thermal hazards

Wear appropriate thermal protective clothing, when necessary.

General hygiene considerations

Handle in accordance with good industrial hygiene and safety practice.

9. Physical and chemical properties**Appearance**

Clear, yellowish liquid.

Physical state

Liquid.

Form

Liquid.

Color

Yellowish.

Odor

Not available.

Odor threshold

Not available.

pH

8.1 Approximate

Melting point/freezing point

Not available.

Initial boiling point and boiling range

Not available.

Flash point

Not available.

Evaporation rate

Not available.

Flammability (solid, gas)

Not applicable.

Upper/lower flammability or explosive limits**Flammability limit - lower (%)**

Not available.

Flammability limit - upper (%)

Not available.

Explosive limit - lower (%)

Not available.

Explosive limit - upper (%)

Not available.

Vapor pressure

Not available.

Vapor density

Not available.

Relative density

1.06

Solubility(ies)**Solubility (water)**

Soluble.

Partition coefficient (n-octanol/water)

Not applicable.

Auto-ignition temperature

Not available.

Decomposition temperature

Not available.

Viscosity

Not available.

10. Stability and reactivity**Reactivity**

The product is stable and non-reactive under normal conditions of use, storage and transport.

Chemical stability

Material is stable under normal conditions.

Possibility of hazardous reactions

Contact with acids liberates toxic gas.

Conditions to avoid	Heat, sparks, flames, elevated temperatures.
Incompatible materials	Strong oxidizing agents.
Hazardous decomposition products	No data available.

11. Toxicological information

Information on likely routes of exposure

Inhalation	Vapors may irritate throat and respiratory system and cause coughing.
Skin contact	May cause skin irritation.
Eye contact	May cause eye irritation.
Ingestion	May cause discomfort if swallowed.

Symptoms related to the physical, chemical and toxicological characteristics
Ingestion may cause irritation and malaise.

Information on toxicological effects

Acute toxicity May cause discomfort if swallowed.

Components	Species	Test Results
Sodium azide (CAS 26628-22-8)		
Acute		
<i>Dermal</i>		
LD50	Rabbit	20 mg/kg
<i>Oral</i>		
LD50	Rat	27 mg/kg

Skin corrosion/irritation May cause skin irritation.

Serious eye damage/eye irritation May cause eye irritation.

Respiratory or skin sensitization

Respiratory sensitization Not classified.

Skin sensitization Not classified.

Germ cell mutagenicity Not classified.

Carcinogenicity Not classifiable as to carcinogenicity to humans.

OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

Not listed.

Reproductive toxicity Not classified.

Specific target organ toxicity - single exposure Not classified.

Specific target organ toxicity - repeated exposure Not classified.

Aspiration hazard Not classified.

Chronic effects No data available.

Further information Contains animal source material of U.S. origin. The human hCG in this preparation was tested by FDA-approved methods and found to be negative for the presence of hepatitis B virus surface antigen (HBsAg), human immunodeficiency virus (HIV) 1 & 2 and hepatitis C virus (HCV). However, because no test method can provide complete assurance that infectious agents are absent, this product should be handled as a potentially biohazardous material in accordance with universal/standard precautions.

12. Ecological information

Ecotoxicity

Components	Species	Test Results
Sodium azide (CAS 26628-22-8)		
Aquatic		
Algae	EC50 Pseudokirchnerella subcapitata	0.35 mg/l, 96 hours

Components	Species	Test Results
Fish	LC50	5.7 mg/l, 96 hours
Persistence and degradability	No data is available on the degradability of this product.	
Bioaccumulative potential	Not available.	
Mobility in soil	Not available.	
Mobility in general	The product is soluble in water.	
Other adverse effects	The product is not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.	

13. Disposal considerations

Disposal instructions	Collect and reclaim or dispose in sealed containers at licensed waste disposal site. Contaminated instruments and surfaces should be disinfected in accordance with your employer's chemical-specific and universal/standard precautions. This preparation contains a small amount of sodium azide which can react with copper, lead, brass or solder in plumbing systems and form potentially explosive metal azides. If preparation enters drain, flush with a large volume of water to prevent azide build-up.
Hazardous waste code	The Waste code should be assigned in discussion between the user, the producer and the waste disposal company.
Waste from residues / unused products	Dispose in accordance with all applicable regulations. Contract with a licensed chemical disposal agency.
Contaminated packaging	Empty containers should be taken to an approved waste handling site for recycling or disposal.

14. Transport information

DOT

Not regulated as dangerous goods.

IATA

Not regulated as dangerous goods.

IMDG

Not regulated as dangerous goods.

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code This substance/mixture is not intended to be transported in bulk.

15. Regulatory information

US federal regulations This product is not hazardous according to OSHA 29CFR 1910.1200. This mixture is a component of an in vitro diagnostic device regulated by the U.S. Food and Drug Administration.

TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D)

Not regulated.

OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

Not listed.

CERCLA Hazardous Substance List (40 CFR 302.4)

Sodium azide (CAS 26628-22-8) LISTED

Superfund Amendments and Reauthorization Act of 1986 (SARA)

Hazard categories Immediate Hazard - No
Delayed Hazard - No
Fire Hazard - No
Pressure Hazard - No
Reactivity Hazard - No

SARA 302 Extremely hazardous substance

Chemical name	CAS number	Reportable quantity	Threshold planning quantity	Threshold planning quantity, lower value	Threshold planning quantity, upper value
---------------	------------	---------------------	-----------------------------	--	--

Sodium azide	26628-22-8	1000	500 lbs		
--------------	------------	------	---------	--	--

SARA 311/312 Hazardous chemical No

SARA 313 (TRI reporting)

Not regulated.

Other federal regulations**Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List**

Not regulated.

Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130)

Not regulated.

Safe Drinking Water Act (SDWA)

Not regulated.

US state regulations

This product does not contain a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

US. Massachusetts RTK - Substance List

Sodium azide (CAS 26628-22-8)

US. New Jersey Worker and Community Right-to-Know Act

Sodium azide (CAS 26628-22-8)

US. Pennsylvania Worker and Community Right-to-Know Law

Sodium azide (CAS 26628-22-8)

US. Rhode Island RTK

Sodium azide (CAS 26628-22-8)

US. California Proposition 65**US - California Proposition 65 - Carcinogens & Reproductive Toxicity (CRT): Listed substance**


Not listed.

International Inventories

Country(s) or region	Inventory name	On inventory (yes/no)*
Australia	Australian Inventory of Chemical Substances (AICS)	Yes
Canada	Domestic Substances List (DSL)	Yes
Canada	Non-Domestic Substances List (NDSL)	No
China	Inventory of Existing Chemical Substances in China (IECSC)	Yes
Europe	European Inventory of Existing Commercial Chemical Substances (EINECS)	Yes
Europe	European List of Notified Chemical Substances (ELINCS)	No
Japan	Inventory of Existing and New Chemical Substances (ENCS)	Yes
Korea	Existing Chemicals List (ECL)	Yes
New Zealand	New Zealand Inventory	Yes
Philippines	Philippine Inventory of Chemicals and Chemical Substances (PICCS)	Yes
United States & Puerto Rico	Toxic Substances Control Act (TSCA) Inventory	Yes

*A "Yes" indicates this product complies with the inventory requirements administered by the governing country(s).

A "No" indicates that one or more components of the product are not listed or exempt from listing on the inventory administered by the governing country(s).

16. Other information, including date of preparation or last revision**Issue date** 13-October-2014**Revision date** -**Version #** 01**NFPA ratings****List of abbreviations**

LD50: Lethal Dose, 50%.

References

ACGIH

HSDB® - Hazardous Substances Data Bank

IARC Monographs. Overall Evaluation of Carcinogenicity

Disclaimer

The information above is provided in good faith. It is believed to be accurate and represents the best information currently available to us. HOWEVER, WE MAKE NO WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OF ANY OTHER TYPE, EXPRESSED OR IMPLIED, WITH RESPECT TO PRODUCTS DESCRIBED OR DATA OR INFORMATION PROVIDED, AND WE ASSUME NO LIABILITY RESULTING FROM THE USE OF SUCH PRODUCTS, DATA OR INFORMATION. Users should make their own investigations to determine the suitability of the information for their particular purposes, and the user assumes all risk arising from their use of the material. The user is required to comply with all laws and regulations relating to the purchase, use, storage and disposal of the material, and must be familiar with and follow generally accepted safe handling procedures. In no event shall Sekisui Diagnostics be liable for any claims, losses, or damages of any individual or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if Sekisui Diagnostics has been advised of the possibility of such damages.

1. Identification

Product identifier	OSOM® hCG Serum Control - Positive
Other means of identification	
Kit number	138
Synonyms	hCG Serum Positive Control
Recommended use	Component of OSOM® hCG Serum Control Set (Catalog # 138). For use in quality control testing with the OSOM® hCG Combo Test. For In Vitro Diagnostic use only.
Recommended restrictions	Use in accordance with supplier's recommendations.
Manufacturer/Importer/Supplier/Distributor information	
Corporate Headquarters	Sekisui Diagnostics, LLC 4 Hartwell Place, Lexington, MA 02421, USA www.sekisuidiagnostics.com Phone: 800-332-1042 Americas 1-760-476-3962
Emergency Telephone Numbers	Europe, Middle East & Africa +1-760-476-3961 Asia Pacific +1-760-476-3960 Access code 333512

2. Hazard(s) identification

Physical hazards	Not classified.
Health Hazards	Not classified.
OSHA defined hazards	Not classified.
Label elements	
Hazard symbol	None.
Signal word	None.
Hazard statement	The mixture does not meet the criteria for classification.
Precautionary statement	
Prevention	None.
Response	None.
Storage	None.
Disposal	None.
Hazard(s) not otherwise classified (HNOC)	None known.

3. Composition/information on ingredients

Mixtures

Chemical name	CAS number	%
Sodium azide	26628-22-8	< 0.2
Horse Serum	N/A	~99.8

Composition comments	All concentrations are in percent by weight unless ingredient is a gas. Gas concentrations are in percent by volume.
-----------------------------	--

4. First-aid measures

Inhalation	Move to fresh air. Call a physician if symptoms develop or persist.
Skin contact	Wash skin thoroughly with soap and water. Get medical attention if irritation develops and persists.
Eye contact	In case of contact, immediately flush eyes with fresh water for at least 15 minutes while holding the eyelids open. Remove contact lenses if worn. Get medical attention if irritation persists.

Ingestion	Rinse mouth thoroughly with water. Never give anything by mouth to an unconscious person. Get medical attention if irritation develops and persists.
Most important symptoms/effects, acute and delayed	Ingestion may cause irritation and malaise.
Indication of immediate medical attention and special treatment needed	Provide general supportive measures and treat symptomatically.
General information	Ensure that medical personnel are aware of the material(s) involved, and take precautions to protect themselves.

5. Fire-fighting measures

Suitable extinguishing media	Extinguish with water spray, carbon dioxide, dry chemical or material appropriate for the surrounding fire.
Unsuitable extinguishing media	None known.
Specific hazards arising from the chemical	When heated to decomposition, may produce hydrazoic acid fumes.
Special protective equipment and precautions for firefighters	Self-contained breathing apparatus and full protective clothing must be worn in case of fire.
Fire fighting equipment/instructions	Use standard firefighting procedures and consider the hazards of other involved materials.
General fire hazards	The product is not flammable.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures	Keep unnecessary personnel away. Wear appropriate protective equipment and clothing during clean-up. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing.
Methods and materials for containment and cleaning up	Absorb spill with vermiculite or other inert material. Dispose of waste in accordance with all applicable federal, state, local and provincial environmental regulations, per Section 13.
Environmental precautions	Do not allow to enter drains, sewers or watercourses. This mixture contains a small amount of sodium azide which can react with copper, lead, brass or solder in plumbing systems and form potentially explosive metal azides. Follow proper disposal procedures.

7. Handling and storage

Precautions for safe handling	Avoid contact with skin and eyes. Wash thoroughly after handling. In case of insufficient ventilation, wear suitable respiratory equipment. Handle and open container with care. Observe good industrial hygiene practices.
Conditions for safe storage, including any incompatibilities	Store at 2-8°C (35-46°F). Store in a closed container away from incompatible materials.

8. Exposure controls/personal protection

Occupational exposure limits

US. ACGIH Threshold Limit Values

Components	Type	Value
Sodium azide (CAS 26628-22-8)	Ceiling	0.29 mg/m3
		0.11 ppm

US. NIOSH: Pocket Guide to Chemical Hazards

Components	Type	Value
Sodium azide (CAS 26628-22-8)	Ceiling	0.3 mg/m3
		0.1 ppm

Biological limit values No biological exposure limits noted for the ingredient(s).

Exposure guidelines Follow standard monitoring procedures.

US - California OELs: Skin designation

Sodium azide (CAS 26628-22-8) Can be absorbed through the skin.

US - Tennessee OELs: Skin designation

Sodium azide (CAS 26628-22-8)

Can be absorbed through the skin.

US. NIOSH: Pocket Guide to Chemical Hazards

Sodium azide (CAS 26628-22-8)

Can be absorbed through the skin.

Appropriate engineering controls

Provide adequate ventilation. Observe Occupational Exposure Limits and minimize the risk of inhalation of vapors.

Individual protection measures, such as personal protective equipment**Eye/face protection**

Wear approved safety glasses or goggles.

Skin protection**Hand protection**

Wear appropriate chemical resistant gloves.

Other

Remove contaminated clothing promptly.

Respiratory protection

In case of inadequate ventilation or risk of inhalation of vapors, use suitable respiratory equipment.

Thermal hazards

Wear appropriate thermal protective clothing, when necessary.

General hygiene considerations

Handle in accordance with good industrial hygiene and safety practice.

9. Physical and chemical properties**Appearance**

Clear, yellowish liquid.

Physical state

Liquid.

Form

Liquid.

Color

Yellowish.

Odor

Not available.

Odor threshold

Not available.

pH

8.1 Approximate

Melting point/freezing point

Not available.

Initial boiling point and boiling range

Not available.

Flash point

Not available.

Evaporation rate

Not available.

Flammability (solid, gas)

Not applicable.

Upper/lower flammability or explosive limits**Flammability limit - lower (%)**

Not available.

Flammability limit - upper (%)

Not available.

Explosive limit - lower (%)

Not available.

Explosive limit - upper (%)

Not available.

Vapor pressure

Not available.

Vapor density

Not available.

Relative density

1.06

Solubility(ies)**Solubility (water)**

Soluble.

Partition coefficient (n-octanol/water)

Not applicable.

Auto-ignition temperature

Not available.

Decomposition temperature

Not available.

Viscosity

Not available.

10. Stability and reactivity**Reactivity**

The product is stable and non-reactive under normal conditions of use, storage and transport.

Chemical stability

Material is stable under normal conditions.

Possibility of hazardous reactions

Contact with acids liberates toxic gas.

Conditions to avoid	Heat, sparks, flames, elevated temperatures.
Incompatible materials	Strong oxidizing agents.
Hazardous decomposition products	No data available.

11. Toxicological information

Information on likely routes of exposure

Inhalation	Vapors may irritate throat and respiratory system and cause coughing.
Skin contact	May cause skin irritation.
Eye contact	May cause eye irritation.
Ingestion	May cause discomfort if swallowed.

Symptoms related to the physical, chemical and toxicological characteristics
Ingestion may cause irritation and malaise.

Information on toxicological effects

Acute toxicity May cause discomfort if swallowed.

Components	Species	Test Results
Sodium azide (CAS 26628-22-8)		
Acute		
<i>Dermal</i>		
LD50	Rabbit	20 mg/kg
<i>Oral</i>		
LD50	Rat	27 mg/kg

Skin corrosion/irritation May cause skin irritation.

Serious eye damage/eye irritation May cause eye irritation.

Respiratory or skin sensitization

Respiratory sensitization Not classified.

Skin sensitization Not classified.

Germ cell mutagenicity Not classified.

Carcinogenicity Not classifiable as to carcinogenicity to humans.

OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

Not listed.

Reproductive toxicity Not classified.

Specific target organ toxicity - single exposure Not classified.

Specific target organ toxicity - repeated exposure Not classified.

Aspiration hazard Not classified.

Chronic effects No data available.

Further information Contains animal source material of U.S. origin. The human hCG in this preparation was tested by FDA-approved methods and found to be negative for the presence of hepatitis B virus surface antigen (HBsAg), human immunodeficiency virus (HIV) 1 & 2 and hepatitis C virus (HCV). However, because no test method can provide complete assurance that infectious agents are absent, this product should be handled as a potentially biohazardous material in accordance with universal/standard precautions.

12. Ecological information

Ecotoxicity

Components	Species	Test Results
Sodium azide (CAS 26628-22-8)		
Aquatic		
Algae	EC50 Pseudokirchnerella subcapitata	0.35 mg/l, 96 hours

Components	Species	Test Results
Fish	LC50	5.7 mg/l, 96 hours
Persistence and degradability	No data is available on the degradability of this product.	
Bioaccumulative potential	Not available.	
Mobility in soil	Not available.	
Mobility in general	The product is soluble in water.	
Other adverse effects	The product is not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.	

13. Disposal considerations

Disposal instructions	Collect and reclaim or dispose in sealed containers at licensed waste disposal site. Contaminated instruments and surfaces should be disinfected in accordance with your employer's chemical-specific and universal/standard precautions. This preparation contains a small amount of sodium azide which can react with copper, lead, brass or solder in plumbing systems and form potentially explosive metal azides. If preparation enters drain, flush with a large volume of water to prevent azide build-up.
Hazardous waste code	The Waste code should be assigned in discussion between the user, the producer and the waste disposal company.
Waste from residues / unused products	Dispose in accordance with all applicable regulations. Contract with a licensed chemical disposal agency.
Contaminated packaging	Empty containers should be taken to an approved waste handling site for recycling or disposal.

14. Transport information

DOT	Not regulated as dangerous goods.
IATA	Not regulated as dangerous goods.
IMDG	Not regulated as dangerous goods.
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	This substance/mixture is not intended to be transported in bulk.

15. Regulatory information

US federal regulations		This product is not hazardous according to OSHA 29CFR 1910.1200. This mixture is a component of an in vitro diagnostic device regulated by the U.S. Food and Drug Administration.			
TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D)					
Not regulated.					
OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)					
Not listed.					
CERCLA Hazardous Substance List (40 CFR 302.4)					
Sodium azide (CAS 26628-22-8)		LISTED			
Superfund Amendments and Reauthorization Act of 1986 (SARA)					
Hazard categories		Immediate Hazard - No Delayed Hazard - No Fire Hazard - No Pressure Hazard - No Reactivity Hazard - No			
SARA 302 Extremely hazardous substance					
Chemical name	CAS number	Reportable quantity	Threshold planning quantity	Threshold planning quantity, lower value	Threshold planning quantity, upper value
Sodium azide	26628-22-8	1000	500 lbs		
SARA 311/312 Hazardous chemical	No				

SARA 313 (TRI reporting)

Not regulated.

Other federal regulations**Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List**

Not regulated.

Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130)

Not regulated.

Safe Drinking Water Act (SDWA)

Not regulated.

US state regulations

This product does not contain a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

US. Massachusetts RTK - Substance List

Sodium azide (CAS 26628-22-8)

US. New Jersey Worker and Community Right-to-Know Act

Sodium azide (CAS 26628-22-8)

US. Pennsylvania Worker and Community Right-to-Know Law

Sodium azide (CAS 26628-22-8)

US. Rhode Island RTK

Sodium azide (CAS 26628-22-8)

US. California Proposition 65**US - California Proposition 65 - Carcinogens & Reproductive Toxicity (CRT): Listed substance**


Not listed.

International Inventories

Country(s) or region	Inventory name	On inventory (yes/no)*
Australia	Australian Inventory of Chemical Substances (AICS)	Yes
Canada	Domestic Substances List (DSL)	Yes
Canada	Non-Domestic Substances List (NDSL)	No
China	Inventory of Existing Chemical Substances in China (IECSC)	Yes
Europe	European Inventory of Existing Commercial Chemical Substances (EINECS)	Yes
Europe	European List of Notified Chemical Substances (ELINCS)	No
Japan	Inventory of Existing and New Chemical Substances (ENCS)	Yes
Korea	Existing Chemicals List (ECL)	Yes
New Zealand	New Zealand Inventory	Yes
Philippines	Philippine Inventory of Chemicals and Chemical Substances (PICCS)	Yes
United States & Puerto Rico	Toxic Substances Control Act (TSCA) Inventory	Yes

*A "Yes" indicates this product complies with the inventory requirements administered by the governing country(s).

A "No" indicates that one or more components of the product are not listed or exempt from listing on the inventory administered by the governing country(s).

16. Other information, including date of preparation or last revision**Issue date** 13-October-2014**Revision date** -**Version #** 01**NFPA ratings****List of abbreviations**

LD50: Lethal Dose, 50%.

References

ACGIH

HSDB® - Hazardous Substances Data Bank

IARC Monographs. Overall Evaluation of Carcinogenicity

Disclaimer

The information above is provided in good faith. It is believed to be accurate and represents the best information currently available to us. HOWEVER, WE MAKE NO WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OF ANY OTHER TYPE, EXPRESSED OR IMPLIED, WITH RESPECT TO PRODUCTS DESCRIBED OR DATA OR INFORMATION PROVIDED, AND WE ASSUME NO LIABILITY RESULTING FROM THE USE OF SUCH PRODUCTS, DATA OR INFORMATION. Users should make their own investigations to determine the suitability of the information for their particular purposes, and the user assumes all risk arising from their use of the material. The user is required to comply with all laws and regulations relating to the purchase, use, storage and disposal of the material, and must be familiar with and follow generally accepted safe handling procedures. In no event shall Sekisui Diagnostics be liable for any claims, losses, or damages of any individual or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if Sekisui Diagnostics has been advised of the possibility of such damages.