


# RAPID NICOTINE TEST CARD



The Rapid Nicotine Test Card is a lateral flow immunochromatographic assay for the qualitative detection of Nicotine metabolite (cotinine) in human urine with 3 cutoff levels (50, 200 and 500 ng/mL).

The Rapid Nicotine Test Card is intended for prescription professional, in-vitro diagnostic use in near-patient (POC), centralized testing locations and for testing in employment settings. The device can consist of any combination of the 3 assays in 1, 2 or 3-test formats.

The test provides a preliminary result only; preliminary positive results should be confirmed using an alternate chemical methodology (GC/MS, GC/MS/MS, LC/MS or LC/MS/MS) if the donor does not admit to nicotine use or anytime your test policies require.

To order, call 888.882.7739