



Certificate of Analysis

Product: BDH Isopropyl Alcohol
Cat: BDH1133-5GL
Lot: 040413B
Date of Manufacture: 04/08/13
Use Before: 10/07/14
Meets ACS Specifications for General Use

<u>Test</u>	<u>Specification</u>		<u>Result</u>	<u>Units</u>
	<u>Min.</u>	<u>Max.</u>		
Water		0.2	0.037	%
Residue after Evaporation		0.001	<0.0001	%
Assay by GC Analysis	99.50		99.93	%
Solubility in Water	Passes test		Passes test	
Color (APHA)		10	<10	
Titration Acid or Base		0.1	0.04	μeq/g
Carbonyl Compounds (as acetone or propionaldehyde)		0.002	<0.002	%

Quality Control Approval

Distributed Exclusively by

