

Nano-Check™ AMI Cardiac Panel Test

Distributed by

CLIAwaived™ Inc.



The Nano-Check™ AMI Cardiac Panel Test is a lateral flow immunoassay for three key cardiac markers (TnI, CK-MB, Myo) in plasma, serum or whole blood. The Nano-Check™ AMI Cardiac Panel Test aids in the early detection of Acute Myocardial Infarction (AMI) in hospitals and physicians' offices. The results can be read visually for qualitative results, or quantitative results can be obtained in conjunction with the Nano-Checker 710 Reader. Quick access to cardiac marker status allows healthcare professionals to quickly initiate treatment to patients experiencing chest pain.

BENEFITS:

- Simple one-step test
- Results in 10-15 minutes
- Intuitive results display
- Fast turn-around
- Easy-to-read controls
- Provides accurate quantitative results
- Assures cost-effective and timely patient care

Ordering Information:


Product Description	Catalog #
Nano-Check™ AMI Cardiac Panel Test - 3 in 1 Whole Blood/Plasma/Serum (20 cassettes)	ND302P
Nano-Check™ AMI Cardiac Panel Test - 2 in 1 Whole Blood/Plasma/Serum (20 cassettes)	ND202P
Nano-Check™ AMI Cardiac Panel Test - cTnI Whole Blood/Plasma/Serum (20 cassettes)	ND7140

Reimbursement Information:


CPT Code #	National Limit Amount
N/A	N/A

All product information and prices are subject to change without notice at any time. Please call our office before ordering to verify current prices.

Nano-Check™ AMI Cardiac Panel Test Instructions



Intuitive results display - visual signal indicates positive; any faint signal on test line, and negative; no signal.



Quantitative results by Nano-Checker 710 are automatically determined and displayed or printed as PDF format report.