



T-Cup[®] Multi-Drug Urine Test Cup



- **CLIA Waived, 510(K) Cleared & Forensic Use Only**
- **Female User Friendly & Customizable**
- **Test up to 18 Drugs per Cup**

CLIA Waived Drugs & Cutoff Levels

AMP AMP AMP AMP AMP	Amphetamine (AMP)	1000ng/mL 500ng/mL	TCA TCA TCA TCA	Tricyclic Antidepressants (TCA)	1000ng/mL
BAR BAR BAR BAR BAR	Barbiturates (BAR)	300ng/mL	MET MET MET MET	Methamphetamine (MET/mAMP)	1000ng/mL 500ng/mL
BZO BZO BZO BZO BZO	Benzodiazepines (BZO)	300ng/mL	MDMA MDMA MDMA MDMA	Methylenedioxymethamphetamine (MDMA)	500ng/mL
BUP BUP BUP BUP BUP	Buprenorphine (BUP)	10ng/mL	MOP MOP MOP MOP	Morphine (MOP/OPI300)	300ng/mL
COC COC COC COC COC	Cocaine (COC)	300ng/mL 150ng/mL	OPI OPI OPI OPI OPI	Opiate (OPI)	2000ng/mL
THC THC THC THC THC	Cannabinoids (THC)	50ng/mL	OXY OXY OXY OXY OXY	Oxycodone (OXY)	100ng/mL
EDDP EDDP EDDP EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300ng/mL	PCP PCP PCP PCP PCP	Phencyclidine (PCP)	25ng/mL
MTD MTD MTD MTD MTD	Methadone (MTD)	300ng/mL	PPX PPX PPX PPX PPX	Propoxyphene (PPX)	300ng/mL

Forensic Use Drugs & Cutoff Levels

6-MAM 6-MAM 6-MAM 6-MAM	Monoacetylmorphine (6-MAM)	10ng/mL	K2 K2 K2 K2 K2	Synthetic Cannabinoid (K2)	50ng/mL
AMP AMP AMP AMP AMP	Amphetamine (AMP)	300ng/mL	KET KET KET KET KET	Ketamine (KET)	300ng/mL 1000ng/mL
BAR BAR BAR BAR BAR	Barbiturates (BAR)	200ng/mL	KRA KRA KRA KRA KRA	Kratom (KRA)	300ng/mL
BZO BZO BZO BZO BZO	Benzodiazepines (BZO)	100ng/mL 200ng/mL	MET MET MET MET MET	Methamphetamine (MET/mAMP)	300ng/mL
COC COC COC COC COC	Cocaine (COC)	100ng/mL	MOP MOP MOP MOP MOP	Morphine (MOP/OPI)	100ng/mL
COT COT COT COT COT	Cotinine (COT)	200ng/mL	MTD MTD MTD MTD MTD	Methadone (MTD)	200ng/mL
EDDP EDDP EDDP EDDP	EDDP	100ng/mL	MQL MQL MQL MQL MQL	Methaqualone (MQL)	300ng/mL
EtG EtG EtG EtG EtG	Ethyl Glucuronide (EtG)	500ng/mL	THC THC THC THC THC	Cannabinoids (THC)	25ng/mL 40ng/mL
FTY FTY FTY FTY FTY	Fentanyl (FTY)	20ng/mL	TRA TRA TRA TRA TRA	Tramadol (TRA)	100ng/mL 200ng/mL 1000ng/mL
GAB GAB GAB GAB GAB	Gabapentin (GAB)	1000ng/mL			

Built-in Adulteration Parameters: Creatinine (CR), Glutaraldehyde (GL), Specific Gravity (S.G.) & pH. Oxidant (OX) & Nitrite (NI) available upon request. Note: Generic stock product adulterants include CR,SG,pH. Other adulterations available through special order.

Order Information

Item Code	Description	Complexity
TDOA-154	5-Panel T-Cup; AMP,COC,OPI,PCP,THC	CLIA WAIVED
TDOA-254	5-Panel T-Cup; AMP,COC,mAMP,OPI,THC	CLIA WAIVED
TDOA-254A3	5-Panel T-Cup; AMP,COC,mAMP,OPI,THC with 3 AD OX,SG,pH	CLIA WAIVED
TDOA-654	5-Panel T-Cup; BZO,COC,mAMP,OPI,THC	CLIA WAIVED
TDOA-264	6-Panel T-Cup; AMP,BZO,COC,mAMP,OPI,THC	CLIA WAIVED
TDOA-264A3	6-Panel T-Cup; AMP,BZO,COC,mAMP,OPI,THC with 3 AD	CLIA WAIVED
TDOA-564	6-Panel T-Cup; BZO,COC,mAMP,OPI,OXY,THC	CLIA WAIVED
TDOA-274	7-Panel T-Cup; AMP,BZO,COC,mAMP,OPI,OXY,THC	CLIA WAIVED
TDOA-3104	10-Panel T-Cup; AMP,BAR,BZO,COC,mAMP,MDMA,MTD,OPI,PCP,THC	CLIA WAIVED
TDOA-3104A3	10-Panel T-Cup; AMP,BAR,BZO,COC,mAMP,MDMA,MTD,OPI,PCP,THC with 3 AD	CLIA WAIVED
TDOA-4104	10-Panel T-Cup; AMP,BAR,BZO,COC,mAMP,MTD,OPI,OXY,PCP,THC	CLIA WAIVED
TDOA-7104	10-Panel T-Cup; BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,TCA,THC	CLIA WAIVED
★ TDOA-7104A3	10-Panel T-Cup; BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,TCA,THC with 3 AD	CLIA WAIVED
TDOA-8104	10-Panel T-Cup; AMP,BUP,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,THC	CLIA WAIVED
TDOA-3124	12-Panel T-Cup; AMP,BAR,BUP,BZO,COC,MDMA,MTD,OPI,OXY,PCP,TCA,THC	CLIA WAIVED
TDOA-4124	12-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,PCP,PPX,THC	CLIA WAIVED
TDOA-6125	12-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC	CLIA WAIVED
TDOA-6125A3	12-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC with 3 AD	CLIA WAIVED
TDOA-6124	12-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,THC	CLIA WAIVED
TDOA-6124A3	12-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,THC with 3 AD	CLIA WAIVED
TDOA-7125	12-Panel T-Cup; AMP,BAR,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,TCA,THC	CLIA WAIVED
TDOA-7125A3	12-Panel T-Cup; AMP,BAR,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,TCA,THC with 3 AD	CLIA WAIVED
TDOA-8135E	13-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC,ETG	Forensic Use Only
TDOA-6135F	13-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC,FTY	Forensic Use Only
TDOA-7135A3EFT	13-Panel T-Cup; AMP,BUP,BZO,COC,mAMP,MDMA,MTD,MOP,OXY,THC,ETG,FTY,TRA with 3AD	Forensic Use Only
TDOA-7135EFT	13-Panel T-Cup; AMP,BUP,BZO,COC,mAMP,MDMA,MTD,MOP,OXY,THC,TRA,FTY,ETG	Forensic Use Only
★ TDOA-3135	13-Panel T-Cup; AMP500,BAR,BUP,BZO,COC150,mAMP500,MDMA,MOP,MTD,OXY,PCP,TCA,THC	CLIA WAIVED
TDOA-3135A3	13-Panel T-Cup; AMP500,BAR,BUP,BZO,COC150,mAMP500,MDMA,MOP,MTD,OXY,PCP,TCA,THC+OX,SG,pH	CLIA WAIVED
TDOA-2135	13-Panel T-Cup;AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,TCA,THC	CLIA WAIVED
TDOA-2135A3	13-Panel T-Cup;AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,TCA,THC with 3 AD	CLIA WAIVED
TDOA-1145A3	14-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,PPX,TCA,THC with 3 AD	CLIA WAIVED
TDOA-8145EF	14-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC,ETG,FTY	Forensic Use Only
TDOA-1144	14-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,PPX,TCA,THC	CLIA WAIVED
TDOA-1144A3	14-Panel T-Cup; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,PPX,TCA,THC with 3 AD	CLIA WAIVED
TDOA-9145EFTK	14-Panel T-Cup; AMP,BUP,BZO,COC,mAMP,MDMA,MTD,MOP,OXY,THC,TRA,FTY,ETG,K2	Forensic Use Only

★ Available through special orders